



Skil-Care™

CORPORATION

Positioning, Safety & Pressure Relief Solutions Since 1978

Healthcare Product Catalog



Wound Management	3
Cushions	14
Wheelchair Accessories	42
Bed Positioning	65
Bed Rail Pads	74
Mattresses	79
Gerl-Chair Accessories	83
Safety	86
Ambulation	104
Transfer Slings	108
Bariatric	109
Physical Restraints	116
Patient Care	122
Incontinent Supplies	124
Recreation/Exercise	126
Sensory Stimulation	130

Laundered Products

Fabrics Include:

(LP1) * **Low Shear I** - A two-way stretch, breathable tricot fabric that protects against skin damage caused by friction and shear force. Fluid-proof, it has a laminated polyurethane backing.

(LP2) * **Polyester** - A knitted fabric that stretches to reduce skin-damaging shear forces.

(LP3) * **Cozy Cloth** - An ultra soft polyester microfiber that wicks away moisture, reduces friction and shear forces.

(LP4) * **Nylon** - Incontinent-proof, durable and laminated with a polyurethane backing.

(LP5) * **Polyester Mesh** - A knit, open weave polyester that is cool to the skin.

(LP6) * **Super Soft Loop Material** - Loop material mainly used in our wound care boots.

(LP7) * **Synthetic Sheepskin** - A fluff polyester that simulates the characteristics of sheepskin. The product helps reduce pressure sores and allows air flow between the fiber.

Cleaning:

- * The fabrics above can be laundered at a temperatures below 180°F (82°C).
- * Use lower temperature during the dry cycle to extend product life.



Wipe Cleanable Products

Fabrics Include:

(WC1) * **Low Shear II** - A two-way stretch, breathable polyurethane fabric with a laminated polyester coating that reduces skin-damaging shear forces.

(WC2) * **Low Shear III** - A four-way stretch, breathable polyurethane fabric with a laminated polyester coating that reduces skin-damaging shear forces.

(WC3) * **Vinyl (UnSupported)** - A heavy duty vinyl that is used for all Skil-Care welded (heat-sealed) products. The material is waterproof and anti-microbial.

(WC4) * **Supported Vinyl (Sewn)** - A vinyl coated fabric that offers superior strength, durability, and cleanability for sewn products.

(WC5) * **Neoprene** - A soft, pliable rubber composite with a knit fabric backing that provides maximum skin protection from abrasions and bruising.

(WC6) * **Non-Skid Fabric** - A fabric coated with a non-skid urethane.

(WC7) * **Nylon** - Incontinent-proof, durable and laminated with a polyurethane backing.

(WC8) * **EVA Foam (EZ Landing Fall Mat)** - A one inch molded floor mat comprised of shock absorbing EVA foam.

(WC9) * **Fiberglass (Smokers Apron)** - A silicon coated woven fiberglass that is fire proof.

(WC10) * **Non-Skid Matting** - A polyester woven fabric with a dipped rubber coating.

Cleaning:

- * Wipe fabrics with either soap and water, liquid disinfectant or approved disinfecting wipes.
- * Rinse thoroughly to remove any residual chemicals or soap.
- * Make sure fabric is completely dry before using.
- * Disinfect blood contamination with a 1:10 dilution of household bleach (5.25% sodium hypochlorite) as recommended by the CDC (Center of Disease Control, US Department of Health and Human Services, February 1989).

Cautions:

- * DO NOT use cleaners that contain Phenol or Benzyl.
- * DO NOT use iodophor type disinfectants, such as Betadine, because they will stain the fabrics listed above.
- * DO NOT use household bleach or concentrated cleaning agents, disinfectants or germicides without diluting the product according to the manufacturer's instructions.
- * DO NOT allow any cleaning agent to remain in contact with the fabrics for an extended time

Foam Only Products

Fabrics Include:

(FO1) * **Vel Foam** - Disposable only

(FO2) * **UnCovered Foam** - At the users risk, foam products sometimes can be laundered at low temperatures (160°F, 71°C) when placed in a mesh laundry bag.

Cleaning:

- * At the users risk, foam products sometimes can be laundered at low temperatures (160°F, 71°C) when placed in a mesh laundry bag.

Table of Contents

<i>Description</i>	<i>Page No.</i>	<i>Description</i>	<i>Page No.</i>	<i>Description</i>	<i>Page No.</i>
Wound Management		Wheelchair Accessories (cont.)		Bariatric (cont.)	
Heelt Off-loading	3	Bags & Holders	60-62	Foot & Leg Supports	111
Heel & Elbow Protection	4-12	Belts	63-64	Accessories	112-113
Compression & Support		Bed Positioning		Cushions	114-115
Garments	13	Lateral Positioning/		Physical Restraints	
Cushions		Roll Control	66-68	Belts	117-118
Cushion Cover Chart	15	Patient Handling	69	Limb Holder	118
Multi-Purpose	16-23	Wedges & Support	70-73	Vests & Jackets	119
- Gel-Foam	16-19	Bed Rail Pads	74-78	Padded Mitts	120-121
- Foam Only	19-22	Mattresses	79-82	Patient Care	122-123
- Air	22-23	Geri-Chair Accessories	83-85	Incontinence	124-125
Slide Control	24-27	Safety		Recreation & Exercise	126-129
Thin-Line	28-30	ChairPro Safety		Sensory Stimulation	130-137
Firm Foundation	31-34	Alarm Sets	87-90	Product Name Index	138-141
Coccyx Cutout	35	BedPro Safety		Numerical Index	142-144
Bariatric	36	Alarm Sets	91-92		
Cushion Covers	37	MultiPro Safety Alarms	92-95		
Pillows	38	FloorPro Safety Alarms	96		
Reminder Cushions	39-41	Magnetic Alarms	97		
Wheelchair Accessories		Alarm Accessories	98-99		
Backrest	43-45	Resident Safety	100-101		
Foot & Leg Supports	52-55	Fall Mats	102-103		
Lap Top Cushion & Trays	46-47	Ambulation			
Lateral Supports	48-51	Transfer Belts	105		
Lumbar Supports	56	Gait Belts	106-107		
Headrests	57	Transfer Slings	108		
Armrests	58-59	Bariatric			
Rain Capes & Aprons	59-60	Heel Off-loading	110		

Satisfaction Guaranteed


The best way we can show our appreciation for the trust you place in Skil-Care is to pledge that we will always offer the highest quality products at the fairest prices. All Skil-Care products composed of fabric are guaranteed to be free from defects in materials and workmanship for a period of six months following purchase. Our Smokers Apron is replacement guaranteed for one year from date of purchase. Skil-Care wheelchair cushions are sold with a one-year guarantee except those cushions specifically noted as having an 18-month guarantee. All alarm electronic units are sold with a one-year guarantee except the Econo-Alarm, Personal Alarm and Door/Window Alarm, which are guaranteed for six months. All Skil-Care mattress guarantees are stated on pages 79-81.

WOUND MANAGEMENT



Heel wounds account for 30% of pressure ulcers. Skil-Care understands the need to address this issue and offers a wide selection of heel protection devices for every budget. Implementing these items into your pressure ulcer prevention program will allow your facility to meet the day to day challenges of risk, prevention and compliance. Skil-Care offers other Skin related items such as elbow protectors, skin sleeves and pads.

HEELS-OFF



Suspends both heels off the bed for total heel pressure relief

- Aids in healing existing ulcers
- Elevates legs to relax spine
- Low-Shear II cover protects against skin-damaging friction
- Wipes clean for easy care
- Optional gel adds superior comfort and cooling

ReOrder #	Description	Size	U/M
503151	Heels-Off Small/Medium	13"L x 24"W x 5.5"H	Ea.
503150	Heels-Off Medium/Large	16"L x 24"W x 7"H	Ea.
503166	Heels-Off w/Gel Pad Small/Medium	13"L x 24"W x 5.5"H	Ea.
503165	Heels-Off w/Gel Pad Medium/Large	16"L x 24"W x 7"H	Ea.
503167	Heels-Off w/Gel Pad Bariatric	16"L x 26"W x 7"H	Ea.
503198	Heels-Off Replacement Cover	Small/Medium	Ea.
503199	Heels-Off Replacement Cover	Medium/Large	Ea.

FOOT/HEEL ELEVATOR



Sheepskin Interior (models 503101 and 503141)

Eliminates heel pressure

- Provides necessary air circulation for accelerated tissue granulation and healing of existing pressure sores
- Resilient polyfoam covered with perspiration-absorbing cozy cloth fabric or optional sheepskin inner padding.

ReOrder #	Description	Size	U/M
503100	Foot/Heel Elevator	4"W	Ea.
503101	Foot/Heel Elevator w/Sheepskin Padding	4"W	Ea.
503140	Foot/Heel Elevator	6"W	Ea.
503141	Foot/Heel Elevator w/Sheepskin Padding	6"W	Ea.

TRI-FUNCTION HEELS OFF-LOADER











The Heel Off-Loader can be used without detachable foot positioner

- Suspends heels off the bed
- Places the knee in a 5° to 10° flexion position
- Elevates legs to relax the spine
- Distributes the weight of the leg along the calf
- Two position foot drop block
- Reduces the risk of pressure ulcer development
- Detachable straps that adjust for any leg length
- Low-Shear II cover reduces skin-damaging friction and wipes clean for easy care

ReOrder #	Description	Size	U/M
503188	Tri-Function Heel Off-Loader System Short	18" Long	Ea.
503189	Heel Off-Loader for 503188 System	18"L x 20"W x 7"H	Ea.
503185	Tri-Function Heel Off-Loader System Long	22" Long	Ea.
503186	Heel Off-Loader for 503185	22"L x 20"W x 7"H	Ea.
503187	Foot Positioner for 503185 & 503188	7"L x 20"W x 12"H	Ea.

visit our website at www.skil-care.com

Pick Your Heel-Float™ Feature

Heel-Float Models	Zero heel pressure (Off-loading)	Gel insert	Low friction inner fabric	Anti-foot drop design	Adjustable strap location	Extended toe	LSII wipe-clean outer fabric	Ambulation appropriate
 Heel-Float Original (see pg. 5)	✓	✓	✓	✓	✓	✗	✗	✗
 Heel-Float II (see pg. 5)	✓	✗	✓	✓	✓	✓	✗	✗
 Heel-Float Plus (see pg. 6)	✓	✓	✓	✓	✓	✓	✗	✗
 Wipe-Clean Heel-Float (see pg. 6)	✓	✓	✓	✓	✓	✗	✓	✗
 Adjustable Heel-Float (see pg. 7)	✓	✓	✓	✓	✗	✓	✗	✗
 Heel-Float Walker Boot (see pg. 7)	✓	✓	✓	✓	✗	✓	✓	✓
 Heel-Float Adjustable Walker Boot (see pg. 7)	✓	✓	✓	✓	✗	✓	✓	✓
 Heel-Float Air Boot (see pg. 8)	✓	✗	✓	✗	✗	✗	✗	✗

Heel-Float™ ORIGINAL



Note: (Cutaway not on actual product)

NEW closure system provides choice of two ways to secure the boot:

1. One strap across the ankle and the other over the foot to prevent foot from moving back in the boot.
2. Both straps across the ankle.



Off-loads the heel over an air cavity to eliminate pressure

- Floats the heel over an air cavity for zero pressure
- Provides air circulation to aid healing of existing ulcers
- Water-based gel insert lowers pressure beneath ankle
- Covered with a soft launderable Cozy Cloth fabric
- Anti-foot-drop construction controls dorsal flexion
- Launderable and reusable (with removable gel pack)
- Prevents plantar flexion

ReOrder #	Description	Size	U/M
503034	Heel-Float - Small	3" wide	Ea.
503035	Heel-Float - Medium	4" wide	Ea.
503036	Heel-Float - Large/Bariatric	5" wide	Ea.
503095		Small	Ea.
503097	Replacement Gel Pack for Heel Float Series	Medium	Ea.
503099		Large	Ea.
503430/503098	Bilateral Anti-Rotation Gripper - Foam	Universal	Ea.



503098 Bilateral Anti-Rotation Gripper

Heel-Float™ II



Bilateral Anti-Rotation Gripper

The ultimate product in comfort, heel off-loading, and cost saving

- Off-loads the heel to eliminate pressure
- Low shear covering to reduce friction and pressure sores
- Anti-foot drop construction controls dorsal flexion
- Two versatile straps that can be multi-positioned
- Prevents plantar flexion
- Extended foot support keeps blanket and sheet off toes
- Reusable and completely launderable

ReOrder #	Description	Size	U/M
503134	Heel-Float II - Small	3" wide	Ea.
503135	Heel-Float II - Medium	4" wide	Ea.
503136	Heel-Float II - Large	5" wide	Ea.
503430/503098	Bilateral Anti-Rotation Gripper - Foam	Universal	Ea.

Heel-Float™ PLUS



Original
Skil-Care
CORPORATION
Design

Bariatric
Available



503098 Bilateral
Anti-Rotation Gripper

Floats the heel over an air cavity for zero pressure

- Extended foot support keeps blanket and sheet off toes
- Cutouts that eliminate pressure on lateral and medial malleoli
- Foot band that secures below bunions
- Water-based gel insert lowers pressure beneath ankle
- Provides maximum air circulation for healing existing sores
- Covered with soft launderable Cozy Cloth fabric
- Launderable after removing solid insert
- Prevents plantar flexion

ReOrder #	Description	Size	U/M
503044	Heel-Float Plus - Small	3" wide	Ea.
503045	Heel-Float Plus - Medium	4" wide	Ea.
503046	Heel-Float Plus - Large/Bariatric	5" wide	Ea.
503430/503098	Bilateral Anti-Rotation Gripper - Foam	Universal	Ea.

WIPE-CLEAN Heel-Float™



Original
Skil-Care
CORPORATION
Design

Bariatric
Available

Floats the heel over an air cavity for zero pressure

- Wipe-clean and reusable
- Low-Shear II cover reduces skin-damaging friction and wipes clean for easy care
- Suspends heel over air cavity to eliminate pressure
- Water-based gel insert lowers pressure beneath ankle
- Prevents plantar flexion

ReOrder #	Description	Size	U/M
503047	Small	3" wide	Ea.
503048	Medium	4" wide	Ea.
503049	Large/Bariatric	5" wide	Ea.

REPLACEMENT GEL PACKS



Replacement Gel Packs for all the Heel-Float Series, see below:

Heel-Float Original
Heel-Float Plus
Heel-Float Wipe-Clean
Heel-Float Adjustable
Heel-Float Walker Boot
Heel-Float Adjustable Walker Boot

ReOrder #	Size	U/M
503095	Small (1)	Ea.
503097	Medium (2)	Ea.
503099	Large (3)	Ea.

ADJUSTABLE Heel-Float™



Floats the heel over an air cavity for zero pressure

- Adjustable belt to accommodate plantar flexion
- Water-based gel insert lowers pressure beneath ankle
- Cozy Cloth inner lining reduces skin shear and wicks away moisture
- Launderable when you remove gel insert and plastic stabilizer
- Wipe-clean nylon outer cover

ReOrder #	Description	Size	U/M
503152	Small	3" wide	Ea.
503153	Medium	4" wide	Ea.
503154	Large/Bariatric	5" wide	Ea.

Heel-Float™



WALKER BOOT



Permits resident to stand, pivot, and transfer

- Floats the heel over an air cavity for zero pressure
- Cozy Cloth inner liner reduces friction and wicks away moisture
- Water-based gel insert lowers pressure beneath ankle
- Detachable tread pad for infection control
- Launderable when you remove gel insert and detachable tread
- Wipe-clean nylon outer cover

ReOrder #	Description	Size	U/M
503142	Small	3" wide	Ea.
503143	Medium	4" wide	Ea.
503144	Large/Bariatric	5" wide	Ea.

Heel-Float™

ADJUSTABLE WALKER BOOT



*Photo shows boot with tread being removed



Permits resident to stand, pivot, and transfer

- Floats the heel over an air cavity for zero pressure
- Water-based gel insert lowers pressure beneath ankle
- Adjustable belt to accommodate foot drop
- Detachable tread pad for infection control
- Cozy Cloth inner liner reduces friction and wicks away moisture
- Launderable when you remove gel insert and detachable tread
- Wipe-clean nylon outer cover

ReOrder #	Description	Size	U/M
503147	Small	3" wide	Ea.
503148	Medium	4" wide	Ea.
503149	Large/Bariatric	5" wide	Ea.

Heel-Float™ AIR BOOT



751580 Hand Pump
(Sold Separately)

Zero pressure on the heel

- Off-loads the heel with hi-rise air pocket
- Provides maximum air circulation for healing or prevention of heel or leg ulcers
- Comfortable inner lining for reduction of shear forces.
- Wipe-clean and reusable


ReOrder #	Description	Size	U/M
503510	Heel-Float Air Boot	Universal	Ea.
751580	Hand Pump	Universal	Ea.

BILATERAL ANTI-ROTATION STABILIZERS



NEW IMPROVED
Bilateral Anti-Rotation
Gripper

503098



503038



Bilateral Stabilizer



Plastic Stabilizer

Universal foot rotation control

- Use with most Skil-Care foot boots to prevent foot rotation
- Recommended for hip-fracture residents

ReOrder #	Description	Size	U/M
503430/503098	Bilateral Anti-Rotation Gripper - Foam	Universal	Ea.
503039	Bilateral Anti-Rotation Plastic	Small/ Medium	Ea.
503038	Bilateral Anti-Rotation Plastic	Large/ Bariatric	Ea.

FOOT DROP BOOT



Model 503171

Adjustable Extended Toe

Model 503174

Rigid plastic boot controls plantar flexion (foot drop)

- Heel is suspended above a cavity for maximum pressure relief
- Supportive plastic inner liner controls plantar flexion (foot drop)
- Synthetic sheepskin liner provides comfort, prevents bruising and minimizes shear force
- Fits either foot and can be laundered after removing solid insert
- Model 503174 is available with extended toe and anti-rotation bar

ReOrder #	Description	Heel to Toe Height	U/M
503170	Small - 3"W	9.5"H	Ea.
503171	Medium - 4"W	11"H	Ea.
503172	Large - 5"W	12.5"H	Ea.
503174	Extended Toe & Lateral Support- Universal - 4"W	11"H	Ea.
503180		Small	Ea.
503181	Replacement Cover	Medium	Ea.
503182		Large	Ea.

SUPER SOFT HEEL PROTECTOR



Heel Protector
w/Bilateral Anti-Rotation Gripper



Bilateral
Anti-Rotation Gripper



This superior poly-filled heel protector combines off-loading of the heel in addition to providing comfort and control of foot drop

- Universally sized poly-filled heel protector designed for treatment and prevention of heel pressure sores
- Open heel design effectively eliminates heel pressure by suspending the heel in the air and also promotes faster healing
- Inner lining is covered with a soft launderable Cozy Cloth fabric
- Extended height prevents sheets and blankets from contacting the toes
- Reduces risk of plantar flexion (foot drop)
- Easy to apply with adjustable straps for custom fit
- Launderable, reusable, and guaranteed for six months
- DVT compatible

ReOrder #	Description	Size	U/M
503410	Super Soft Heel Protector	Universal	Ea.
503420	Super Soft Heel Protector w/Bilateral Anti-Rotation Gripper	Universal	Ea.
503430/503098	Bilateral Anti-Rotation Gripper - Foam	Universal	Ea.


FOAM PRESSURE RELIEVING HEEL PROTECTORS



Convoluted (503450)



Smooth (503400)



Auxiliary Pad

Low Friction
Heel Insert

This economical heel protector comes in smooth and convoluted foam

- Off-loads the heel for decubitus ulcer protection
- Heel insert covered with low friction fabric for reduction of sweating and friction burns
- Two soft, comfortable cozy cloth tie straps
- Includes auxiliary foam pad for extending the length or height
- Multiple ventilation holes for cooling affect
- Reusable and completely launderable

ReOrder #	Description	Size	U/M
503400	Foam Heel Protector - Smooth	Universal	Ea.
503450	Foam Heel Protector - Convoluted	Universal	Ea.

SOFT BOOT



Provides pressure-relief under heel, ankle, and lower calf

- New longer length for upper-calf skin protection
- Made from plush synthetic sheepskin
- Large vent holes promote air circulation to reduce perspiration and moisture build-up
- Secures with Hook & Loop closures
- Completely launderable

ReOrder #	Description	Size	U/M
503160	Soft Boot	Universal	Pair

ULTRA-SOFT FIBER-FILLED HEEL CUSHION



Helps reduce the risk of pressure sores

- Filled with resilient spun fiber
- Covered with a soft launderable Cozy Cloth fabric
- Hook & Loop closure for easy application and secure positioning
- Reusable and completely launderable
- Cover wicks away moisture to eliminate skin-damaging perspiration

ReOrder #	Description	Size	U/M
503030	Ultra-Soft Heel Cushion	Universal	Pair

GEL-FOAM HEEL CUSHION



Helps reduce the risk of pressure sores

- Combines gel and foam for maximum pressure relief
- Water-based gel insert improves pressure relief and cools leg area
- Ample foam padding adds extra cushioning
- Covered with a soft Cozy Cloth fabric
- Launderable and reusable (with removal of gel pack)
- Cover wicks away moisture to eliminate skin damaging perspiration
- Extra-wide Hook & Loop strap secures heel cushion

ReOrder #	Description	Size	U/M
503025	Gel-Foam Heel Cushion	Universal	Pair
503088	Replacement Gel Pack	3.25"W x 6.5"D	Ea.

COVERED FOAM HEEL PROTECTOR



Filled with soft foam padding

Helps reduce the risk of pressure sores


- Covered with a soft launderable Cozy Cloth fabric
- Hook & Loop for easy application to secure positioning
- Reusable and completely launderable
- Cover wicks away moisture to eliminate skin-damaging perspiration

ReOrder #	Description	Size	U/M
503040	Foam Heel Cushion	Universal	Pair

ECONO FOAM HEEL PROTECTOR



Standard Model 503070



Convuluted Model 503080

Oversized design protects the heel

- Made from pressure-relieving polyfoam
- Hook & Loop closures assure a snug, comfortable fit
- Available in standard or convuluted polyfoam

ReOrder #	Description	Size	U/M
503070	Econo Foam Heel Protector	Universal	Pair
503077	Standard	Universal	12 Pair/Cs
503080	Econo Foam Heel Protector	Universal	Pair
503087	Convuluted	Universal	12 Pair/Cs

ECONO SHEEPSKIN HEEL PROTECTOR



- Light weight and economical heel protector
- Synthetic sheepskin is low shear and allows for air circulation
- Stays soft even after repeated launderings



ReOrder #	Description	Size	U/M
503005	Econo Sheepskin Heel Protector	Universal	12 Pair/Cs
503006	Heel Protector	Universal	36 Pair/Cs

MULTI-PLY HEEL PROTECTOR



Triple-Ply Technology

Single-Ply Technology



Provides three layers of sheepskin under sensitive heels

- Maintains resiliency for constant pressure relief
- Side vent holes promote air flow and reduce moisture build-up
- Hook & Loop closures assure a snug, comfortable fit
- Completely launderable

ReOrder #	Description	Size	U/M
503015	Multi-Ply Heel Protector	Universal	Pair
503016	Protector	Universal	36 Pair/Cs

STAY-ON HEEL PROTECTOR



Slide-Lock Buckle Closure



Special design makes it difficult to kick off

- Secure slide-lock buckle closure
- Two layers of plush synthetic sheepskin for superior pressure relief
- Stays soft and resilient even after repeated laundering

ReOrder #	Description	Size	U/M
503050	Stay-On Heel Protector	Universal	Pair
503060		Universal	12 Pair/Cs

HEEL/ELBOW PROTECTOR



Affordable and effectively protects sensitive heels and elbows

- Prevents pressure sores, friction burns & abrasions
- Soft polyfoam pad is specially contoured to comfortably cover resident's heels, elbows, or knees
- Seamless inside to prevent abrasion
- Oversized stockinette sleeve keeps pad securely positioned without constricting arm or ankle

ReOrder #	U/M	ReOrder#	U/M	Size
503203	Pair	503213	36 Pair/Cs	Small
503204		503214		Medium
503205		503215		Large

ELBOW PAD PROTECTOR



Visco-Foam Elbow Pad Protector
(Wipe-Clean Fabric)



Synthetic Sheepskin Elbow Protector (Launderable Fabric)



Provides protection against pressure

- Increases comfort
- Designed to flex naturally with resident's arm

ReOrder #	Description	Size	U/M
503230	Visco-Foam Elbow Pad w/LSII Cover	Universal	Pair
503240		Universal	36 Pair/Cs
503310	Synthetic Sheepskin Elbow Protector	Universal	Pair
503320		Universal	36 Pair/Cs

POUCH ARM SLING

Universal arm support for either arm

- Adjustable shoulder strap
- Covered with a soft launderable Cozy Cloth fabric
- Non-slip buckle

ReOrder #	Description	Size	U/M
902011	Pouch Arm Sling	Small-Medium	Ea.
902012		Medium-Large	Ea.



GERI-SLEEVES




LIGHT TONE



PADDED



DARK TONE



Close up of padding




- Extra-heavy stockinette maintains skin integrity
- Breathable cotton-blend Lycra® protects against skin tears
 - Thumb-hole cutouts keep arm sleeves in place

- Leg sleeves prevent skin shear at the heels and calves
- Protects I.V. and wound sites
- Reusable and launderable

ReOrder #	Description	Size	U/M	ReOrder #	U/M	Stripe Color	Length	Normal Circumference	Fits Maximum Limb Circumference
503355	Light Tone Geri-Sleeves	Small Arm	Pair	503356	25 Pr/Pk	Red	15"	9 1/2" / 24.1cm	20" / 51cm
503350		Medium Arm	Pair	503352	25 Pr/Pk	Green	20"	11 3/4" / 29.8cm	25" / 63.5cm
503357		Bariatric Arm	Pair	503358	25 Pr/Pk	Blue	20"	15 3/4" / 40.0cm	31" / 78.7cm
503360		Universal Leg	Pair	503361	25 Pr/Pk	Green	16"	9 1/2" / 24.1cm	15 3/4" / 40.0cm
503364	Dark Tone Geri-Sleeves	Small Arm	Pair	503365	25 Pr/Pk	Red	15"	9 1/2" / 24.1cm	20" / 51cm
503362		Medium Arm	Pair	503363	25 Pr/Pk	Green	20"	11 3/4" / 29.8cm	25" / 63.5cm
503366		Bariatric Arm	Pair	503367	25 Pr/Pk	Blue	20"	15 3/4" / 40.0cm	31" / 78.7cm
503368		Universal Leg	Pair	503369	25 Pr/Pk	Green	16"	9 1/2" / 24.1cm	15 3/4" / 40.0cm
503392	Padded Geri-Sleeves	Small	Pair				9 1/2" / 24.1cm		20" / 51cm
503390		Medium	Pair				11 3/4" / 29.8cm		25" / 63.5cm
503394		Large/Bariatric	Pair				15 3/4" / 40.0cm		31" / 78.7cm

Note: For padded geri-sleeve circumference see regular geri-sleeve sizing above

CUSHIONS




Comfort

753178



Premium

751046



Advanced

751545



Factors to consider when selecting a cushion are pressure ulcer risk, skin condition, positioning requirements and weight. Our wide selection of cushions are specifically designed to meet the needs of today's residents. Whether that need is Comfort, Pelvic Positioning, Forward Sliding, Coccyx Relief or Bariatric support you can find the right cushion to meet your needs. The cushion materials range from Foam, Gel over Foam, Gel within Foam, All Gel and Air.

Foam cushions provide a comfortable seating surface and basic pressure redistribution.

Gel Cushions provide premium pressure redistribution because the gel cradles to the individual therefore providing better pressure redistribution. The gel helps cool the seat surface by reducing heat buildup and provides a customized surface for optimum pressure redistribution. The gel chambers keep the gel in the area of highest pressure and provide greater support and stability.

Air Cushions provide an advanced level of pressure redistribution properties because the weight is redistributed evenly throughout the entire cushion

Review the Skil-Care Cushion selection chart to help find the ideal cushion.

WHAT IS YOUR CUSHION TYPE?


Comfort

Premium

Advanced

Level of Pressure Redistribution	Cushion Material	Cushion Type	Page	Comfort	Forward Sliding Control	Coccyx Relief (option)	Bariatric Support (option)	Cooling Gel	Thin-Line
Comfort	Foam	Pressure Check Foam Cushion	20	•					
		Conform	27	•	•	•			•
		Contour	27	•	•				
		Anti-Thrust	24	•	•				
		Comfort	19	•					
Premium	Gel-Foam	Econo-Gel	16	•				•	•
		Gel-Foam	17	•		•	•	•	•
		Super-Soft	16	•				•	
		Ultra Cushion	31	•				•	•
		Position Plus	26	•	•		•	•	•
		Anti-Thrust with Gel	24	•	•			•	
		X-Gel	18	•	•	•		•	•
		Slide Guard	25	•	•			•	
	Gel Only	X-Gel Overlay	30	•				•	
		Gel Lift	19	•				•	
Advanced	Air	Cushion Air	22	•					
		Foam Air	22	•		•			•
		Air Lift	23	•					
		Air Lock	23	•					

ECONO-GEL CUSHION



Two-part construction for greater stability



Original Skil-Care Corporation Design



751162-751164
Econo-Gel Cushion with Sheepskin Cover




Bimini Blue Vinyl (Model 751180 & 751181)

Pressure relief and redistribution at a budget price

- Two chamber gel pad design for greater stability
- Comfort foam on top and bottom
- Gel Pad between foam layers provides pressure-relieving floatation
- Two toned cover serves as a reminder to reposition resident
- Available polyester outer cover reduces skin shear
- Soft synthetic sheepskin cover enhances comfort and helps prevent decubitus ulcers (see models 751162, 751163 & 751164)

ReOrder #	Description	Size	U/M
751148	Vinyl w/Poly Cvr	18"W x 18"D x 2"H	Ea.
751156	Vinyl Only		Ea.
751149	Vinyl w/Poly Cvr	20"W x 18"D x 2"H	Ea.
751154	Vinyl Only		Ea.
751180	Bimini Blue Vinyl	16"W x 16"D x 2"H	Ea.
751150	Vinyl w/Poly Cvr		Ea.
751162	Vinyl with Sheepskin Cvr - Not Reminder		Ea.
751155	Vinyl Only		Ea.
751181	Bimini Blue Vinyl	18"W x 16"D x 2"H	Ea.
751151	Vinyl w/Poly Cvr		Ea.
751163	Vinyl with Sheepskin Cvr - Not Reminder		Ea.
751157	Vinyl Only		Ea.
751153	Vinyl w/Poly Cvr	20"W x 16"D x 2"H	Ea.
751164	Vinyl with Sheepskin Cvr - Not Reminder		Ea.

SUPER SOFT ECONO GEL CUSHION




Reminder

This cushion is designed to relieve pressure and increase comfort on any seating surface.

- Voted Best Value
- Super soft surface for maximum comfort
- Provides pressure relief and redistribution for all body types
- Can be used on either side
- Two toned cover serves as a reminder to reposition resident
- Low-Shear II stretch cover surrounding a gel-foam composition
- Water-resistant cover is wipe clean

ReOrder #	Description	Size	U/M
751166		16"W x 16"D x 2"H	Ea.
751167	Super Soft Econo-Gel Cushion w/ LSII Cover	18"W x 16"D x 2"H	Ea.
751168		20"W x 16"D x 2"H	Ea.

GEL-FOAM CUSHION (CLASSIC)



Two-layered design for maximum pressure relief and redistribution

- Water-based gel lowers interface pressure and cools to minimize perspiration
- High-density foam prevents bottoming-out
- Low-Shear I launderable outer cover that reduces skin damaging friction
- Vinyl inner cover wipes clean for easy care and prevents urine from entering cushion

ReOrder #	Description	Size	U/M	
751018	Gel-Foam Cushion Vinyl Cover Only	16"W x 16"D x 2.5"H	Ea.	
751010	Gel-Foam Vinyl Cushion w/LSI Cover		Ea.	
751021	Gel-Foam Cushion Vinyl Cover Only	18"W x 18"D x 2.5"H	Ea.	
751013	Gel-Foam Vinyl Cushion w/LSI Cover		Ea.	
751019	Gel-Foam Cushion Vinyl Cover Only	18"W x 16"D x 2.5"H	Ea.	
751011	Gel-Foam Vinyl Cushion w/LSI Cover		Ea.	
751020	Gel-Foam Cushion Vinyl Cover Only	20"W x 16"D x 2.5"H	Ea.	
751012	Gel-Foam Vinyl Cushion w/LSI Cover		Ea.	
751014	Gel-Foam Vinyl Cushion w/LSI Cover	20"W x 18"D x 2.5"H	Ea.	
ReOrder #	Description	Coccyx Cutout	Size	U/M
751025	Gel-Foam Vinyl Cushion w/ Coccyx Cutout & LSI Cover		18"W x 16"D x 2.5"H	Ea.
754528	Gel-Foam Vinyl Cushion w/ Coccyx Cutout & LSI Cover		18"W x 16"D x 3.5"H	Ea.
751024	Gel-Foam Vinyl Cushion w/ Coccyx Cutout & LSI Cover		20"W x 16"D x 2.5"H	Ea.

SUPER SOFT GEL-FOAM CUSHION


- Two-Chamber water based gel distributes weight to relieve pressure.
- High-density base prevents bottoming-out and increases sitting comfort
- 4-Way stretch cover provides an industry top choice of 4-way stretch fabric that provides less moisture, more comfort, and cooler feel than a standard fabric cover.
- Non-slip bottom keeps cushion securely in place for extra safety.
- Water-resistant cover is wipe clean



ReOrder #	Description	Size	U/M
751050		16"W x 16"D x 2.5"H	Ea.
751051	Super Soft Gel-Foam Cushion w/ 4-Way Stretch Cover	18"W x 16"D x 2.5"H	Ea.
751052		20"W x 16"D x 2.5"H	Ea.

SITTIN' PRETTY™ GEL-FOAM CUSHION

- Water-based gel lowers interface pressure and cools to minimize perspiration
- High-density comfort foam that reduces pressure
- Vinyl cover wipes clean for easy care and prevents urine from entering cushion



ReOrder #	Description	Size	U/M
751200		16"W x 16"D x 2.5"H	Ea.
751201	Bimini Blue Gel-Foam Cushion	18"W x 16"D x 2.5"H	Ea.
751202		18"W x 18"D x 2.5"H	Ea.

Sittin' Pretty™

THIN-LINE GEL-FOAM CUSHION




Maximum pressure relief and redistribution with a reduced height

- Gel flotation provides pressure relief and redistribution
- Water-based gel cools to reduce skin damaging perspiration
- High-density foam prevents bottoming out
- Incontinence-proof vinyl inner cover wipes clean
- Available in Low-Shear I launderable outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
751038	Thin-Line Gel-Foam Cushion - Vinyl Cover Only	16"W x 16"D x 1.5"H	Ea.
751039	Thin-Line Gel-Foam Cushion - Vinyl w/LSI Cover	16"W x 16"D x 1.5"H	Ea.
914481	Ultra Thin Gel-Foam Cushion - Vinyl w/LSI Cover	16"W x 16"D x 1"H	Ea.
751041	Thin-Line Gel-Foam Cushion - Vinyl Cover Only	18"W x 16"D x 1.5"H	Ea.
751040	Thin-Line Gel-Foam Cushion - Vinyl w/LSI Cover	18"W x 16"D x 1.5"H	Ea.
914482	Ultra Thin Gel-Foam Cushion - Vinyl w/LSI Cover	18"W x 16"D x 1"H	Ea.
751043	Thin-Line Gel-Foam Cushion - Vinyl Cover Only	20"W x 16"D x 1.5"H	Ea.
751042	Thin-Line Gel-Foam Cushion - Vinyl w/LSI Cover	20"W x 16"D x 1.5"H	Ea.
914483	Ultra Thin Gel-Foam Cushion - Vinyl w/LSI Cover	20"W x 16"D x 1"H	Ea.

STARRY-NIGHT GEL-FOAM CUSHION



☆☆☆
Starry Night

ReOrder #	Description	Size	U/M
751196		16"W x 16"D x 2.5"H	Ea.
751197	Starry Night Gel-Foam Cushion	18"W x 16"D x 2.5"H	Ea.
751198		18"W x 18"D x 2.5"H	Ea.

GEL-FOAM X-CUSHION



Gel-Foam X-Cushion



Gel-Foam X-Cushion w/Coccyx Cutout



X-tra positioning, stability and pressure relief

- Water-based gel relieves pressure and cools the skin to prevent breakdown caused by perspiration and heat build-up
- High-density foam adds comfort and prevents bottoming out
- Available in a Low-Shear I launderable outer cover that reduces skin damaging friction
- Incontinent-proof vinyl inner cover prevents urine from entering cushion
- Wipes clean for easy care

ReOrder #	Description	Size	U/M
751030		16"W x 16"D x 2.5"H	Ea.
751026		18"W x 16"D x 2.5"H	Ea.
751028	GFX-Cushion Vinyl Only	18"W x 18"D x 2.5"H	Ea.
751032		20"W x 16"D x 2.5"H	Ea.
751031		16"W x 16"D x 2.5"H	Ea.
751027		18"W x 16"D x 2.5"H	Ea.
751029	GFX-Cushion Vinyl w/LSI Cover	18"W x 18"D x 2.5"H	Ea.
751033		20"W x 16"D x 2.5"H	Ea.
751037		20"W x 18"D x 2.5"H	Ea.

Coccyx Cutout

751034		16"W x 16"D x 2.5"H	Ea.
751035	GFX Coccyx Cutout Vinyl w/LSI Cover	18"W x 16"D x 2.5"H	Ea.
751036		20"W x 16"D x 2.5"H	Ea.

TOILET SEAT CUSHION



Gel-Foam Toilet Seat

Foam Toilet Seat

Layered with gel and foam to provide comfort and maximum pressure relief and redistribution


- Gel-foam or Foam only cushion provides comfort and maximum pressure relief and redistribution
- Foam base on both models protects against bottoming out
- Horseshoe shape fits most commercial toilet seats. Enlarged size (14" W x 18" L x 1.5" H) assures that there's no exposed toilet seat surfaces
- Four removable straps hold the cushion in place
- Vinyl surface can be cleaned with a spray and wipe anti-bacterial cleaner

ReOrder #	Description	Size	U/M
915325	Gel-Foam Toilet Seat Cushion	14" W x 18" D x 1.5" H	Ea.
915330	Foam Only Toilet Seat Cushion		
909397	Four replacement straps/pack	Universal	24 pcs (6 pks)

GEL-LIFT CUSHION

This water-based gel surface provides pressure-relieving flotation and cooling.

- Provides added comfort
- Pressure relieving cooling gel
- Wipe-clean anti-bacterial vinyl




ReOrder #	Description	Size	U/M
751100	Gel-Lift Cushion w/ Safety Ties - Green	16"W x 16"D x 1"H	Ea.
751101	Gel-Lift Cushion w/ Safety Ties - Blue	18"W x 16"D x 1"H	Ea.
781030	Optional LSII Cover Wipe-clean	18"W x 16"D x 1"-2"H	Ea.

FOAM ONLY CUSHIONS

COMFORT FOAM WHEELCHAIR CUSHION



Reminder



Bimini Blue Cushion Models 751188 & 751189



Comfort and quality at reasonable prices

- High-density foam provides comfort and relieves pressure
- Two-color surfaces act as a reminder to reposition resident
- Extra-high E-Z Rise models support post-hip surgery care
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction

Note: Foam only products are guaranteed for 6 months.

ReOrder #	Description	ReOrder #	Description	ReOrder #	Description	Size	U/M
753232	Foam w/two color vinyl	-	-	-	-	16"W x 16"D x 2"H	Ea.
753272	Foam w/two color vinyl	-	-	-	-	16"W x 16"D x 3"H	Ea.
753230	Foam w/two color vinyl	753210	Foam w/two color vinyl and LSII Cover (as shown)	753240	Foam Only	18"W x 16"D x 2"H	Ea.
753270	Foam w/two color vinyl	753250	Foam w/two color vinyl and LSII Cover	753280	Foam Only	18"W x 16"D x 3"H	Ea.
753271	Foam w/two color vinyl	753251	Foam w/two color vinyl and LSII Cover	-	-	18"W x 18"D x 3"H	Ea.
753290	E-Z Rise w/ Vinyl Cover	753285	E-Z Rise w/LSII Cover	753295	E-Z Rise Foam Only	18"W x 16"D x 4"H	Ea.
707052	E-Z Rise w/ Vinyl Cover	-	-	-	-	18"W x 16"D x 6"H	Ea.
751188	Bimini Blue Vinyl Foam Cushion					18"W x 16"D x 2"H	Ea.
751189	Bimini Blue Vinyl Foam Cushion					18"W x 16"D x 3"H	Ea.

SOLID FOAM CUSHION W/ SHEEPSKIN COVER



Economical pressure relieving surface

- Soft foam provides pressure relief and redistribution
- Stable surface for residents who have difficulty maintaining upper trunk control
- Lightweight and can be used under all temperature conditions
- Synthetic sheepskin cover for cooling effect and low shear surface

ReOrder #	Description	Size	U/M
753220	Solid Foam Cushion with Synthetic Sheepskin Cover	16"W x 16"D x 2"H	Ea.
753260	Solid Foam Cushion with Synthetic Sheepskin Cover	18"W x 16"D x 2"H	Ea.

PRESSURE-CHECK LIGHT/HEAVY FOAM CUSHION



- **Dual-sided** checkered cushion provides foam cells that adjust to the resident's anatomy
- Side with larger pods designed for heavier residents (150lbs or more)
- Side with smaller pods designed for lighter residents (149lbs or less)
- Both sides provide maximum weight distribution for pressure relief
- Available in a Low-Shear II wipe-clean outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
753184	Pressure-Check Light/Heavy Foam Cushion (Foam Only)	18"W x 16"D x 3"H	Ea.
753186	Pressure-Check Light/Heavy Foam Cushion w/LSII Cover	18"W x 16"D x 3"H	Ea.

PRESSURE-CHECK FOAM CUSHION




Reduces pressure and increases comfort

- Foam sections adjust to meet individual needs
- Provides maximum weight distribution for pressure relief and redistribution
- Available in a Low-Shear II wipe-clean outer cover that reduces skin damaging friction
- Available in a urethane waterproof clear film cover
- Soft synthetic sheepskin cover enhances comfort and helps prevent decubitus ulcers (*for models 753182 and 753183*)

ReOrder #	Description	Size	U/M
753178	Pressure-Check Foam w/Urethane Waterproof Clear Film Cover	18"W x 16"D x 3"H	Ea.
753179	Pressure-Check Foam w/Urethane Waterproof Clear Film Cover	18"W x 16"D x 4"H	Ea.
753180	Pressure-Check Foam w/LSII Cover	18"W x 16"D x 3"H	Ea.
753181	Pressure-Check Foam w/LSII Cover	18"W x 16"D x 4"H	Ea.
753182	Pressure-Check Foam w/Sheepskin Cover	18"W x 16"D x 3"H	Ea.
753183	Pressure-Check Foam w/Sheepskin Cover	18"W x 16"D x 4"H	Ea.

Note: Urethane waterproof clear film covering is fragile and for more durable cover refer to the Low Shear II cover.


GEL-INFUSED VISCO CUSHION TOPPER



- Covered Gel-Infused Visco foam topper upgrades any sub-level 18"x16" cushion
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
914441	Gel-Infused Visco Cushion Topper	18"W x 16"D x 1"H	Ea.

WEDGED FOAM WHEELCHAIR CUSHION



- Controls sliding and provides pressure relief and redistribution
- High-density foam provides comfort and reduces pressure
 - Incontinent-proof, wipe-clean vinyl inner cover
 - Two-color surfaces act as a reminder to reposition resident
 - Available polyester cover removes for easy laundering

ReOrder #	Description	ReOrder #	Description	ReOrder #	Description	Size	U/M
754016	Soft Foundation Vinyl Only	-	-	-	-	16"W x 16"D x 6"H x 3"H	Ea.
754039	Soft Foundation Vinyl Only	-	-	-	-	16"W x 16"D x 4"H x 2"H	Ea.
754020	Soft Foundation Vinyl Only	754010	Soft Foundation Vinyl w/LSII	754021	Soft Foundation Vinyl w/Polyester	18"W x 16"D x 6"H x 3"H	Ea.
754036	Firm Foundation Vinyl Only	754037	Firm Foundation Vinyl w/LSII	754035	Firm Foundation Vinyl w/Polyester	18"W x 16"D x 6"H x 3"H	Ea.
754041	Soft Foundation Vinyl Only	754040	Soft Foundation Vinyl w/LSII	754048	Soft Foundation Vinyl w/Polyester	18"W x 16"D x 4"H x 2"H	Ea.
754051	Firm Foundation Vinyl Only	754050	Firm Foundation Vinyl w/LSII	-	-	18"W x 16"D x 4"H x 2"H	Ea.
754043	Soft Foundation Vinyl Only	-	-	-	-	20"W x 16"D x 4"H x 2"H	Ea.

Note: Firm foundation cushions have a rigid foam bottom to minimize the sling-seat effect

Wedged GEL-FOAM cushion

ReOrder #	Description	Size	U/M
754032	Gel-Foam Soft Foundation Vinyl w/Polyester	18"W x 16"D x 6"D x 3"H	Ea.
754044	Gel-Foam Soft Foundation Vinyl w/LSI	18"W x 16"D x 4"D x 2"H	Ea.
754054	Gel-Foam Firm Foundation Vinyl w/Polyester	18"W x 16"D x 4"D x 2"H	Ea.

EZ DRY FOAM CUSHION




Reticulated foam cushion that can be hand washed

- Soft foam provides pressure relief and redistribution
- Stable surface for residents who have difficulty in maintaining upper trunk control
- Light weight
- Two-color surfaces act as a reminder to reposition resident
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
753400	16"W w/LSII Cover		
753402	18"W w/LSII Cover	16"D x 2"H	Ea.
753404	20"W w/LSII Cover		
753406	16"W w/LSII Cover		
753408	18"W w/LSII Cover	16"D x 3"H	Ea.
753410	20"W w/LSII Cover		

FOAM ONLY CUSHION

CONVOLUTED
FOAM CUSHION




- Economical foam cushion for pressure relief and redistribution
- Soft foam provides pressure relief and redistribution
 - Stable surface for residents who have difficulty in maintaining upper trunk control
 - Available in two heights to meet a range of needs
 - Light weight and can be used under all temperature conditions
 - Available in a Low-Shear II wipe-clean outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
910115	Convoluted Foam Cushion	18"W x 16"D x 3"H	Ea.
910110		18"W x 16"D x 4"H	Ea.

ReOrder #	Description	Size	U/M
910117	Convoluted Foam Cushion w/LSII Cover	18"W x 16"D x 3"H	Ea.
910112		18"W x 16"D x 4"H	Ea.

CUSHION AIR



Shown with LSII Cover

- A proven design for pressure relief and positioning
- Position equalization results from the variable action of individual air cells
 - Low surface tension reduces friction and skin shear
 - Easy-to-use air pump is included with all models
 - Adjustable air pressure capability
 - All models come with a Low-Shear II wipe-clean outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
751545	Cushion Air	16"W x 16"D x 2"H	Ea.
751560		18"W x 16"D x 2"H	Ea.
751564		20"W x 16"D x 2"H	Ea.
751567		18"W x 18"D x 2"H	Ea.
751569		20"W x 18"D x 2"H	Ea.
751540	Cushion Air	16"W x 16"D x 4"H	Ea.
751550		18"W x 16"D x 4"H	Ea.
751562		20"W x 16"D x 4"H	Ea.
751566		18"W x 18"D x 4"H	Ea.
751568		20"W x 18"D x 4"H	Ea.
751580	Hand Pump		

FOAM AIR CUSHION



751580 Hand Pump Included

- Air chamber seating surface is easily adjusted to user's needs
- Air-quilt surface inflates or deflates to meet resident's pressure reduction and comfort requirements
 - Pre-inflated
 - High-density foam base provides stability and prevents "bottoming-out"
 - Incontinent-proof vinyl inner cover wipes clean for easy care
 - Low-Shear II wipe-clean outer cover that reduces skin damaging friction
 - Comes with safety straps
 - 1.5yr guarantee


ReOrder #	Description	Size	U/M
914251	Foam Air Cushion	16"W x 16"D x 2.5"H	Ea.
914250		18"W x 16"D x 2.5"H	Ea.
914252		20"W x 16"D x 2.5"H	Ea.

Coccyx Cutout

914247	Foam Air Cushion w/Coccyx Cutout	16"W x 16"D x 2.5"H	Ea.
914248		18"W x 16"D x 2.5"H	Ea.
914249		20"W x 16"D x 2.5"H	Ea.

751580	Hand Pump		
--------	-----------	--	--

AIR LIFT CUSHION



751580 Hand Pump
(Sold Separately)



Pre-inflated economical cushion

- Provides decubitus prevention and comfort
- Pressure relieving air flotation surface
- Easy to inflate or deflate (shipped pre-inflated)
- Wipe-clean anti-bacterial heat-sealed vinyl
- Available in a Low-Shear II wipe-clean outer cover that reduces skin damaging friction
- Comes with safety straps
- Meets flammability requirements of CAL 117

ReOrder #	Description	Size	Cushion Color	U/M
756210		17"W x 17"D	Blue	Ea.
756211		17"W x 17"D	Blue	12/Cs
756225	Air Lift Seat Cushion Only	19"W x 19"D	Grey	Ea.
756226		19"W x 19"D	Grey	12/Cs
756235		22"W x 28"D	Green	Ea.
756236		22"W x 28"D	Green	6/Cs
751580	Hand Pump			
756215		17"W x 17"D	Blue	Ea.
756230	Air Lift Seat Cushion w/LSII Cover	19"W x 19"D	Grey	Ea.
756237		22"W x 28"D	Green	Ea.

Cushion Covers

756219		17"W x 17"D	Ea.
756234	Available Two-Color Cushion Cover w/Straps	19"W x 19"D	Ea.
756239		22"W x 28"D	Ea.

AIR LOCK CUSHION



This soft pre-inflated anti-bacterial vinyl cushion makes any wheelchair seat more comfortable and also provides pressure relief and redistribution.

- Decubitus prevention
- Air flotation (pre-inflated)
- Pressure relief and redistribution
- Vinyl, heat-sealed construction
- Anti-bacterial wipe clean vinyl
- Safety straps
- No valves to interfere with pressure relief
- Meets flammability requirements of CAL 117

ReOrder #	Description	Size	U/M
756240	Air Lock Seat Cushion	19"W x 19"D	Ea.
756241	Air Lock Seat Cushion	19"W x 19"D	12/Cs
756245	Air Lock Seat Cushion w/LSII Cover	19"W x 19"D	Ea.
ReOrder #	Description	Size	U/M
756250	Air Lock Seat Cushion	22"W x 28"D	Ea.
756251	Air Lock Seat Cushion	22"W x 28"D	6/Cs
756255	Air Lock Seat Cushion w/LSII Cover	22"W x 28"D	Ea.

ANTI-THRUST CUSHION



Original
Skil-Care
CORPORATION
Design



Control sliding with this comfortable alternative to a pommel cushion

- Anti-thrust barrier eliminates forward and sacral sitting problems
- Helps maintain pelvic neutrality
- High-density foam increases comfort and relieves pressure
- Contoured base eliminates sling-seat hammocking and minimizes seat rise
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction
- Optional water-based gel pad provides pressure relief and redistribution in the coccyx area

ReOrder #	Description	Gel-Foam	Size	U/M
757108	Gel-Foam - Soft Base		16"W x 16"D	Ea.
757109	Gel-Foam - Firm Base*		16"W x 16"D	Ea.
757115	Gel-Foam - Soft Base (as shown)		18"W x 16"D	Ea.
757125	Gel-Foam - Firm Base*		18"W x 16"D	Ea.
757117	Gel-Foam Soft Base		20"W x 16"D	Ea.
757127	Gel-Foam Firm Base*		20"W x 16"D	Ea.

Foam Only

757110	Foam Only - Soft Base		18"W x 16"D	Ea.
757120	Foam Only - Firm Base*		18"W x 16"D	Ea.
757112	Foam Only - Soft Base		20"W x 16"D	Ea.
757122	Foam Only - Firm Base*		20"W x 16"D	Ea.
757130	Gel-Foam Cushion for Geri-Chair		20"W x 18"D	Ea.

***Note:** The firm base includes a firm, closed cell foam that helps eliminated sling-seat (hammocking)

E-Z TRANSFER SLIDER POMMEL SYSTEM



Original
Skil-Care
CORPORATION
Design

Removable pommel controls sliding and facilitates transfers

- Comfortable foam-padded pommel maintains pelvic neutrality
- Keeps legs in abduction
- Attached Hook & Loop pommel (6"H) removes for easy transfer
- May be used with or without cushion
- Effective in any chair
- Low-Shear II launderable outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
706362	E-Z Transfer Slider Pommel System	18"W x 16"D x 1"H	Ea.
706364	E-Z Transfer Slider Pommel System	20"W x 16"D x 1"H	Ea.
706363	E-Z Transfer Slider Pommel System for Geri-Chair	20"W x 18"D x 1"H	Ea.

SLIDE GUARD LEVEL & WEDGED CUSHION



For slide control and improved sitting posture

- High-density foam increases comfort and lowers pressure
- Available water-based gel lowers interface pressure and reduces perspiration
- 4" foam pommel controls sliding by eliminating posterior pelvic tilt
- Available in a Low-Shear I launderable outer cover that reduces skin damaging friction
- Incontinent-proof, vinyl inner cover prevents urine from entering cushion

Level Cushions w/LSI Cover - Gel-Foam

ReOrder #	Description	Size	U/M
751210	Level - Gel-Foam, Firm Foundation w/Vinyl and LSI Cover	18"W x 16"D x 3"H	Ea.
751230	Level - Gel-Foam, Soft Foundation w/Vinyl and LSI Cover	18"W x 16"D x 2.5"H	Ea.

Wedged Cushions w/LSI Cover - Gel-Foam

751316		16"W x 16"D x 4"- 2"H	Ea.
751320	Wedge - Gel-Foam, Soft Foundation w/Vinyl and LSI Cover	18"W x 16"D x 4"- 2"H	Ea.
751324		20"W x 16"D x 4"- 2"H	Ea.

Level Cushions w/Vinyl Cover - Foam Only

751265		18"W x 16"D x 2"H	Ea.
751269	Level - Foam Only, Soft Foundation w/Vinyl Cover	20"W x 16"D x 2"H	Ea.
751260	Level - Foam Only, Soft Foundation w/Vinyl and LSI Cover	18"W x 16"D x 2"H	Ea.

Wedged Cushions w/LSI Cover - Foam Only

751370	Wedge - Foam Only, Soft Foundation w/Vinyl and LSI Cover	18"W x 16"D x 4"- 2"H	Ea.
--------	--	-----------------------	-----

LEG ABDUCTOR CUSHION



Comfortably keeps legs in abduction and prevents hip rotation

- **Firm, Contoured Foundation** eliminates sling-seat hammocking
- 5"H contoured foam pommel provides abduction without compromising comfort
- High-density foam guarantees comfort and stability
- Fluid-proof, wipe-clean vinyl inner cover protects against incontinence
- Available in a Low-Shear I launderable outer cover that reduces skin damaging friction
- Wedged base available for additional slide control


Level Cushions

ReOrder #	Description	Size	U/M
758010	Level - Foam Cushion w/Vinyl Cover	18"W x 16"D x 3.5"H	Ea.
758015	Level - Gel-Foam Cushion w/Vinyl and Low Shear I Cover	18"W x 16"D x 4"H	Ea.

Wedged Cushions

758020	Wedged Foam w/Vinyl Cover	18"W x 16"D x 4"- 1.5"H	Ea.
758025	Wedged Gel-Foam w/Vinyl and Low Shear I Cover	18"W x 16"D x 4.5"- 2"H	Ea.

POSITION PLUS - LEVEL & WEDGED CUSHION



Four separate gel-filled chambers provide stability, pressure relief and redistribution

- Water-based gel relieves pressure and cools the skin to prevent breakdown caused by perspiration and heat build-up
- Side gel pads improve lateral stability by controlling left-right leaning
- Rear gel pad provides pressure relief in the coccyx area
- Front gel pads controls forward sliding by eliminating posterior pelvic rotation
- High-density foam adds comfort and prevents bottoming out
- Incontinent-proof vinyl inner cover prevents urine from entering cushion
- Low-Shear I launderable outer cover reduces shear forces
- Inner vinyl cover wipes clean for easy care
- Firm Foundation
- Ideal for residents who have compromised upper trunk stability

Level Cushions

ReOrder #	Description	Size	U/M
751420		16"W x 16"D x 2.5"H	Ea.
751435	Position Plus Contour Cushion - Level	18"W x 16"D x 2.5"H	Ea.
751450		20"W x 16"D x 2.5"H	Ea.

Wedged Cushions

751465		16"W x 16"D x 3.5"- 1"H	Ea.
751480	Position Plus Contour Cushion - Wedged	18"W x 16"D x 3.5"- 1"H	Ea.
751495		20"W x 16"D x 3.5"- 1"H	Ea.

STABILITY PLUS GEL-FOAM CUSHION



This is the ultimate cushion that combines stability, comfort, pressure relief and positioning

- Four dynamic chambers to assure stability
- Raised pommel to control sliding
- Water-based gel to provide pressure relief, redistribution, and therapeutic cooling
- High-density foam base to prevent bottoming out
- Incontinent-proof vinyl inner cover prevents urine from entering cushion
- Low-Shear I launderable outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
751046		16"W x 16"D x 2.5"H	Ea.
751047	Stability Plus Cushion w/LSI Cover	18"W x 16"D x 2.5"H	Ea.
751048		20"W x 16"D x 2.5"H	Ea.

FOAM CONTOUR CUSHIONS

ANTI-THRUST VISCO FOAM CUSHION



Gel-Infused Visco Foam



Visco foam in a wedged design

- Raised sides for lateral stabilization
- Raised pommel and wedge design for improved slide control
- Low-Shear II wipe-clean outer cover reduces skin damaging friction
- Rear strap for cushion stabilization

ReOrder #	Description	Size	U/M
753160	Anti-Thrust Visco Foam Cushion	18"W x 16"D x 2"H	Ea.

VISCO CONFORM CUSHION



Temperature-sensitive foam provides uniform pressure relief

- Visco foam molds to user's anatomy
- Contour leg troughs promote thigh alignment for postural symmetry
- Low-rise barrier prevents sliding by controlling forward thrust
- Wedged and coccyx cutout models are available
- Low-Shear II wipe-clean outer cover reduces skin damaging friction

Level Cushions


ReOrder #	Description	Size	U/M
753150	ConForm Cushion	18"W x 16"D x 2"H*	Ea.
753155	ConForm Cushion w/Coccyx Cutout	18"W x 16"D x 2"H*	Ea.

Wedged Cushions

753157	ConForm Wedge Cushion	18"W x 16"D x 4"- 2"H*	Ea.
--------	-----------------------	------------------------	-----

*Front Low Point

CONTOUR CUSHION



Superior positioning, pressure relief and redistribution

- Raised front prevents forward thrust
- Maintains pelvic neutrality for slide control
- Leg troughs promote thigh alignment and postural symmetry
- Elevated sides improve lateral stability
- Visco foam or gel pad available in coccyx area
- High-density foam increases comfort
- Low-Shear II wipe-clean cover reduces skin damaging friction

Visco Pads

ReOrder #	Description	Size	U/M
751610		16"W x 16"D x 4"H	Ea.
751615	Visco Pad	18"W x 16"D x 4"H	Ea.
751620		20"W x 16"D x 4"H	Ea.

Gel Pads

751625		16"W x 16"D x 4"H	Ea.
751630		18"W x 16"D x 4"H	Ea.
751635	Gel Pad	20"W x 16"D x 4"H	Ea.
751636		22"W x 18"D x 4"H	Ea.
751637		24"W x 18"D x 4"H	Ea.

CONTOUR PLUS CUSHION



Contour design plus anti-hammocking base

- Low-rise barrier controls forward thrust to prevent sliding
- Lateral support sides control left-right leaning
- Leg troughs promote proper thigh alignment and postural symmetry
- Convex base eliminates sling-seat hammocking
- Water-based gel pad relieves coccyx pressure and cools
- High-density foam provides pressure relief, comfort, and stability
- Available wedged base offers additional slide control
- Low-Shear II wipe-clean cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
751675	Level - Contour Plus Cushion w/Gel Pad	18"W x 16"D x 3"H	Ea.
751680	Wedged - Contour Plus Cushion w/Gel Pad	18"W x 16"D x 5"D x 2"H	Ea.

THIN-LINE GEL-FOAM CUSHION



- Maximum pressure relief and redistribution with a reduced height
- Gel flotation provides pressure relief and redistribution
 - Water-based gel cools to reduce skin damaging perspiration
 - High-density foam prevents bottoming out
 - Incontinence-proof vinyl inner cover wipes clean
 - Available in a Low-Shear I launderable outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
751038	Thin-Line Gel-Foam Cushion - Vinyl Cover Only	16"W x 16"D x 1.5"H	Ea.
751039	Thin-Line Gel-Foam Cushion - Vinyl w/LSI Cover	16"W x 16"D x 1.5"H	Ea.
914481	Ultra Thin Gel-Foam Cushion - Vinyl w/LSI Cover	16"W x 16"D x 1"H	Ea.
751041	Thin-Line Gel-Foam Cushion - Vinyl Cover Only	18"W x 16"D x 1.5"H	Ea.
751045	Thin-Line Gel-Foam Cushion - Vinyl Cover Only	18"W x 18"D x 1.5"H	Ea.
751040	Thin-Line Gel-Foam Cushion - Vinyl w/LSI Cover	18"W x 16"D x 1.5"H	Ea.
751044	Thin-Line Gel-Foam Cushion - Vinyl w/LSI Cover	18"W x 18"D x 1.5"H	Ea.
914482	Ultra Thin Gel-Foam Cushion - Vinyl w/LSI Cover	18"W x 16"D x 1"H	Ea.
751043	Thin-Line Gel-Foam Cushion - Vinyl Cover Only	20"W x 16"D x 1.5"H	Ea.
751042	Thin-Line Gel-Foam Cushion - Vinyl w/LSI Cover	20"W x 16"D x 1.5"H	Ea.
751049	Thin-Line Gel-Foam Cushion - Vinyl w/LSI Cover	20"W x 18"D x 1.5"H	Ea.
914483	Ultra Thin Gel-Foam Cushion - Vinyl w/LSI Cover	20"W x 16"D x 1"H	Ea.


THIN-LINE GEL-FOAM X-CUSHION



- Improved positioning and stability while providing pressure relief
- Unique X-design controls lateral leaning and forward sliding
 - Water-based gel surface provides pressure-relieving flotation
 - Cooling action reduces skin-damaging perspiration and heat build-up
 - High-density foam base prevents bottoming-out
 - Incontinent-proof vinyl inner cover wipes clean for easy care
 - Available in a Low-Shear I launderable outer cover reduces skin damaging friction
 - Available with coccyx cutout models

ReOrder #	Description	Size	U/M
751030	Thin-Line GFX-Cushion - Vinyl Only	16"W x 16"D x 2.5"H	Ea.
751026	Thin-Line GFX-Cushion - Vinyl Only	18"W x 16"D x 2.5"H	Ea.
751032	Thin-Line GFX-Cushion - Vinyl Only	20"W x 16"D x 2.5"H	Ea.
751031	Thin-Line GFX-Cushion	16"W x 16"D x 2.5"H	Ea.
751027	Thin-Line GFX-Cushion Vinyl w/LSI Cover	18"W x 16"D x 2.5"H	Ea.
751033	Thin-Line GFX-Cushion Vinyl w/LSI Cover	20"W x 16"D x 2.5"H	Ea.
751034	Thin-Line GFX-Cushion w/Coccyx Cutout	16"W x 16"D x 2.5"H	Ea.
751035	Thin-Line GFX-Cushion w/Coccyx Cutout Vinyl w/LSI Cover	18"W x 16"D x 2.5"H	Ea.
751036	Thin-Line GFX-Cushion w/Coccyx Cutout Vinyl w/LSI Cover	20"W x 16"D x 2.5"H	Ea.


THIN-LINE GEL-INFUSED VISCO CUSHION TOPPER



- Covered Gel-Infused Visco foam topper upgrades any sub-level 18"x16" cushion
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction
- Two-tone reminder cover

ReOrder #	Description	Size	U/M
914441	Gel-Infused Visco Cushion Topper	18"W x 16"D x 1"H	Ea.

THIN-LINE™ ECONO-GEL CUSHION



Two-part construction for greater stability




Low-rise pressure relief and redistribution at a budget price

- Top and bottom layers of foam provide comfort
- Two part construction relieves pressure and increases stability
- Incontinence-proof vinyl inner cover wipes clean for easy care
- Available launderable polyester outer cover reduces skin shear
- Two-tone covers serve as a reminder to reposition resident

ReOrder #	Description	Size	U/M
751172	Vinyl Only	16"W x 16"D x 1.5"H	Ea.
751170	Vinyl w/Polyester	16"W x 16"D x 1.5"H	Ea.
751173	Vinyl Only	18"W x 16"D x 1.5"H	Ea.
751176	Vinyl Only	18"W x 18"D x 1.5"H	Ea.
751171	Vinyl w/Polyester	18"W x 16"D x 1.5"H	Ea.
751178	Vinyl w/Polyester	18"W x 18"D x 1.5"H	Ea.
751175	Vinyl Only	20"W x 16"D x 1.5"H	Ea.
751177	Vinyl Only	20"W x 18"D x 1.5"H	Ea.
751174	Vinyl w/Polyester	20"W x 16"D x 1.5"H	Ea.
751179	Vinyl w/Polyester	20"W x 18"D x 1.5"H	Ea.

THIN-LINE™ ULTRA CUSHION



Firm Foundation



The ultimate low-rise firm base pressure-relief cushion

- **Firm Foundation** eliminates sling-seat hammocking
- Water-based gel relieves interface pressure
- High-density foam adds comfort and stability
- Incontinence-proof vinyl inner cover wipes clean for easy care
- Low-Shear I launderable outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
754339	Thin-Line Ultra	16"W x 16"D x 2.5"H	Ea.
754340	Cushion	18"W x 16"D x 2.5"H	Ea.

THIN-LINE™ POSITION PLUS CUSHION



Firm Foundation


Four low-profile gel pads provide lateral and forward stability

- Water-based gel relieves pressure and cools skin for effective skin protection
- Front gel pad controls sliding; side pods reduce leaning
- Rear gel pad lowers coccyx pressure
- Low-Shear I launderable outer cover reduces skin damaging friction
- Incontinence-proof vinyl inner cover wipes clean for easy care
- Firm Foundation

ReOrder #	Description	Size	U/M
751420		16"W x 16"D x 2.5"H	Ea.
751435	Position Plus - Level	18"W x 16"D x 2.5"H	Ea.
751450		20"W x 16"D x 2.5"H	Ea.

***Note: The Thin-Line Cushions are helpful for residents who leg propel the wheelchair and require a lower seat height.**

THIN-LINE™ FOAM AIR CUSHION



Original Skil-Care CORPORATION Design

Hand Pump Included

Air chamber seating surface is easily adjusted to user's needs

- Air-quilt surface inflates or deflates to meet resident's pressure relief and comfort requirements
- Pre-inflated
- High-density foam base provides stability and prevents bottoming-out
- Incontinent-proof vinyl inner cover wipes clean for easy care
- Low-Shear II wipe-clean outer cover reduces skin damaging friction
- Comes with safety straps
- 1.5yr guarantee
- Available in coccyx cutout models

ReOrder #	Description	Size	U/M
914251		16"W x 16"D x 2.5"H	Ea.
914250	Foam Air Cushion	18"W x 16"D x 2.5"H	Ea.
914252		20"W x 16"D x 2.5"H	Ea.

Coccyx Cutout

914247		16"W x 16"D x 2.5"H	Ea.
914248	Foam Air Cushion w/Coccyx Cutout	18"W x 16"D x 2.5"H	Ea.
914249		20"W x 16"D x 2.5"H	Ea.

THIN-LINE™ X-GEL OVERLAY



Reminder

SIX MONTH GUARANTEE

Original Skil-Care CORPORATION Design

Adds gel flotation to any foam cushion

- Gel topper enhances foam cushions for more effective pressure relief and redistribution
- May be used over any similar sized cushion
- Water-based gel cools to prevent skin-damaging perspiration
- Optional Low-Shear II wipe-clean outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
914429	X-Gel Overlay	16"W x 16"D x 1"H	Ea.
914430	Vinyl Only	18"W x 16"D x 1"H	Ea.
781066	Low Shear II Cover for 914429, 16" x 16"	16"W x 16"D	Ea.
781065	Low Shear II Cover for 914430, 18" x 16"	18"W x 16"D	Ea.

THIN-LINE™ VISCO CONFORM CUSHION



Visco Foam

ONE YEAR GUARANTEE


Low-rise, temperature-sensitive foam cushion

- Visco foam molds to user's anatomy
- Maximizes pressure relief and redistribution
- Leg troughs promote thigh alignment and postural symmetry
- Raised front controls sliding
- Safety strap keeps cushion in place
- Low-Shear II wipe-clean outer cover reduces skin damaging friction
- Available in coccyx cutout models

ReOrder #	Description	Size	U/M
753150	Thin-Line Visco ConForm Cushion	18"W x 16"D x 2"H*	Ea.
753155	Thin-Line Visco ConForm Cushion w/ Coccyx Cutout	18"W x 16"D x 2"H*	Ea.

*Front Low Point

ULTRA-CUSHIONS



Gel-Foam Wedge Cushion



Foam Wedge Cushion



The ultimate in pressure relief, redistribution, comfort, and support

- **Firm foundation** eliminates pressure caused by sling-seats
- High-density foam increases comfort and instability
- Contoured base rests below the seat frame so user sits lower in wheelchair
- For slide control, choose the Ultra-Wedge (below) available with all the same features.

- Available Low-Shear I launderable outer cover reduces skin damaging friction
- Optional water-based gel relieves interface pressure and cools to prevent skin damage
- Incontinent-proof, inner vinyl cover prevents urine from entering cushion

LEVEL

ReOrder #	Description	Size	U/M	ReOrder #	Description	Size	U/M
754312	Foam w/Vinyl Cover	16"W x 16"D x 3.5"H	Ea.	754332	Gel-Foam w/Vinyl & LSI Cover	16"W x 16"D x 4"H	Ea.
754310		18"W x 16"D x 3.5"H	Ea.	754330		18"W x 16"D x 4"H	Ea.
754322	Foam w/Vinyl & LSI Cover	16"W x 16"D x 3.5"H	Ea.	754334	Thin-Line Gel-Foam w/Vinyl & LSI Cover	20"W x 16"D x 4"H	Ea.
754320		18"W x 16"D x 3.5"H	Ea.	754339		16"W x 16"D x 2.5"H	Ea.
754324		20"W x 16"D x 3.5"H	Ea.	754340		18"W x 16"D x 2.5"H	Ea.

WEDGED

ReOrder #	Description	Size	U/M	ReOrder #	Description	Size	U/M
754110	Foam Ultra-Wedge w/Vinyl Cover	18"W x 16"D x 3.5"-1.5"H	Ea.	754114	Foam Ultra-Wedge w/Vinyl Cover	20"W x 16"D x 3.5"-1.5"H	Ea.
754132		16"W x 16"D x 4"-1.5"H	Ea.	754124		20"W x 16"D x 3.5"-1.5"H	Ea.
754130	Gel-Foam Ultra-Wedge w/Vinyl and LSI Cover	18"W x 16"D x 4"-1.5"H	Ea.	754174	Foam Ultra-Wedge w/Vinyl and LSI Cover	20"W x 16"D x 4.5"-1.5"H	Ea.
754180		18"W x 16"D x 5"-1.5"H	Ea.	754164		20"W x 16"D x 4.5"-1.5"H	Ea.
754160	Foam Ultra-Wedge w/Vinyl Cover	18"W x 16"D x 4.5"-1.5"H	Ea.	754134	Gel-Foam Ultra-Wedge w/Vinyl and LSI Cover	20"W x 16"D x 4"-1.5"H	Ea.
754170		18"W x 16"D x 4.5"-1.5"H	Ea.	754184		20"W x 16"D x 5"-1.5"H	Ea.
754120	Foam Ultra-Wedge w/Vinyl and LSI Cover	18"W x 16"D x 3.5"D-1.5"H	Ea.				

SLIDE GUARD LEVEL & WEDGED CUSHIONS



For slide control and improved sitting posture

- High-density foam increases comfort and lowers pressure
- **Firm foundation** eliminates pressure caused by sling-seats
- Available water-based gel lowers interface pressure
- 4" foamommel controls sliding by eliminating posterior pelvic tilt
- Low-Shear I launderable outer cover reduces skin damaging friction
- Incontinent-proof, inner vinyl cover prevents liquids from entering

Level Cushion

ReOrder #	Description	Size	U/M
751210	Level, Gel-Foam, Firm Foundation w/Vinyl and LSI Cover	18"W x 16"D x 3"H	Ea.
751280	Level, Foam Only, Firm Foundation w/Vinyl and LSI Cover	18"W x 16"D x 3"H	Ea.

Wedged Cushion

ReOrder #	Description	Size	U/M
751314	Wedged - Gel-Foam, Firm Foundation w/Vinyl and LSI Cover	16"W x 16"D x 4"-1.5"H	Ea.
751310	Wedged - Gel-Foam, Firm Foundation w/Vinyl and LSI Cover	18"W x 16"D x 4"- 2"H	Ea.
751350	Wedged - Foam Only, Firm Foundation w/Vinyl and LSI Cover	18"W x 16"D x 4"- 2"H	Ea.

LEG ABDUCTOR CUSHION



Comfortably keeps legs in abduction and prevents hip rotation

- **Firm, Contoured Foundation** eliminates sling-seat hammocking
- 5"H contoured foamommel provides abduction without compromising comfort
- High-density foam guarantees comfort and stability
- Optional water-based gel relieves pressure and cools preventing skin-damage
- Fluid-proof, wipe-clean vinyl inner cover protects against incontinence
- Available in a Low-Shear I launderable outer cover that reduces skin damaging friction
- Wedge base available for additional slide control

LEVELED

ReOrder #	Description	Size	U/M
758010	Level - Foam Cushion w/Vinyl Cover	18"W x 16"D x 3.5"H	Ea.
758015	Level - Gel-Foam Cushion w/Vinyl and Low Shear I Cover	18"W x 16"D x 4"H	Ea.

WEDGED

ReOrder #	Description	Size	U/M
758020	Wedged Foam w/ Vinyl Cover	18"W x 16"D x 4"- 1.5"H	Ea.
758025	Wedged Gel-Foam w/Vinyl and Low Shear I Cover	18"W x 16"D x 4.5"- 2"H	Ea.

FOAM WHEELCHAIR WEDGE CUSHION



Controls sliding and provides pressure relief and redistribution

- High-density foam provides comfort and reduces pressure
- **Firm foundation** eliminates pressure caused by sling-seats
- Incontinent-proof, wipe-clean vinyl cover
- Polyester cover removes for easy laundering

ReOrder #	Description	Size	U/M
754050	Foam Firm Foundation w/Vinyl & Polyester Outer Cover	18"W x 16"D x 4"- 2"H	Ea.
754051	Foam Firm Foundation w/Vinyl Cover	18"W x 16"D x 4"- 2"H	Ea.

TRI-FOAM/ GEL-INFUSED VISCO CUSHION



Firm Foundation



Gel-Infused Visco Foam

- Top layer responds to body temperature
- Temperature-sensitive Gel-Infused Visco upper surface molds to user for optimum pressure relief and redistribution
 - Middle layer of high-density foam adds comfort and prevents bottoming-out
 - **Firm foundation** eliminates sling-seat hammocking and distributes weight
 - Low-Shear II wipe-clean outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
754750		16"W x 16"D x 3"H	Ea.
754755		18"W x 16"D x 3"H	Ea.
754757	Tri-Foam Gel-Infused Visco Cushion	18"W x 18"D x 3"H	Ea.
754810		20"W x 16"D x 3"H	Ea.
754812		20"W x 18"D x 3"H	Ea.

LATERAL POSITIONING CUSHION



Gel-Infused Visco Foam

Firm Foundation


Controls Lateral leaning

- Gel-Infused Visco foam surface molds to user for reduction of high pressure points
- High-density foam core adds comfort and prevents bottoming-out
- **Firm foundation** eliminates sling-seat hammocking
- Elevated side corrects left or right leaning
- Low-Shear II wipe-clean outer cover reduces skin damaging friction
- Can be oriented to right or left side

Note: The cover must be oriented with the straps to the back of the cushion

ReOrder #	Description	Size	U/M
756010	Lateral Positioning Cushion	18"W x 16"D x 3.5"H x 5"H	Ea.

WHEELCHAIR WEDGED BASE



Firm Foundation

Slopes any flat cushion

- Elevated rigid base for eliminating hammocked surface
- Elevates cushions to help residents from sliding forward
- Maintains 90-degree hip alignment when in inclined back wheelchairs
- Wipe-clean strong ABS plastic
- Wheelchair cushion required, sold separately

ReOrder #	Description	Size	U/M
914471		18"W x 16"D x 2"H	Ea.
914472	Wheelchair Wedge Base	18"W x 16"D x 3"H	Ea.
914473		18"W x 16"D x 4"H	Ea.

J-HOOK DROP SEAT



*Shown without a Cushion



Firm Foundation

Adjustable drop-seat can be lowered up to four inches

- Base can be set to level or inclined position to meet user's positioning needs
- Safety brackets keep the back of the seat down when resident leans forward
- Conveniently lifts off when chair must be folded
- Drop seat weight limit:
 - 16"-18" Wheelchairs 200lbs. max
 - 20" Wheelchairs 350lbs. max

ReOrder #	Description	Size	U/M
704311		Fits 16" Wheelchairs	Ea.
704310	Drop Seat Base Only	Fits 18" Wheelchairs	Ea.
704312		Fits 20" Wheelchairs	Ea.
704315		Fits 16" Wheelchairs	Set
704320	Drop Seat w/Gel Cushion	Fits 18" Wheelchairs	Set
704321		Fits 20" Wheelchairs	Set

Replacement Cushions

704325		Fits 16" Wheelchairs	Ea.
704330	Gel-Foam Cushion	Fits 18" Wheelchairs	Ea.
704331		Fits 20" Wheelchairs	Ea.
704348	Replacement J-Hooks Hardware	Universal	Sets

SOLID SEAT PLATFORM



Eliminates wheelchair sling-seat hammocking

- Vinyl-covered hardboard platform levels sagging wheelchair seat
- Provides a firm foundation for any wheelchair cushion
- Fits wheelchairs from 18 to 20 inches wide
- Vinyl cover is incontinent-proof
- Wipes clean for easy care
- Comes with safety straps
- Wheelchair cushion required, sold separately

ReOrder #	Description	Size	U/M
703082	Solid Seat Platform	18"W x 16"D x .5"H	Ea.

WHEELCHAIR LEVELING PAD



757015



Firm Foundation

Eliminates Sling-Seat Hammocking



The perfect foundation for all therapeutic cushions


- Makes wheelchair sling seats flat and level with chair frame
- Eliminates sling seat hammocking
- Solid, closed cell foam insert
- Available with or without vinyl cover
- Wheelchair cushion required, sold separately

ReOrder #	Description	Size	U/M
757010	Wheelchair Leveling Pad No Cover	18"W x 16"D	Ea.
757015	Wheelchair Leveling Pad Vinyl Cover	18"W x 16"D	Ea.
757017	Wheelchair Leveling Pad Vinyl Cover	20"W x 16"D	Ea.

GEL-FOAM CUSHION (CLASSIC)



Original Skil-Care Corporation Design



Two-layered design for maximum pressure relief and redistribution

- Water-based gel lowers interface pressure and cools to prevent pressure sores
- High-density foam prevents bottoming-out
- Vinyl inner cover wipes clean for easy care and keeps fluids out

ReOrder #	Description	Size	U/M
751025		18"W x 16"D x 2.5"H	Ea.
754528	Gel-Foam Vinyl Cushion w/ Coccyx Cutout & LSI Cover	18"W x 16"D x 3.5"H	Ea.
751024		20"W x 16"D x 2.5"H	Ea.

GEL-FOAM X-CUSHION

Original Skil-Care Corporation Design




X-tra positioning, stability and pressure relief and redistribution

- Water-based gel relieves pressure and cools the skin to prevent breakdown
- High-density foam adds comfort and prevents bottoming out
- Low-Shear I launderable outer cover reduces skin damaging friction
- Incontinent-proof vinyl inner cover prevents fluids from entering cushion

ReOrder #	Description	Size	U/M
751034		16"W x 16"D x 2.5"H	Ea.
751035	GFX w/Coccyx Cutout w/LSI Cover	18"W x 16"D x 2.5"H	Ea.
751036		20"W x 16"D x 2.5"H	Ea.

FOAM AIR CUSHION



Air chamber seating surface is easily adjusted to user's needs


- Air-quilt surface inflates or deflates to meet resident's pressure reduction and comfort requirements (Pre-inflated)
- High-density foam base provides stability and prevents "bottoming-out"
- Incontinent-proof vinyl inner cover wipes clean for easy care
- Low-Shear II wipe-clean outer cover reduces skin damaging friction
- 1.5yr guarantee

ReOrder #	Description	Size	U/M
914247		16"W x 16"D x 2.5"H	Ea.
914248	Foam Air Cushion w/Coccyx Cutout	18"W x 16"D x 2.5"H	Ea.
914249		20"W x 16"D x 2.5"H	Ea.

VISCO CONFORM CUSHION

Temperature-sensitive foam provides uniform pressure relief

- Visco foam molds to user's anatomy
- Contour leg troughs promote thigh alignment for postural symmetry
- Low-rise barrier prevents sliding by controlling forward thrust
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction



Visco Foam



Original Skil-Care Corporation Design

ReOrder #	Description	Size	U/M
753155	Visco ConForm Cushion w/Coccyx Cutout	18"W x 16"D x 2"H*	Ea.

*Front Low Point

visit our website at www.skil-care.com

BARIATRIC CUSHIONS

TRI-FOAM GEL-INFUSED VISCO BARIATRIC CUSHION



Weight Tested to 550 lbs

Maximum comfort and stability

- Top layer of Gel-Infused Visco foam provides pressure relief
- Middle layer of high-density foam adds comfort and cushioning
- A **Firm Foundation** controls sling-seat hammocking
- Three-part design prevents bottoming-out
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction


Contoured Bottom - For Sling Seat Wheelchairs

ReOrder #	Description	U/M	ReOrder #	Description	U/M
754810	20"W x 16"D x 3"H	Ea.	754821	24"W x 20"D x 4"H	Ea.
754812	20"W x 18"D x 3"H	Ea.	754823	26"W x 18"D x 4"H	Ea.
754814	22"W x 16"D x 4"H	Ea.	754825	26"W x 20"D x 4"H	Ea.
754816	22"W x 18"D x 4"H	Ea.	754827	28"W x 18"D x 4"H	Ea.
754818	24"W x 18"D x 4"H	Ea.	754829	28"W x 20"D x 4"H	Ea.

Flat Bottom - For Solid Seat Wheelchairs

754833	20"W x 16"D x 3"H	Ea.	754843	24"W x 20"D x 4"H	Ea.
754835	20"W x 18"D x 3"H	Ea.	754845	26"W x 18"D x 4"H	Ea.
754837	22"W x 16"D x 4"H	Ea.	754847	26"W x 20"D x 4"H	Ea.
754839	22"W x 18"D x 4"H	Ea.	754849	28"W x 18"D x 4"H	Ea.
754841	24"W x 18"D x 4"H	Ea.	754851	28"W x 20"D x 4"H	Ea.

GEL-FOAM CUSHION



Weight Tested to 550 lbs

Double-chamber construction relieves pressure for heavy residents

- High-density foam prevents bottoming-out
- Low-Shear I launderable outer cover reduces skin damaging friction
- Incontinent-proof vinyl inner cover

ReOrder #	Description	Size	U/M
915105	Bariatric Thin-Line Gel-Foam	18"W x 18"D x 2"H	Ea.
915110	Bariatric Gel-Foam	20"W x 18"D x 3"H	Ea.
915115	Bariatric Thin-Line Gel-Foam	20"W x 18"D x 2"H	Ea.
915120	Bariatric Gel-Foam	22"W x 18"D x 3"H	Ea.
915125	Bariatric Thin-Line Gel-Foam	22"W x 18"D x 2"H	Ea.
915130	Bariatric Gel-Foam	24"W x 18"D x 3"H	Ea.
915131	Bariatric Thin-Line Gel-Foam	24"W x 18"D x 2"H	Ea.
915135		26"W x 18"D x 3"H	Ea.
915133		26"W x 20"D x 3"H	Ea.
915134	Bariatric Gel-Foam	28"W x 18"D x 3"H	Ea.
915137		28"W x 20"D x 3"H	Ea.
915139		30"W x 20"D x 3"H	Ea.

POSITION PLUS CUSTOMIZED BARIATRIC CUSHION



Weight Tested to 650 lbs

A total seating system for bariatric residents

- Available in sizes up to 28" x 28" (see insert on website for specific sizes)
- Custom-made to user's special seating requirements
- Gel flotation relieves interface pressure
- Reduces skin-damaging shear and friction
- Controls sliding and lateral leaning
- Incontinent-proof construction

All custom cushions and wedges are made to user specification and measurement. Call Skil-Care customer service at 800-431-2972 for complete ordering information.

CUSHIONS COVERS



LOW-SHEAR I



LOW-SHEAR II



SYNTHETIC SHEEPSKIN



TWO-TONE LOW-SHEAR II



TWO-TONE POLYESTER

WEDGED

Comes with a Non-Slip Bottom

Size	Low-Shear II (LSII)	UM
16"W x 16"D x 4" - 2"H	781040	Ea.
16"W x 16"D x 6" - 3"H	781041	Ea.
18"W x 16"D x 4" - 2"H	781042	Ea.
18"W x 16"D x 6" - 3"H	781043	Ea.

LEVEL

Size	Comes with a Non-Slip Bottom			Two-Tone Reminder Cover		UM
	Low-Shear I (LSI)	Low-Shear II (LSII)	Synthetic Sheepskin	LSII	Polyester	
16"W x 16"D x 1" - 2"H	781014	781028	781034	781054	781044	Ea.
16"W x 16"D x 3" - 4"H	781015	781029	781035	781055	781045	Ea.
18"W x 16"D x 1" - 2"H	781016	781030	781036	781056	781046	Ea.
18"W x 16"D x 3" - 4"H	781017	781031	781037	781057	781047	Ea.
20"W x 16"D x 1" - 2"H	781018	781032	781038	781058	781048	Ea.
20"W x 16"D x 3" - 4"H	781019	781033	781039	781059	781049	Ea.

SUPER SOFT HEAD PILLOW



This super soft pillow is made from a top of the line soft foam material. It can be used under the head, knees, ankles or anywhere you need to feel an elevated comfort.

- Super soft foam for added comfort, maintain support and prevent bottoming out
- Sold with either LSII stretched cover or vinyl heat-sealed cover
- Wipes clean for easy care
- Optional Cozy Cloth pillow case is available and sold separately



905023 Optional Cover

ReOrder #	Description	Size	U/M
905020	Super Soft Head Pillow w/LSII Cover	24"W x 18"D	Ea.
905021	Super Soft Head Pillow w/Vinyl Cover	24"W x 18"D	Ea.
905022	Optional Blue Cozy Cloth Pillow Case	24"W x 18"D	Ea.
905023	Optional Pink Cozy Cloth Pillow Case	24"W x 18"D	Ea.

INFINITY PILLOW



Optimum product for off-loading the heel, especially when the pressure and rubbing of a boot cannot be tolerated.

- Off loads both heels
- Prevents formation of heel pressure sores
- Allows air circulation to aid healing of existing ulcers
- Supports the lower legs
- Universally sized and fits most residents
- Visco-elastic foam redistributes pressure, therefore not causing high pressure points
- Low Shear II cover protects against skin damaging friction
- Wipes clean for easy care
- Can be used as a positioning pillow

ReOrder #	Description	Size	U/M
905025	Infinity Pillow	15" x 22" x 4"	Ea.

SKIL-CARE'S ALL NEW **Reminder** CUSHIONS

Reminder cushions comes with a two-color cover to help remind the resident and the caregiver to reposition. The cushions can be used on both sides.

SUPER SOFT ECONO GEL CUSHION




This cushion is designed to relieve pressure and increase comfort on any seating surface.

- Voted Best Value
- Super soft surface for maximum comfort
- Provides pressure relief and redistribution for all body types
- Can be used on either side
- Two toned cover serves as a reminder to reposition resident
- Low-Shear II stretch cover surrounding a gel-foam composition
- Water-resistant cover is wipe clean

ReOrder #	Description	Size	U/M
751166		16"W x 16"D x 2"H	Ea.
751167	Super Soft Econo-Gel Cushion w/ LSII Cover	18"W x 16"D x 2"H	Ea.
751168		20"W x 16"D x 2"H	Ea.

ECONO-GEL CUSHION




Pressure relief and redistribution at a budget price

- Two chamber gel pad designed for greater stability
- Comfort foam on top and bottom
- Gel pad between foam layers provides pressure-relieving floatation
- Available polyester outer cover reduces skin shear

ReOrder #	Description	Size	U/M
751154	Vinyl Only	16"W x 16"D x 2"H	Ea.
751150	Vinyl w/Poly Cvr		Ea.
751155	Vinyl Only	18"W x 16"D x 2"H	Ea.
751151	Vinyl w/Poly Cvr		Ea.
751157	Vinyl Only	20"W x 16"D x 2"H	Ea.
751153	Vinyl w/Poly Cvr		Ea.

GEL-INFUSED VISCO CUSHION TOPPER



- Gel-Infused Visco foam topper provides additional comfort
- Low-Shear II wipe-clean outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
914441	Gel-Infused Visco Cushion Topper	18"W x 16"D x 1"H	Ea.

SKIL-CARE'S ALL NEW **Reminder** CUSHIONS

CUSHIONS

REMINDER CUSHIONS

EZ DRY FOAM CUSHION




Reticulated foam cushion that can be hand washed

- Soft foam provides pressure relief and redistribution
- Stable surface for residents who have difficulty in maintaining upper trunk control
- Light weight
- Low-Shear II wipe-clean outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
753400	16"W w/LSII Cover		
753402	18"W w/LSII Cover	16"D x 2"H	Ea.
753404	20"W w/LSII Cover		
753406	16"W w/LSII Cover		
753408	18"W w/LSII Cover	16"D x 3"H	Ea.
753410	20"W w/LSII Cover		

COMFORT FOAM WHEELCHAIR CUSHION




Comfort and quality at reasonable prices

- High-density foam provides comfort and relieves pressure
- Extra-high E-Z Rise models support post-hip surgery care
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction

Note: Foam only products are guaranteed for 6 months.

ReOrder #	Description	ReOrder #	Description	Size	U/M
753232	Foam w/two color vinyl	-	-	16"W x 16"D x 2"H	Ea.
753272	Foam w/two color vinyl	-	-	16"W x 16"D x 3"H	Ea.
753230	Foam w/two color vinyl	753210	Foam w/two color vinyl and LSII Cover	18"W x 16"D x 2"H	Ea.
753270	Foam w/two color vinyl	753250	Foam w/two color vinyl and LSII Cover	18"W x 16"D x 3"H	Ea.
753271	Foam w/two color vinyl	753251	Foam w/two color vinyl and LSII Cover	18"W x 18"D x 3"H	Ea.
753290	E-Z Rise w/ Vinyl Cover	753285	E-Z Rise w/ LSII Cover	18"W x 16"D x 4"H	Ea.
707052	E-Z Rise w/ Vinyl Cover	-	-	18"W x 16"D x 6"H	Ea.
751188	Bimini Blue Vinyl Foam Cushion			18"W x 16"D x 2"H	Ea.
751189	Bimini Blue Vinyl Foam Cushion			18"W x 16"D x 3"H	Ea.

WEDGED FOAM WHEELCHAIR CUSHION



Controls sliding and provides pressure relief and redistribution

- High-density foam provides comfort and reduces pressure
- Incontinent-proof, wipe-clean vinyl inner cover
- Available polyester cover removes for easy laundering

ReOrder #	Description	Size	U/M
754016	Soft Foundation Vinyl Only	16"W x 16"D x 6"H - 3"H	Ea.
754039	Soft Foundation Vinyl Only	16"W x 16"D x 4"H - 2"H	Ea.
754020	Soft Foundation Vinyl Only		Ea.
754010	Soft Foundation Vinyl w/LSII	18"W x 16"D x 6"H - 3"H	Ea.
754021	Soft Foundation Vinyl w/Polyester		Ea.
754041	Soft Foundation Vinyl Only		Ea.
754040	Soft Foundation Vinyl w/LSII	18"W x 16"D x 4"H - 2"H	Ea.
754048	Soft Foundation Vinyl w/Polyester		Ea.
754043	Soft Foundation Vinyl Only	20"W x 16"D x 4"H - 2"H	Ea.

SKIL-CARE'S ALL NEW **Reminder** CUSHIONS

AIR LOCK CUSHION



This soft pre-inflated anti-bacterial vinyl cushion makes any wheelchair seat more comfortable and also provides pressure relief and redistribution.

- Decubitus prevention
- Air flotation (pre-inflated)
- Pressure relief and redistribution
- Vinyl, heat-sealed construction
- Anti-bacterial wipe clean vinyl
- Safety straps
- No valves to interfere with pressure relief
- Meets flammability requirements of CAL 117

ReOrder #	Description	Size	U/M
756245	Air Lock Seat Cushion w/ LSII Cover	19"W x 19"D	Ea.
756255		22"W x 28"D	Ea.

AIR LIFT CUSHION



Pre-inflated economical cushion

- Provides decubitus prevention and comfort
- Pressure relieving air flotation surface
- Easy to inflate or deflate (shipped pre-inflated)
- Wipe-clean anti-bacterial heat-sealed vinyl
- Optional Low-Shear II wipe-clean outer cover reduces skin damaging friction
- Comes with safety straps
- Meets flammability requirements of CAL 117

ReOrder #	Description	Size	Cushion Color	U/M
756215	Air Lift Seat Cushion	17"W x 17"D	Blue	Ea.
756230	w/LSII Cover	19"W x 19"D	Grey	Ea.
756237		22"W x 28"D	Green	Ea.

Cushion Covers

756219	Optional Two-Color Cushion Cover w/Straps	17"W x 17"D	Ea.
756234		19"W x 19"D	Ea.
756239		22"W x 28"D	Ea.

X-GEL OVERLAY



Adds gel flotation to any foam cushion


- Gel topper enhances foam cushions for more effective pressure relief and redistribution
- May be used over any similar sized cushion
- Water-based gel cools to prevent skin-damaging perspiration

ReOrder #	Description	Size	U/M
914429	X-Gel Overlay Vinyl Only	16"W x 16"D x 1"H	Ea.
914430		18"W x 16"D x 1"H	Ea.
781066	Low Shear II Cover for 914429, 16" x 16"	16"W x 16"D	Ea.
781065	Low Shear II Cover for 914430, 18" x 16"	18"W x 16"D	Ea.

WHEELCHAIR ACCESSORIES



703470



706050



703121



SkiL-Care offers a wide selection of wheelchair accessories that provide stability, positioning, support and comfort. Categories include forward leaning, lateral leaning, lumbar support, leg positioning, head control and skin protection. Additional accessories include bags & holders, trays, ponchos and aprons.

RECLINING WHEELCHAIR BACKREST



Economically converts any wheelchair to a recliner

- Provides up to 35-degrees of tilt
- Available in three heights (19", 24", 33") for 16", 18" and 20" wheelchairs
- Covered in durable, wipe-clean vinyl
- Can be power-washed with chair
- Easily customizable with six pockets (Foam supplied)

All sizes of Skil-Care's Reclining Wheelchair Backrests have six foam-padded compartments that can be customized to meet the resident's back alignment needs. These pockets may also be used to hold hot or cold gel packs.

ReOrder #	Description	Size	U/M
703111		16"W x 19"H	Ea.
703112		16"W x 24"H	Ea.
703113		16"W x 33"H	Ea.
703101		18"W x 19"H	Ea.
703102	Reclining WC Backrest	18"W x 24"H	Ea.
703103		18"W x 33"H	Ea.
703105		20"W x 19"H	Ea.
703106		20"W x 24"H	Ea.
703107		20"W x 33"H	Ea.
703199	Replacement Strap for Reclining WC Backrest	Universal	Ea.



ReOrder #	Description	Size	U/M
703130	Extended Full Arm - Armrest Pads	22"L	Pair
703131	Extended Half Arm - Armrest Pads	18"L	Pair

Optional armrest available to reduce pressure on sensitive arms and elbows. These pads are available for full and half-arm wheelchairs. Extended Armrest are recommended for patient comfort.

BACKREST SEAT COMBO



914373 & 914375 Only

- Provides added comfort and cushioning to any chair
- All models come with a fixed foam back pad that is 16"H x 1.5" thick
 - Low-Shear I launderable outer cover reduces skin damaging friction
 - Models 914366, 914368, 914443 & 914444 comes with an X-Gel cushion in seating area for improved pressure relief, cooling and comfort
 - Models 914373 and 914375 both have foam for the back and seat with a carrying handle

ReOrder #	Description	Seat Cushion Size	U/M
914366	Backrest Seat Combo w/Gel-Foam Vinyl	16"W x 16"D x 2.5"H	Ea.
914368	X-Cushion	18"W x 16"D x 2.5"H	Ea.
914443	Backrest Seat Combo w/Gel-Foam Vinyl	16"W x 16"D x 2.5"H	Ea.
914444	X-Cushion with Coccyx Cutout	18"W x 16"D x 2.5"H	Ea.
914372	Backrest w/pocket for 16" Cushion	16"W x 16"D	Ea.
914374	Backrest w/pocket for 18" Cushion	18"W x 16"D	Ea.
914373	Backrest Seat Combo w/Foam Seat Cushion	16"W x 16"D x 1.5"H	Ea.
914375	and carrying handle	18"W x 16"D x 1.5"H	Ea.

LATERAL SUPPORT ORTHOSIS



Helps lateral-leaning residents sit erect in wheelchair

- Contoured plastic frame provides firm support under axilla
- Foam backrest may be customized to meet resident's needs
- Restraint-free design uses no restrictive straps or belts
- Covered in durable, wipe-clean nylon
- For use with 18" wheelchairs



ReOrder #	Description	Height	Inner Width	U/M
702210		19"H Backrest	15"W	Ea.
702220	Lateral Support Orthosis	24"H Backrest	15"W	Ea.
702230		33"H Backrest	15"W	Ea.

AIR LIFT BACKREST CUSHION



751580 Hand Pump (Sold Separately)

Pre-inflated economical cushion

- Provides decubitus prevention and comfort
- Pressure relieving air flotation surface
- Pre-inflated, easy to inflate or deflate
- Wipe-clean anti-bacterial heat-sealed vinyl
- Available in a Low-Shear II wipe-clean outer cover reduces skin damaging friction
- Easily ties to handle bars
- Meets CAL 117 flammability requirements

ReOrder #	Description	Size	Color	U/M
756210	Air Lift Backrest Cushion Only	17"W x 17"D	Blue	Ea.
756215	Air Lift Backrest Cushion w/LSII Cover			
756219	Optional Cushion Cover w/Straps	19"W x 19"D	Grey	Ea.
756225	Air Lift Backrest Cushion Only			
756230	Air Lift Backrest Cushion w/LSII Cover			
756234	Optional Cushion Cover w/Straps			
751580	Hand Pump			Ea.

WHEELCHAIR BACKREST PAD



A portable foam pad that cushions the back and any seat for a 16" or 18" wheelchair

- Cushions the back or seat with 1.5" of foam
- Covered with vinyl
- Fits any 16" or 18" wide wheelchair

ReOrder #	Description	Size	U/M
703063	Wheelchair Backrest Pad	Fits 16" - 18"W Wheelchairs	Ea.

SEAT & BACK AIR CUSHION



This combination seat and back air cushion is foldable and comes pre-inflated

- Provides decubitus prevention and comfort
- Pressure relieving air flotation surface
- Pre-inflated, easy to inflate or deflate for portability
- Easily ties to handle bars
- Wipe-clean anti-bacterial heat-sealed vinyl

ReOrder #	Description	Size	U/M
703300	Seat & Back Air Cushion	34"L x 17"W	Ea.

WHEELCHAIR COZY SEAT




Foam padding provides seating comfort

- Protects against skin tears and relieves pressure
- Covered with a soft, launderable cozy cloth fabric
- Models secure to either 16"x16" or 18"x16" wheelchairs
- Easy installation with loops and hook & loop straps

ReOrder #	Description	Size	U/M
703005	Wheelchair Cozy Seat	Fits 16"W Wheelchairs	Ea.
703006	Wheelchair Cozy Seat	Fits 18"W Wheelchairs	Ea.

LIFT-AWAY WHEELCHAIR TRAY



912410



914504



Hook & Loop Strips

Armrest Pad

Hook & Loop Attachment and Armrest Pads



705014 LSII Cover

Easily removed by most residents

- Helps achieve proper body position, balance and alignment
- Padded to provide a comfortable, low-pressure surface for arms/elbows
- Water-based gel surface provides tactile stimulation and the stars create visual appeal (model 912410)
- Hook & Loop strips secure tray to wheelchair armrests armrest pads provided
- Covered in durable, wipe-clean vinyl

ReOrder #	Description	Size	U/M
705015	SofTop Lift-Away Wheelchair Tray w/ Vinyl Cover	Fits 16" - 18" Wheelchairs	Ea.
705025	SofTop Lift-Away Wheelchair Tray w/ Super Soft LSII Cover	Fits 20" - 22" Wheelchairs	Ea.
705014	SofTop Lift-Away Wheelchair Tray w/ Super Soft LSII Cover	Fits 16" - 18" Wheelchairs	Ea.
705026	SofTop Lift-Away Wheelchair Tray w/ Super Soft LSII Cover	Fits 20" - 22" Wheelchairs	Ea.
912410	Gel Top Sensory Tray, Black Vinyl, Lift-Away	Fits 16" - 18" Wheelchairs	Ea.
914504	Gel Top Sensory Tray, White Vinyl, Lift-Away	Fits 16" - 18" Wheelchairs	Ea.
705016	Replacement Armrest Sleeves	Fits 16" - 18" Wheelchairs	Ea.

WHEELCHAIR LAP TRAY WITH STRAPS



Gel Top Tray 912415



Easy Release Hook & Loop Strap

- Foam-padded surface is gentle on arms and elbows
- Helps achieve proper body position, balance, and alignment
- Padded to provide a comfortable, low-pressure surface for arms/elbows
- Water-based gel surface provides tactile stimulation and the stars create visual appeal (model 912415)
- Around-the-chair strap prevents resident from sliding forward
- Easy-open Hook & Loop closures
- Larger size fits wheelchairs from 24" to 28"
- Available with new nylon or standard vinyl cover



705018 LSII Cover

ReOrder #	Description	Size	U/M
705021	SofTop Wheelchair Tray w/Vinyl Cover	Fits 14" Wheelchairs	Ea.
705018	SofTop Wheelchair Tray w/Super Soft LSII Cover	Fits 16" - 18" Wheelchairs	Ea.
705027	SofTop Wheelchair Tray w/Super Soft LSII Cover	Fits 24" - 28" Wheelchairs	Ea.
705020	SofTop Wheelchair Tray w/Vinyl Cover (as shown)	Fits 16" - 18" Wheelchairs	Ea.
705022	SofTop Wheelchair Tray w/Vinyl Cover	Fits 20" - 22" Wheelchairs	Ea.
705019	SofTop Wheelchair Tray w/Nylon Cover	Fits 16" - 18" Wheelchairs	Ea.
705024	SofTop Wheelchair Tray w/Nylon Cover	Fits 24" - 28" Wheelchairs	Ea.
912415	Gel Top Sensory Tray, Black Vinyl, Hook & Loop Strap	Fits 16" - 18" Wheelchairs	Ea.
705091	Replacement Straps for Wheelchair Tray	—	Set

WHEELCHAIR FLIP TRAY



Skil-Care
CONSUMER DESIGN

Padding



E-Z exiting arm tray

- Wheelchair tray flips out of the way for easy exiting
- Hard, stable, soil and scratch-resistant plastic
- Side and armrest padded when in the down or flipped position
- Two cup holders for drinks, eye glasses or small items

ReOrder #	Description	Size	U/M
705031	Right	23"L x 19.5"W x 0.25"D	Ea.
705032	Left	0.25"D	Ea.

LAPTOP CUSHION



Provides safety, positioning, and slide control

- Polyfoam padded
- Covered in durable, wipe-clean nylon
- Available for both full and half-arm wheelchairs

For 16"-18"W Full Arm Wheelchairs

ReOrder #	Description	U/M
307012	2.5" Thick, LapTop w/Cutouts - Black	Ea.
307014	3.25" Thick, LapTop w/Cutouts - Black	Ea.
307010	4" Thick, LapTop w/Cutouts - Black	Ea.
307020	4" Thick, LapTop w/o Cutouts - Black	Ea.

For 16"W Half Arm Wheelchairs

ReOrder #	Description	U/M
307032	2.5" Thick, LapTop w/Cutouts - Black	Ea.
307034	3.25" Thick, LapTop w/Cutouts - Black	Ea.
307030	4" Thick, LapTop w/Cutouts - Black	Ea.

For 18"W Half Arm Wheelchairs

ReOrder #	Description	U/M
307042	2.5" Thick, LapTop w/Cutouts - Black	Ea.
307044	3.25" Thick, LapTop w/Cutouts - Black	Ea.
307040	4" Thick, LapTop w/Cutouts - Black	Ea.

LIFT-OFF LAP CUSHION



Note: Lift-Off Lap Cushion is 3.5"H

Hook & Loop closures with red loops enable release by sight or hand-impaired residents.

Provides restraint-free safety and positioning

- Releases easily as resident stands
- Held in place by easy-to-apply Hook & Loop closures
- Foam padded for resident comfort
- Covered in durable, wipe-clean nylon

ReOrder #	Size	U/M
307050	Fits 16"-18" Wheelchairs	Ea.
307051	Fits 20"-22" Wheelchairs	Ea.

WHEELCHAIR POSTURE SUPPORT



Promotes correct sitting posture

- Limits forward and lateral leaning
- Secures with easy-open Hook & Loop closures
- Available Hook & Loop with red finger grip closures for hand-impaired residents

Hook & Loop

ReOrder #	Description	U/M
610111	Small-Medium, Hook & Loop	Ea.
610112	Medium-Large, Hook & Loop	Ea.

Hook & Loop w/Red Loop

610115	Small-Medium, Hook & Loop w/Red finger grip	Ea.
610116	Medium-Large, Hook & Loop w/Red finger grip	Ea.

WHEELCHAIR SIDE WINGS



Comfortable support for leaning resident

- Foam-padded, rigid side wings attach easily to most wheelchairs
- Covered with durable, wipe-clean vinyl
- X-tra Support Side Wings include armrest pads for additional stability and comfort

ReOrder #	Description	U/M
706320	Wheelchair Side Wings - Standard Support	Pair
706322	Wheelchair Side Wings - X-Tra Support (as shown)	Pair

ADJUSTABLE LATERAL SUPPORT



Easy adjustment to accommodate a range of body widths.



Maximum support is provided by width-adjustable side panels

- Minimizes lateral leaning in wheelchairs and most other chairs
- Foam-padded sides adjust laterally for a custom fit
- Will not interfere with resident's ability to reach wheelchair hand-rims
- Wipes clean for easy care
- Available in three sizes for petite to bariatric residents
- All models are 7"H x 9"D

ReOrder #	Description	Adjustment Range	U/M
706050	Adjustable Lateral Support Small/Petite	Width adjusts from 10"-15"	Ea.
706055	Adjustable Lateral Support Medium	Width adjusts from 12"-19"	Ea.
706059	Adjustable Lateral Support Large/Bariatric	Width adjusts from 18"-30"	Ea.

ONE-WAY GLIDE SYSTEM



Continuous Loop



Single Layer w/Handles



The Skil-Care One-Way Glide systems are designed to restrict patient movement in one direction. The Unique material helps to maintain proper seating position and allows for easy repositioning of the patient.

- Promotes proper posture and provides additional comfort
- Reduces staff time and helps reduce caregiver back strain while repositioning
- Can be used on any seating surface
- Eliminates the need for disposable anti-slip material

Available in two Styles: Continuous loop style and Single layer with handles

ReOrder #	Description	Size	U/M
556080	Continuous Loop	17"L x 17"W	Ea.
556081	Single Layer w/Handles	50"L x 15.5"W	Ea.

CHAIR GLIDE SYSTEM



Chair Glide makes it easier to safely glide the patient to an upright seat position without lifting.

- Promotes proper posture and provides additional comfort
- Comes with disposable overlay sheet which is moisture absorbing and breathable
- Reduces staff time for repeat positioning
- One-Way glide restricts forward movement
- Patient glides into position without lifting
- Glide pad is wipe-clean or launderable

ReOrder #	Description	Size	U/M
556084	Chair Glide System (Glide Pad & Overlay Sheet)	17"L x 17"W	Set
556085	Chair Glide Pad Only	17"L x 17"W	Ea.
556086	Disposable Cozy Cloth Overlay Sheet	22"L x 22"W	12/pk

LATERAL/LUMBAR SUPPORT



Replacement Gel Packs



Gel Pads Positioned for Lumbar Support

Comfortable gel packs control lateral leaning or provide lumbar support

- Position two gel packs vertically to reduce lateral leaning
- Position one gel pad horizontally to provide lumbar support
- Firm back support distributes pressure caused by backrest hammocking
- Gel packs and backrest are covered in wipe-clean vinyl
- Secures to any 16-inch or 18-inch wheelchair, geri-chair or utility chair

ReOrder #	Description	Size	U/M
706350	Lateral Support w/Vinyl Cover (w/2 8"L Gel Pads)	Fits 16"-18" WC	Ea.
706359	Optional Low Shear II Cover (Fits over entire backrest)	Fits 16"-18" WC	Ea.
706351	Replacement Gel Pack	8"L x 4"W	Ea.
706352		4"L x 4"W	Ea.

SNUG SUPPORT



- Controls lateral leaning in any armrest chair
- Fills the gap between resident and armrest
 - Effective for promoting proper sitting posture
 - Foam padded for resident comfort
 - Covered in durable, wipe-clean nylon
 - Provides hip guide for pelvic alignment
 - Can be used on either side, the wheelchair
 - Can be used in all wheelchair widths

ReOrder #	Description	Size	U/M
706309	Snug Support	16"L x 6"H (front) x	Ea.
706310		15"H (back) x 1.5" Thick	2/pk

UNIVERSAL CHAIR SIDE WINGS



- Comfortable lateral positioning for all chairs
- Improves sitting posture when chair is in upright or reclined position
 - Prevents sideways movement when chair is reclined
 - Foam padded for comfort
 - Low-Shear II wipe-clean outer cover that reduces skin damaging friction
 - Secures to all wheelchairs, geri-chairs and lounge chairs

ReOrder #	Description	Size	U/M
706319	Universal Chair Side Wings	18.5"W x 10"D x 10"H	Ea.

LATERAL BODY SUPPORT



Controls lateral leaning in wheelchair or geri-chair

- Vertical support stabilizes torso
- Foam-cushioned armrest trough supports arm
- Arm trough has a Hook & Loop strap which steadies arm for support
- May be used on left or right side
- Rigid model numbers have a plastic insert to maintain rigidity



ReOrder #	Description	Size	U/M
706010	Body Support Rigid Insert w/Sheepskin Cover	12 3/4"L	Ea.
706011	Body Support Rigid Insert w/Polyester Cover	12 3/4"L	Ea.
706012	Body Support Rigid Insert w/LSII Cover (as shown)	12 3/4"L	Ea.
914210	Body Support Rigid Insert w/LSII Cover/Gel Pad	12 3/4"L	Ea.
706017	Body Support Rigid Insert w/Nylon Cover	14"L	Ea.
706015	Body Support Rigid Insert w/Sheepskin Cover (as shown)	14"L	Ea.
706016	Body Support Rigid Insert w/Polyester Cover	14"L	Ea.
706020	Body Support Rigid Insert w/Sheepskin Cover Right Arm (Elevated)	12 3/4"L	Ea.
706021	Body Support Rigid Insert w/Sheepskin Cover Left Arm (Elevated)	12 3/4"L	Ea.

ReOrder #	Description	Size	U/M
914394	Lateral Support Not Rigid w/Sheepskin	14"L	Ea.
914395	Lateral Support Not Rigid w/Nylon & Sheepskin Cover	14"L	Ea.

EZ ON LATERAL SUPPORT




- Convenient slip-on design
- Controls lateral leaning by stabilizing torso
 - Foam cushioned for comfort and pressure relief
 - **All foam** arm trough with Hook & Loop straps steadies arm for support
 - May be used on left or right side
 - Covered in durable, wipe-clean nylon
 - Universally sized

ReOrder #	Description	U/M
706005	EZ On Lateral Support with Nylon Cover	Ea.

LATERAL STABILIZER ARM TROUGH



706140 Stabilizer w/Polyester Cover



914220 Arm Tray with Gel



Converts a standard armrest to a stable platform

- Arm trough stabilizes arm when resting on the wheelchair armrest
- Foam cushioned for comfort and pressure relief
- Helps stabilize upper trunk through added support
- Can be used on either side of the wheelchair
- Conveniently slips onto the armrest without straps
- Available with a gel pad for additional pressure relief

ReOrder #	Description	Size	U/M
706140	Lateral Stabilizer Arm Trough - Polyester	16"L x 5"H x 6"W	Ea.
706141	Lateral Stabilizer Arm Trough - Synthetic Sheepskin Cover	16"L x 5"H x 6"W	Ea.
706142	Lateral Stabilizer Arm Trough - LSII Cover	16"L x 5"H x 6"W	Ea.
914220	Lateral Stabilizer Arm Trough w/Gel Pad LSII Cover	16"L x 5"H x 6"W	Ea.

LATERAL STABILIZER ARMREST BOLSTERS



Shown with two Bolsters



Reduces lateral leaning

- Cushions wheelchair armrests
- Available for full-arm and half-arm wheelchairs in 5" or 6" diameters
- Covered in a polyester, nylon or synthetic sheepskin fabric

ReOrder #	Description	Arm	U/M
706110	5" w/Polyester Cover		Ea.
706111	5" w/Sheepskin Cover		Ea.
706112	5" w/Nylon Cover	Full Arm	Ea.
706120	6" w/Polyester Cover		Ea.
706121	6" w/Sheepskin Cover		Ea.
706122	6" w/Nylon Cover		Ea.

ReOrder #	Description	Arm	U/M
706114	5" w/Polyester Cover		Ea.
706115	5" w/Sheepskin Cover		Ea.
706116	5" w/Nylon Cover	Half Arm	Ea.
706124	6" w/Polyester Cover		Ea.
706125	6" w/Sheepskin Cover		Ea.
706126	6" w/Nylon Cover		Ea.

SWING-AWAY FOOT SUPPORT



Original Skil-Care CORPORATION Design

Bariatric Available

- Swings away for easy wheelchair transfers
- Cushioned side and back panels control foot positioning
 - Covered in durable, wipe-clean nylon
 - Easy installation
 - Available for standard size and bariatric wheelchairs
 - Specify left, right, or pair when ordering
 - Fits most wheelchairs

ReOrder #	Description	Size	U/M
703470	Left		Ea.
703471	Right	Fits 16" - 18"W Wheelchair	Ea.
703472	Right and Left		Pair
703475	Left		Ea.
703476	Right	Fits 22" - 30"W Wheelchair	Ea.
703477	Right and Left		Pair


WHEELCHAIR FOOTREST EXTENDER W/ LEG SEPARATION



- Helps control foot drop
- Prevents legs from slipping behind the footrest
- Protects feet from contacting wheelchair metal surfaces
- Separates the legs for optimal hip alignment with either a 1.5" or 3" leg separator
- Removable foot separator that converts the Footrest Extender to a standard footrest
- Covered in durable wipe-clean vinyl
- Two quick release buckle straps allow for easy attachment to most 16", 18" or 20" wide wheelchairs


ReOrder #	Description	Size	U/M
703262	Footrest Extender with 1.5" separator	19"W x 11"D x 9.5"H	Ea.
703264	Footrest Extender with 3" separator	19"W x 11"D x 9.5"H	Ea.
703267	Optional Separator 1.5" Thick	8"H x 7"D x 1.5" Thick	Ea.
703269	Optional Separator 3" Thick	8"H x 7"D x 3" Thick	Ea.

FOOT CRADLE



- Controls flexion contracture of the legs
- Accommodates existing flexion contracture
 - Prevents foot drop
 - Cushioned back and side panels prevent feet from slipping off footrests
 - Includes one set of 1-inch-thick pads
 - Covered in durable, wipe-clean vinyl
 - Available Contracture Accommodation Kit addresses a wide range of flexion contracture needs

Contracture Accommodation Kit Includes:
 (2 ea.) 1" x 7" x 7" pads
 (1 ea.) 2" x 7" x 7" pad
 (1 ea.) 3" x 7" x 7" pad



ReOrder #	Description	Size	U/M
703410	Foot Cradle	Fits 16"-18"W W/C	Ea.
703416	Contracture Accommodation Kit	Fits 16"-18"W W/C	Set (Kit)

TWO-PIECE FOOTREST EXTENDER



Effective for lower extremity positioning

- Raises footrest for improved posture and seating alignment
- Prevents foot drop and feet from slipping behind footrests
- Cushioned back protects lower leg against contact with wheelchair
- Secures easily with quick-release straps
- Covered in durable, wipe-clean vinyl
- Removable foot pad for flexible height selection

ReOrder #	Description	Wheelchair Size	Foot Pad Thickness	U/M
703280		16" - 18"W	1 Inch	Ea.
703281		16" - 18"W	2 Inch	Ea.
703282		16" - 18"W	3 Inch	Ea.
703283		16" - 18"W	4 Inch	Ea.
703210		18" - 20"W	1 Inch	Ea.
703220		18" - 20"W	2 Inch	Ea.
703230	Footrest Extender	18" - 20"W	3 Inch	Ea.
703240		18" - 20"W	4 Inch	Ea.
703250		20" - 22"W	1 Inch	Ea.
703251		20" - 22"W	2 Inch	Ea.
703252		20" - 22"W	3 Inch	Ea.
703253		20" - 22"W	4 Inch	Ea.

TWO-PIECE FOOTREST EXTENDER REPLACEMENT PADS

ReOrder #	Wheelchair Size	Foot Pad Thickness	U/M
703291		1 Inch	Ea.
703292	16" - 18"W/C	2 Inch	Ea.
703293		3 Inch	Ea.
703295		4 Inch	Ea.
703271		1 Inch	Ea.
703272	18" - 20"W/C	2 Inch	Ea.
703273		3 Inch	Ea.
703274		4 Inch	Ea.
703256		1 Inch	Ea.
703257	20" - 22"W/C	2 Inch	Ea.
703258		3 Inch	Ea.
703259		4 Inch	Ea.

SIDE-KICK



Side-Kick shown with Footrest Extender

Prevents feet from slipping off sides

- Attaches to all Skil-Care two-piece footrest extender devices*
- Accepts footrest pad shown on right
- Comfortable, foam-padded side panels
- Covered in durable, wipe-clean vinyl

ReOrder #	Description	Size	U/M
703455	Side-Kick	Fits 16"-18"W WC	Ea.

* Does not fit on the One-Piece models

ONE-PIECE FOOTREST EXTENDER



Effective for lower extremity positioning

- Raises footrest for improved posture and seating alignment
- Prevents foot drop and feet from slipping behind footrests
- Cushioned back protects lower leg against contact with wheelchair
- Secures easily with quick-release straps
- Covered in durable, wipe-clean vinyl

ReOrder #	WC Size	Foot Pad Thickness	U/M
703286	Fits	1 Inch	Ea.
703288	16" - 18" W/C	2 Inch	Ea.
703294	Fits	1 Inch	Ea.
703296	20" - 24" W/C	2 Inch	Ea.
703297		3 Inch	Ea.

WHEELCHAIR LEG PAD



Protects legs from contact with metal wheelchair parts

- New E-Z clip buckles and straps simplify attachment
- Covered in durable, wipe-clean vinyl
- Prevents resident's feet from slipping behind footrests
- Foam-padded, rigid-backed leg pad

ReOrder #	Description	Size	U/M
703070		Fits 16" - 18"W Wheelchairs	Ea.
703071	Wheelchair Leg Pad	Fits 18" - 20"W Wheelchairs	Ea.
703072		Fits 20" - 24"W Wheelchairs	Ea.

ABDUCTION WEDGE



Comfortable foam wedge controls abduction contracture

- Held securely in place with adjustable Hook & Loop straps (each strap is 25"L)
- Wipe-clean vinyl cover
- Promotes proper thigh alignment

ReOrder #	Description	Size	U/M
703075	Abduction Wedge	7"W x 8"H x 4"D	Ea.

GULL WING LEG ABDUCTOR



- Accommodates abduction contracture
- Promotes proper thigh/leg alignment
- Held securely in place while seated with Hook & Loop straps
- Wipe clean plastic separator
- Can remove straps and hand wash

ReOrder #	Description	Size	U/M
703077	Gull Wing Leg Abductor	Universal	Ea.

WHEELCHAIR LEGREST W/PADDED SIDES

- Prevents the legs from falling behind the footrests
- Controls the legs from falling off the sides
- Foam padded and prevents bruising from contact with the wheelchair
- Using side release buckles, it easily attaches to all wheelchairs
- Covered with durable, wipe-clean vinyl



Note: Fit most 16", 18" and 20" wide wheelchairs.

ReOrder #	Description	Size	U/M
703460	Legrest w/Padded Sides	20"W x 11"D x 8"H	Ea.

SYNTHETIC SHEEPSKIN SLIPPERS



Non-Skid Bottom



Wheelchair Attachment Slipper

Positions legs on footrest and provides a non-slip walking slipper

- Two models available: one attaches to the footrest and the other doesn't
- Protects legs while resting on wheelchair footrest
- Positions foot on leg rest for improved posture
- Soft synthetic sheepskin that decreases restlessness through added comfort
- Non-detachable model serves as a non-skid slipper

ReOrder #	Description	Size	U/M
703350	Relief Slipper for footrest	Small-Medium	Pair
703351		Medium-Large	Pair
703352	Relief Slipper for Walking	Small-Medium	Pair
703353		Medium-Large	Pair

CALF PAD COVER



Original Skil-Care Corporation Design



Reduces pressure and skin tears

- Temperature-sensitive foam padding lowers pressure on calves
- Protects against tears and abrasions
- Low-Shear II wipe-clean outer cover reduces skin damaging friction
- Hook & Loop straps secure to all wheelchair calf pads

ReOrder #	Description	U/M
703068	Calf Pad Cover	Pair

SKIN-GUARD LEG PROTECTOR



Protects legs against skin tears and bruises

- Soft, durable neoprene pads wrap around wheelchair footrest-support bars
- Hook & Loop tabs provide easy and secure application
- Designed to fit all wheelchairs
- Can be used on products with 1" tubular surfaces

ReOrder #	Description	U/M
703073	Skin-Guard Leg Protector	Pair

LEG BOLSTERS



Skil-Care Leg Bolsters are designed to guard against skin tears due to sharp and/or protruding edges around footrest/legrest hanger brackets.

- Low-Shear II wipe-clean outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
703060	Leg Bolster	4"D x 6.5"H	Pair
703061	Leg Bolster	4"D x 6.5"H	16 Pairs/ Cs



LATERAL/LUMBAR SUPPORT



Gel Pads Positioned for Lateral Support



Gel Pads Positioned for Lumbar Support




Comfortable gel packs control lateral leaning or provide lumbar support

- Position two gel packs vertically to reduce lateral leaning
- Position one gel pad horizontally to provide lumbar support
- Can be moved vertically to adjust height of lumbar support
- Soft back support distributes pressure caused by backrest hammocking and adds to patients comfort
- Gel packs and backrest are covered in wipe-clean vinyl
- Secures to any 16-inch or 18-inch wheelchair, geri-chair or utility chair

ReOrder #	Description	Size	U/M
706350	Lateral Support w/Vinyl Cover (w/2 8"L Gel Pads)	Fits 16"-18" Wheelchairs	Ea.
706359	Optional Low Shear II Cover (Fits over entire backrest)	Fits 16"-18" Wheelchairs	Ea.
706351	Replacement Gel Pack	Full (8"L x 4"W)	Ea.
706352	Replacement Gel Pack	Half (4"L x 4"W)	Ea.

ADJUSTABLE THICKNESS LUMBAR SUPPORT



Easy to Remove Foam Insert



Universal, adjustable thickness lumbar support

- Removable foam pad that converts the lumbar support to either a 2" or 3" thickness
- Easily attaches to any chair with a quick release belt attachment
- Full 17" long and 8" wide
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction
- Non-skid back to ensure proper positioning

ReOrder #	Description	Size	U/M
706365	Modular Lumbar Support - Half Round	17"L x 8"W x 2" thick (w/o foam insert) or 3" thick (w/foam insert)	Ea.

FULL ROUND LUMBAR SUPPORT

Universal, round lumbar support

- Easily attaches to any chair with a quick release belt attachment
- Full 16" long and 5" diameter
- Covered in wipe-clean vinyl



ReOrder #	Description	Size	U/M
703066	Full Round Lumbar Support	16"L x 5"Dia.	Ea.

HEAD CONTROL

STANDARD HEADREST

Cushions and positions resident's head

- Can be used for all high-back wheelchairs, geri-chairs, and recliners
- Foam support provides cushioned comfort
- Polyester cover removes for easy laundering
- Available in three depths

ReOrder #	Description	Size	U/M
703115		3" Depth	Ea.
703116	Standard Headrest	3.75" Depth	Ea.
703117		4.75" Depth	Ea.



ADJUSTABLE HEAD POSITIONER W/GEL

Prevents hyperextension of the neck

- Head support adjusts in width to customize head placement
- Water-based gel packs relieve pressure on ears
- Can be used for all high-back wheelchairs, geri-chairs, and recliners
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction

ReOrder #	Description	U/M
914364	Adjustable Head Positioner w/Gel	Ea.
503099	Replacement Gel Pack	Ea.



Gel Inserts



TRIANGULAR HEAD POSITIONER

Wedge shape provides comfortable lateral head support

- Designed to prevent hyperextension of the neck
- Can be used for all high-back wheelchairs, geri-chairs and recliners
- Cushioned with resilient foam
- Launderable polyester outer cover

ReOrder #	Description	U/M
703121	Triangular Head Positioner	Ea.



ARMREST WHEELCHAIR CUSHION WITH STORAGE POUCH



Dual purpose armrest pad comes with storage pouch

- Prevents arm pressure sores and skin abrasion
- Storage pouch (1) on one side only for personal items
- Armrest pad with pouch can be used on either side of wheelchair

ReOrder #	Description	Size	U/M
703018	Armrest Storage Pouch with two arm pads (One pouch per set)	16"L	Pair
703019	Armrest Storage Pouch with two arm pads (One pouch per set)	11"L	Pair

MOBILE ARM SUPPORT



Level w/Gel Pad (Model 914230 and 914235)



- Rotational arm support for contracted-arm residents
- Especially effective for contracted-arms
 - Adjustable-angle capability to adapt to resident's specific needs
 - Easily installs to any half-arm or full-arm wheelchair
 - Level or elevated-arm models allow for full hand extension
 - Low-Shear II wipe-clean outer cover reduces skin damaging friction
 - Fits most wheelchairs

ReOrder #	Description	Arm	U/M
706210	Level		Ea.
706220	Elevated	Right	Ea.
914230	Level w/Gel Pad		Ea.
706215	Level		Ea.
706225	Elevated	Left	Ea.
914235	Level w/Gel Pad		Ea.

ARMREST PADS



Full Arm



Extended Full Arm



Sheepskin Cover




Half Arm

Comfort and pressure relief protection for arms and elbows

- Foam-padded and covered in nylon or synthetic sheepskin
- Hook & Loop fasteners for easy application
- Available for full-arm and half-arm wheelchairs
- For added protection, choose the extended full-arm model
- Fits most standard open-arm chairs

ReOrder #	Description	Size	U/M
703010	Full Arm - Sheepskin Cover	15" L	Pair
703015	Full Arm - Nylon Cover	15" L	Pair
703016	Half Arm - Nylon Cover	11" L	Pair
703130	Extended Full Arm - Nylon Cover	22" L	Pair
703131	Extended Half Arm - Nylon Cover	18" L	Pair

SHEEPSKIN COVERINGS



Added comfort and skin protection for wheelchair-bound residents

- High-pile, non-allergenic synthetic sheepskin
- Adds a touch of luxury to any wheelchair
- Stands up to frequent institutional laundering
- Fits 16"-18" wheelchair sizes

ReOrder #	Description	U/M
703020	Wheelchair Leg Pads	Ea.
703010	Wheelchair Armrest Pads	Pair
703040	Wheelchair Footrest Pads	Pair
703050	Wheelchair Seat/Backrest Pads	Set

LAP BLANKET W/ HAND WARMER



Synthetic Sheepskin Lined Bottom

Traps heat in for warmer hands and feet

- Covers entire lap and is universally sized
- Provides warmth to lap, hands and feet
- Two hand warmer pockets and one center utility pocket
- Soft cozy cloth fabric and strategically lined with synthetic sheepskin

ReOrder #	Description	Size	U/M
914761	Lap Blanket w/ Hand Warmer	Universal	Ea.

RAIN CAPE FOR WHEELCHAIR



The rain cape carrying case attaches easily to the wheelchair handles.

Keeps residents and wheelchair dry

- Easy-to-use hooded cape covers resident and wheelchair
- Made from durable, lightweight, waterproof nylon
- Universal size fits small to large users
- Dries quickly, wipes clean or is launderable for easy care

ReOrder #	Description	U/M
703080	Rain Cape for Wheelchair with Carrying Case	Ea.

SHOWER PONCHO



Keeps resident warm before and after the shower

- Provides modesty and warmth for added comfort
- Use before and after bathing or showering
- Recommended for use when washing and styling resident's hair
- Durable and completely launderable
- Cozy cloth on the inside

ReOrder #	Description	Size	U/M
909120	Short Back - No Hood	34" long (front) 23" long (back)	Ea.
909130	*Long Back - No Hood	34" long (front) 34" long (back)	Ea.
909140	*Long Back - With Hood	34" long (front) 34" long (back)	Ea.
909150	Short Back - With Hood	34" long (front) 23" long (back)	Ea.

* Recommended, especially when resident is ambulated



MODESTY APRON




Provides warmth and dignity to wheelchair residents

- Fits over residents legs to provide warmth and a modesty covering for female residents
- Extra wide length and width to fit most residents
- Two convenient pockets to hold personal items
- Easily attaches with waist tie straps and can be used on any chair
- Super-soft cozy cloth, launderable micro-fiber fabric

ReOrder #	Description	Size	U/M
707030	Modesty Apron	38"W x 33"L	Ea.

ACTIVITY AIDS



Activity Aids come in three styles: Apron, Vest and Overlay. These items are designed for patients with decreased cognitive function and the need for tactile stimulation.

- Helps redevelop fine motor skills
- Improves hand-eye coordination
- Helps improve dressing skills, dexterity and attention span
- Provides color and textures stimulation
- Attaches to tray, lap or worn as a vest depending on the model type.
- Pockets can hold small photos

ReOrder #	Description	Size	U/M
912461	Activity Apron	Small-Medium	Ea.
912460	Activity Apron	Large	Ea.
912462	Activity Apron w/Hook & Loop Pocket	Large	Ea.
912464	Activity Vest	Small-Medium	Ea.
912465	Activity Vest	Large	Ea.
912467	Activity Overlay	28"W x 11"D	Ea.

OXYGEN CYLINDER HOLDER



- Fits all push-handle wheelchairs
- Durable nylon for long-lasting quality
 - Easy to insert and remove tank
 - Adjustable straps accommodate most wheelchair widths
 - Converts to hold "D" and "E" size cylinders
 - Convenient carry handle

ReOrder #	Description	Size	U/M
707027	Single	Fits 16" - 24" Wheelchair	Ea.
707026		Fits 16" - 24" Wheelchair	10/pk
707028	Double	Fits 16" - 24" Wheelchair	Ea.
707029		Fit all Geri-Chairs with handles	Ea.

WHEELCHAIR FOOTREST BAG



- Prevents misplacement of wheelchair footrests
- Fits most footrests and wheelchairs
 - Easily hangs off wheelchair
 - Covered in durable, wipe-clean vinyl

ReOrder #	Description	Size	U/M
914362		22"H x 14"W	Ea.
914161	Footrest Bag	22"H x 16"W	Ea.
914363		22"H x 18"W	Ea.

CHART HOLDER



- 16"W x 14"H bag that houses critical medical charts or other flat items
- Two loop-straps to hang off the wheelchair
- 5.5" x 6.5" clear window to view the bag contents

ReOrder #	Description	Size	U/M
914341	Chart Holder	14"W x 16"H	Ea.

UNIVERSAL STORAGE BAG



707010



914393

These bags are designed to attach to the walker or wheelchair with adjustable straps that have side release buckles. The pockets will house both small and larger personal items.

- Pockets for housing personal items.
- Attaches to the wheelchair or walker with 1" webbing straps and plastic side release buckles
- Wipe clean vinyl

ReOrder #	Description	Size	U/M
707010	Universal Storage Bag w/Three pockets	11.5"W x 3.5"D x 11"H	Ea.
914393	Universal Storage Bag w/One pocket	13.5"W x 3.5"D x 11"H	Ea.

WHEELCHAIR/ WALKER HANDY BAG



Multi-pocket bag that allows residents to carry personal items

- Assorted sized pockets
- Includes name tag and key ring holder
- Easy installation on wheelchair or walker

ReOrder #	Description	Size	U/M
914396	Handy Bag	18"W x 11"L	Ea.

CHAIRPACK



A convenient storage bag for your personal items

- Stores books, newspapers and personal items
- ChairPack has a foam back for added comfort
- Easily fits on back of any chair
- Has easy carrying handle
- Wipes-clean
- Clear window for a name tag
- Pencil, pen holder closes with a zipper

ReOrder #	Description	Size	U/M
707035	ChairPack	14"W x 18"H	Ea.

CARPACK



A convenient storage bag for your car rides

- Stores books, newspapers and personal items
- CarPack conveniently attaches to headrest
- Wipes-clean
- Clear window for a name tag
- Pencil, pen holder closes with a zipper

ReOrder #	Description	Size	U/M
707038	Car Pack	14"W x 18"H	Ea.

RESIDENT-RELEASE CLOSURES


Skil-Care™ offers four options in resident-release closures to suit your resident's needs and capabilities. Choose from the following:



Hook & Loop



Hook & Loop with Red Finger Grip



Side-Release Buckle



Resident-Friendly Buckle

RESIDENT-RELEASE SOFT BELTS



- Promotes proper seating and positioning
- Provides wheelchair safety
 - Available with four closures
 - Universally sized
 - Completely launderable
 - Conveniently attaches to the wheelchair

ReOrder #	Description	Size	U/M
301250	Hook & Loop		Ea.
301255	Hook & Loop with Red Finger Grip	26"L x 5"W (Front Pad w/ Buckle)	Ea.
301270	Side-Release Buckle	42"L Loop Straps	Ea.
301275	Resident-Friendly Buckle (as shown)		Ea.

RESIDENT-RELEASE NYLON BELTS



- Promotes proper seating and positioning
- Completely size adjustable
 - Easily secures to wheelchair
 - Available with four closures

ReOrder #	Description	Size	U/M
701010	Side-Release Buckle	42"L x 2"W	Ea.
701011	Side-Release Buckle	54"L x 2"W	Ea.
701015	Resident-Friendly Buckle	42"L x 2"W	Ea.
701030	Hook & Loop (as shown)	37"L x 2"W	Ea.
701031	Hook & Loop (as shown)	60"L x 2"W	Ea.
701035	Hook & Loop with Red Finger Grip	37"L x 2"W	Ea.

RESIDENT-RELEASE SLIDER BELTS



Cross-over design prevents forward sliding

- Positions hips to prevent pelvic thrust
- 2" W webbing adjusts from 34" - 42"
- Easily attaches to the wheelchair

ReOrder #	Description	Size	U/M
610310	Hook & Loop	Adjusts from 35"L - 42"L	Ea.
610315	Hook & Loop with Red Finger Grip (as shown)	Adjusts from 35"L - 42"L	Ea.

RESIDENT-RELEASE SLIDE CONTROL BELTS



This non-restrictive slide control belt provides pelvic support

- Comes in two lengths for wheelchair and geri-chair
- Launderable
- Universal size

ReOrder #	Description	Size	U/M
610220	Slide Control Wheelchair Belt	Universal	Ea.
610230	Slide Control Geri-Chair Belt	Universal	Ea.

RESIDENT-RELEASE UNIVERSAL SAFETY BELTS



Multi-purpose belt for ambulation assistance


- 60" or 72" continuous loop belt for positioning assistance
- Wipe-clean nylon
- Convenient 2" side release buckle

ReOrder #	Description	Size	U/M
620010	Universal Safety Belt - Nylon	72" L, 2"H	Ea.
620020	Belt - Nylon	60" L, 2"H	Ea.

RESIDENT-RELEASE ECONO SEAT BELT

This nylon belts promotes proper seating, positioning and can be easily opened by the resident.

- Adjustable width
- Easily secures to wheelchair
- Press button closure
- Easy to release and is not considered to be a restraint when used by residents who are unable to release them.
- Wipe-clean fabric




ReOrder #	Description	Size	U/M
701055	Resident-Release Econo Seat Belt	50"L x 2"W	Ea.

BED POSITIONING



556010



556150




555055



BED POSITIONING

Residents who are bed-bound are at greater risk for developing pressure ulcers. Skil-Care understands the need to address this issue and offers a wide selection of Turning and Positioning Systems, Wedges, Pillows, Transfer System and Slide Sheets. Turning and Positioning Systems are designed to aid the process by making turning easier, safer and more secure. Wedges can be used to provide numerous areas of support, such as side-lying, trunk stabilization, hip-spine alignment, prevention of ulcers, and comfort.

30° POSITIONING WEDGES



30-Degree Wedge
w/Convoluted Surface

Provides optimum positioning for pressure relief

- Provides optimal 30-degree angle of incline for relieving pressure on hipbone and coccyx
- Wedge is comfortable, high-density convoluted foam
- Available in a Low-Shear II wipe-clean outer cover that reduces skin damaging friction
- Can be disinfected for infection control
- Non-slip bottom keeps wedge securely in place
- Available in Standard, Convoluted or Hi-Density Foams

ReOrder #	Description	Size	U/M
554022		16"L x 12"D x 7"H	Pair
554025	Standard Foam w/LSII Cover	24"L x 12"D x 7"H	Pair
554026		34"L x 12"D x 7"H	Pair
914451	Convoluted surface	24"L x 12"D x 7"H	Ea.
914452	w/LSII Cover	34"L x 12"D x 7"H	Ea.
554030	Hi-Density	24"L x 12"D x 7"H	Pair
554028	w/LSII Cover	34"L x 12"D x 7"H	Pair

45° POSITIONING WEDGES



Promotes proper hip-spine alignment

- Reduces sacral pressure to prevent decubitus ulcers
- Provides 45° angle of inclination for relieving pressure on hipbone and coccyx
- Available in a Low-Shear II cover with Anti-Slip Bottom (model 554010 only)

ReOrder #	Description	Size	U/M
554010	Wedge w/LSII Cover	17"L x 8"D x 8"H - 45°	Pair
554020	Wedge No Cover	17"L x 8"D x 8"H - 45°	Pair

IN-BED RESIDENT POSITIONING SYSTEM



Note: For more information on TLC Positioning Pad see pg. 108

TLC Positioning Pad

Stabilizes and positions bed residents


- Helps position resident in bed
- Reduces caregiver exposure to back injury
- Allows one caregiver to turn and position heavier residents
- Two wedged bolsters support the side-lying resident
- Hook & Loop strips on bolsters hold TLC positioning pad securely in place
- TLC Positioning pad's fluid-proof barrier protects linens and mattress
- Launderable pad also serves as transfer aid
- Bolsters are covered in durable, wipe-clean vinyl

ReOrder #	Description	Size	U/M
555036	Positioning System comes with One (1) 555015 TLC Positioning Pad & Two (2) 555035 34" Positioning Wedges	P - 48"L x 40"W W - 34"L x 8"D x 8"H	Set
555015	TLC Positioning Pad with Hook & Loop compatible bottom	48"L x 40"W	Ea.
555018		40"L x 36"W	Ea.
555025	45° Replacement Positioning Wedge	17"L x 8"D x 8"H	Ea.
555035		34"L x 8"D x 8"H	Ea.

ROLL CONTROL BOLSTERS



Original Skil-Care Design




Protects against falls by transforming any mattress to a raised edge mattress

- Lowers the risk of entrapment hazards of metal rails
- Provides secure positioning for side-lying residents
- Separation width between the bolsters are adjustable
Inside width expands to 29"W
Outside width expands to 44"W
- Can rotate resident for coccyx relief
- Attaches to the bed with straps
- Foam-padded for comfort
- Covered in durable, wipe-clean vinyl
- The optional TLC pad aids in resident positioning, see model 555012

ReOrder #	Description	Size	U/M
556010	Roll-Control Bolsters	34"L x 7"H x 8"D	Pair (as shown)
556012	Roll-Control Bolster	34"L x 7"H x 8"D	Single Sided (Ea.)
555012	TLC Positioning Pad	48"L x 40"W	Ea.

BED SUPPORT SYSTEM W/TWO ATTACHED 30° BOLSTERS



This bed bolster system has two (2) 30 degree wedges that are attached to the nylon sheet

- Launderable absorbent pad attaches to the system for added comfort and incontinence protection
- Absorbent pad has a non-skid bottom, can be used alone as a bed overlay and is launderable
- Sold with or without absorbent pad

ReOrder #	Description	Size	U/M
556020	Bed Support Bolster System w/Pad	32"L x 34"W x 6"H	Set (as shown)
556021	Bed Support Bolster System with Mesh* without Pad	32"L x 34"W x 6"H	Ea.
556025	Bed Support Bolster System without Pad	32"L x 34"W x 6"H	Ea.
556030	Optional Pad for Bed Support System	32"L x 20"W	Ea.

* For low air-loss beds

30° FULL BODY BED SUPPORT SYSTEM



This Full Body four bolster positioning system provides a stable platform for side lying residents when coccyx off-loading is needed.

- Covers full length of the bed (74"L x 34"W)
- 30 degree angle is the optimum angle of incline for relieving pressure on the hipbone and coccyx
- Comfortable, wipe clean, launderable nylon covering
- Easily attaches to the bed with six quick release straps
- Foam wedges are removable for laundering of cover

ReOrder #	Description	Size	U/M
556038	30° Full Body Bed Support System	74"L x 34"W	Set

30° BED SYSTEM WITH SLIDER SHEET AND TWO 16" OR 24" WEDGES



Provides coccyx pressure relief and a safe repositioning system for the caregiver

- Pressure relieving wedge off-loads the patients coccyx area eliminating the need for pillow positioning
- Slider sheet is waterproof and comes with continuous handles placed along the sides to enable the caregiver to conveniently grab at all positions
- Dual Wedges mate with sheet with high friction surface for optimum positioning

ReOrder #	Description	Size	U/M
556150	Bed Positioning Nylon Sheet w/Two 30 Degree 16" Wedges	16" Wedge, 50"L x 38"W slider sheet	Set
556151	16"L 30 Degree Wedges (2) for 556150 System	16"L x 12"D x 7"H	Pair
556152	Nylon Positioning Sheet for 556150 System	50"L x 38"W	Ea.
556153	24"L 30 Degree Wedges (2) for 556150 and 556154 System	24"L x 12"D x 7"H	Pair
556157	Bed Positioning Mesh Sheet w/Two 30 Degree 16" Wedges	16" Wedge, 50"L x 38"W slider sheet	Set
556158	Mesh Positioning Sheet for 556157 System	50"L x 38"W	Ea.
556154	Bed Positioning Nylon Sheet w/Two 30 Degree 16" Wedges	16" Wedge, 50"L x 48"W slider sheet	Set
556155	Nylon Positioning Sheet for 556154 System	50"L x 48"W	Ea.

UNIVERSAL GLIDE SHEET



Universal Glide Sheet provides a low friction surface

- Can be used with any Skil-Care transfer sheet or systems

ReOrder #	Description	Size	U/M
556065	Standard Mattress	72"L x 38"W	Ea.
556066	For Mattress 38"W - 44"W	72"L x 44"W	Ea.
556067	For Mattress 44"W - 56"W	72"L x 56"W	Ea.

BED GLIDE SYSTEM



The Bed Glide system is designed to stay under the patient for an extended period and is always ready to assist in boosting and repositioning.

- Reduces the risk of injury during boosting & repositioning
- Comes with disposable overlay sheet which is moisture absorbing and breathable.
- Allows resident to sit on it for an extended period
- Four pull handles for easy gripping
- Covers entire bed
- Glide pad is wipe-clean or launderable

ReOrder #	Description	Size	U/M
556087	Bed Glide System (Glide Pad & Overlay Sheet)	62"L x 30"W	Set
556088	Bed Glide Pad Only	62"L x 30"W	Ea.
556089	Disposable Cozy Cloth Overlay Sheet	53"L x 40"W	12/pk

ROLL SHEET



This Skil-Care sheet helps minimize the risk of back injury for the healthcare workers. The sheet makes it possible for nurses and staff to achieve compliance to your turning protocol while providing the best care and minimizing additional stress on the patient.

- Transfers and repositioning is the leading cause of caregiver injury, stop turning nurses into patients.
- Facilitates the transfer and repositioning of bed lying residents.
- Reduces time, caregiver exertion and back injury.
- Low-friction fabric helps to prevent shear and friction forces on resident's skin.
- Stands up to repeated laundering

ReOrder #	Description	Size	U/M
556068	Skil-Care Roll Sheet	46"L x 38"W	Ea.
556070		50"L x 42"W	Ea.

FOURHANDLE SLIDE SHEET



This Skil-Care nylon sheet helps minimize the risk of back injury for the healthcare workers. The sheet makes it possible for nurses and staff to achieve compliance to your turning protocol while providing the best care and minimizing additional stress on the patient.

- Transfers and repositioning is the leading cause of caregiver injury, stop turning nurses into patients.
- Facilitates the transfer and repositioning of bed lying residents with four handles on each corner.
- Reduces time, caregiver exertion and back injury.
- Low-friction single layer nylon fabric helps to prevent shear and friction forces on resident's skin.
- Stands up to repeated laundering
- 300lbs max. capacity

ReOrder #	Description	Size	U/M
556069	Nylon Slide Sheet	48"L x 40"W	Ea.

SLIDER SHEET WITH CONTINUOUS HANDLES




Double layer laminate with continuous handles

- Transfers and repositioning is the leading cause of caregiver injury, stop turning nurses into patients.
- Facilitates the transfer and repositioning of bed lying residents with continuous handles on each corner.
- Reduces time, caregiver exertion and back injury.
- Low-friction fabric helps to prevent shear and friction forces on resident's skin.
- Stands up to repeated laundering
- Heavy-duty reinforced laminate material

ReOrder #	Description	Size	U/M
556152	Nylon Positioning Sheet for 556150 System	50"L x 38"W	Ea.
556158	Mesh Positioning Sheet for 556157 System	50"L x 38"W	Ea.
556155	Nylon Positioning Sheet for 556154 System	50"L x 48"W	Ea.

LEG/THIGH ELEVATOR



Aligns the legs while simultaneously relaxing the spine

- Comfort foam off-loads the heels for prevention of pressure sores
- Aligns the leg and thighs
- Covered with a soft launderable Cozy Cloth fabric

ReOrder #	Description	Size	U/M
555040	Leg Elevator w/ Scoop Cutouts - Small/Medium	17"L x 7"H	Ea.
555041	Thigh Elevator w/ Scoop Cutouts - Medium/Large	20"L x 7"H	Ea.

KNEE ELEVATOR



Decreases leg and lower back pain

- Elevates the legs to relax the spine
- Covered in durable, wipe-clean anti-microbial vinyl
- Improves circulation
- Foam padded for comfort

ReOrder #	Description	Size	U/M
555043	Elevator Half Round Small	17"L x 8"D x 5"H	Ea.
555045	Elevator Half Round Medium	17"L x 8"D x 7"H	Ea.
555046	Elevator Half Round Bariatric	26"L x 8"D x 7"H	Ea.

PILLOW PROP



*As shown under the arm



Comfortable arm/head/leg support while lying in the bed

- Elevates head for added comfort
- Elevates and supports arms for residents with edema
- Effective for separating legs while sleeping
- Covered in fluid-proof, wipe-clean anti-bacterial nylon

ReOrder #	Description	Size	U/M
905010		13.5"L x 13.5"W	Ea.
905015	Pillow Prop	20"L x 15"W	Ea.

ELEVATING ARM SUPPORT




This soft foam elevator with recessed cavity creates a trough for resting the arm in a comfortable and stable position.

- Fits either the left or right arm
- Can be placed on any leveled surface
- Ideal for reducing swelling and edema
- Optional Cozy Cloth cover (905032)

ReOrder #	Description	Size	U/M
905030	Elevating Arm Support		Each
905031	Elevating Arm Support w/Cozy Cloth Cover	20"L x 6"W x 5.5"H	Each
905032	Optional Cozy Cloth Cover for elevating Arm Support	One Size	Each

BODY ALIGNER



- Economical foam body positioner
- Provides 45° angle of inclination for relieving pressure on hipbone and coccyx
 - Extended 8" height
 - Non-skid foam base that minimizes sliding

ReOrder #	Description	Size	U/M
911182	Body Aligner	20"L x 9"D x 8"H	Pair
911184			6/Cs

HIP/LEG ABDUCTION PILLOW



The Hip/Leg Abduction Pillow separates the legs of patients who are recovering from procedures or who are at risk of skin breakdown caused from friction or pressure.

- Polyfoam pillow with concave sides and two foam straps
- Product is intended to be used as a single-patient-use item
- Hook and loop fasteners on each side to comfortably secure the product in place between your legs

ReOrder #	Description	Size	U/M
710010	Hip/Leg Abduction Pillow	Small	Pair
710015	Hip/Leg Abduction Pillow	Medium/Large	Pair


ABDUCTOR/ CONTRACTURE CUSHION



- Prevents leg contracture and abduction in bed
- Designed to allow resident to sleep comfortably on back or sides
 - Secures behind knees with Hook & Loop straps
 - Low-Shear I outer cover is launderable and reduces skin damaging friction
 - The inner cover is made out of durable wipe clean nylon

ReOrder #	Description	Size	U/M
555060	Abductor/ Contracture Cushion	10" x 14" x 12"	Ea.

ADJUSTABLE ANTI-REFLUX BED WEDGE




2" Foam

- Doctors recommend upper-body elevation for GERD relief
- Specially designed wedge decreases night time reflux and heartburn
 - Allows user to sleep on back or side
 - Low-Shear II wipe-clean outer cover reduces skin damaging friction
 - Non-slip bottom keeps wedge in place
 - Height adjustable - 5" or 7", 8" or 10"

ReOrder #	Description	Size	U/M
555055	Anti-Reflux Bed Wedge w/LSII Cover Adjust to 5" or 7"H	24"L x 24"W	Ea.
555056	Anti-Reflux Bed Wedge w/LSII Cover Adjust to 8" or 10"H	24"L x 24"W	Ea.

BED-FOOT SUPPORT



- Provides comfortable foot positioning
- Helps prevent and controls foot drop
 - Keeps sheets and blankets off resident's feet
 - Securing straps keep the foot support in place
 - Foam padded and covered in wipe-clean nylon
 - Reversible for availability of 90° or 70° angle (as shown)

ReOrder #	Description	Size	U/M
554040	Bed-Foot Support	24"W x 13"H x 10"D	Ea.

ELEVATING BED WEDGES



Cozy Cloth Cover

Vinyl Cover



Comfortably elevates resident's head or feet

- Provides secure, effective positioning
- Available with a soft launderable Cozy Cloth fabric or wipe clean vinyl cover
- Choice of four heights (7.5", 10", 12" or 27")

ReOrder #	Description	Size	U/M
553010	Wedge w/Cozy Cloth Cover - 15 Degrees	24"W x 27"L x 7.5"H	Pair
553011	Wedge w/Vinyl Cover - 15 Degrees	24"W x 27"L x 7.5"H	Ea.
553030	Wedge w/Cozy Cloth Cover - 20 Degrees	24"W x 27"L x 10"H	Pair
553031	Wedge w/Vinyl Cover - 20 Degrees	24"W x 27"L x 10"H	Ea.
553050	Wedge w/Cozy Cloth Cover - 25 Degrees	24"W x 27"L x 12"H	Pair
553051	Wedge w/Vinyl Cover - 25 Degrees	24"W x 27"L x 12"H	Ea.
553061	Wedge w/Vinyl Cover - 45 Degrees	24"W x 27"L x 27"H	Ea.

ULTRA-SOFT BOLSTER



Original Skil-Care
comfort Design



Model# 554120 as shown above



Model# 554110 as shown above



Model# 554120 as shown above

Versatile bolsters provide positioning, pressure relief, and comfort

- Covered with a soft launderable Cozy Cloth fabric
- Inner cover wipes clean with disinfectant
- Elevates knees to reduce leg and lower back pain
- Provides leg/knee abduction to prevent pressure
- Filled with soft spun fiber
- Covered with two covers:
 - Inner-water proof wipe-clean
 - Outer-soft, launderable micro-fiber

ReOrder #	Description	Size	U/M
554110	Ultra-Soft Bolster	8" Diameter, 24" Long	Ea.
554120		6" Diameter, 18" Long	Ea.

BED RAIL PADS



Skil-Care developed an effective line of bedside safety products including bed rail pads, mattress extenders, gap guards, positioning wedges and bed barrier surfaces.

Bed rail pads serve the purpose of cushioning the bed rails and filling the gap between the rails and the mattress and outer rail. Skil-Care offers various types, lengths and materials. Custom sizes can be made per your specification.

VINYL BED RAIL PADS (CLASSIC)




15-inch-high pads provide maximum cushioning

- Pads extend below and between rail and mattress for added protection against entrapment
- Hook & Loop closures fasten pads securely to bed rails
- Covered in durable, reinforced, anti-microbial wipe-clean vinyl
- Available in 1" and 2" thickness's

ReOrder #	Description	Size	U/M
401010		60"L x 15"H x 1"D	Pair
401040		72"L x 15"H x 1"D	Pair
401050		80"L x 15"H x 1"D	Pair
401100	Vinyl Bed Rail Pads	60"L x 15"H x 2"D	Pair
401140		72"L x 15"H x 2"D	Pair
401150		80"L x 15"H x 2"D	Pair

BED RAIL WEDGE PAD



16-inch-high pads provide maximum cushioning

- Pads extend below and between rail and mattress for added protection
- Helps reduce entrapment in rails or mattress
- Hook & Loop closures fasten pads securely to bed rails
- Covered in durable, reinforced, anti-microbial wipe-clean vinyl
- Rises 10.5" from the bed surface

ReOrder #	Description	Size	U/M
401235	Bed Rail Wedge Pad - Full	70" Long	Pair
401240	Bed Rail Wedge Pad - Half	35" Long	Pair

SPLIT-RAIL VINYL BED RAIL PADS



Sized to fit four-section bed rails

- Cushioned with 1"-thick foam padding
- Prevents entrapment in rails
- Hook & Loop closures fasten pads securely to bed rails
- Covered in durable, reinforced, anti-microbial wipe-clean vinyl

ReOrder #	Description	Size	U/M
401080	Split-Rail Vinyl Bed Rail Pads	28"L x 9"H x 1"D	Set of 4

HALF RAIL VINYL BED PADS



Fits most half-rails, assist rails, and new anti-entrapment rails

- Protects against contact with metal rails
- Cushions top of rail for added protection
- Prevents entrapment in rails or mattress
- Hook & Loop closures secure pads easily to rail
- Covered in durable, reinforced, anti-microbial wipe-clean vinyl

ReOrder #	Description	Size	U/M
401090	Half-Size Vinyl Bed Rail Pads	36"L x 14"H x 1"D	Pair

CUSHION TOP VINYL BED RAIL PADS




Top and side cushioning provides added protection

- Pads extend between and below mattress and bed rails
- Over the top foam padded for comfort and extends 4" over the rail
- Covered in durable, reinforced, anti-microbial wipe-clean vinyl
- Available in three lengths, 1" thickness

ReOrder #	Description	Size	U/M
401210		60"L x 16"H x 4" fold over	Pair
401220	Cushion Top Vinyl Bed Rail Pads	72"L x 16"H x 4" fold over	Pair
401230		80"L x 16"H x 4" fold over	Pair

GAP GUARD



Gap Guard is 2 3/4" wide at top and tapers to 2" at bottom, allowing for a tight fit.



Closes space between mattress and headboard, footboard or side railst

- Prevents arm and leg entrapment
- Wedge shape fills gaps up to three inches
- Works when side rails are raised or lowered
- Durably constructed of polyfoam
- Covered in incontinence-proof, wipe-clean nylon
- Available for full and half length bed rails

ReOrder #	Description	Size	U/M
401310	Gap Guard - Full Length Side Rail	70" L, 2 3/4" max thickness	Ea.
401320	Gap Guard - Half Length Side Rail, for use with Headboard & Footboard	35" L, 2 3/4" max thickness	Pair

THRU-VIEW VINYL BED RAIL PADS



Standard Window
(13.5"W x 6.5"H)



Large Window (35"W x 5"H)



Choose standard or large clear vinyl window

- Reduces resident's feelings of isolation
- Makes monitoring resident easier
- Protects against bed rail related injuries
- 15"-high pads fit well below rail and mattress for added protection
- Straps with side-release buckles fasten pad securely to rails
- Prevents entrapment in rails
- Covered in durable, reinforced, anti-microbial wipe-clean vinyl

ReOrder #	Description	Size	U/M
401015	Standard Window	60" Long	Pair
401017	Large Window	60" Long	Pair
401045	Standard Window	72" Long	Pair
401047	Large Window	72" Long	Pair
912220	Full View - Half Rail	29" Long	Set of 4

SYNTHETIC SHEEPSKIN BED RAIL PADS



Soft, high-pile pads provide full-rail protection

- Covers resident's side and top of bed rail
- Prevents entrapment in rails
- Health grade sheepskin is non-allergenic
- Completely launderable
- Attached with easy release snaps



ReOrder #	Description	Size	U/M
402010	Synthetic Sheepskin Bed Rail Pads	60"L x 18"W	Pair

BED/WALL PROTECTOR



Protects the wall from costly repairs resulting from the bed bumping into the wall

- Easy to clean nylon cover
- Non-slip backing
- Fits all beds from 36"-42"W

ReOrder #	Description	Size	U/M
401165	Bed/Wall Protector	38"L x 6"H x 3"D	Ea.

CUSHION FOR FOOT AND HEADBOARD




Soft vinyl covered foam pads that protect the residents head and/or feet from banging into the foot or headboard


- Headboard cushion is a full 2" x 13" x 36" that protects the head from banging into the headboard
- Footboard cushion is a full 2" x 9" x 36" that protects the feet from banging into the footboard.
- Wipe clean vinyl
- Elastic straps that wrap around the head or footboard.
- Can be placed between the mattress and the frame to "fill the gap"
- Can be universally fitted to all medical beds

ReOrder #	Description	Size	U/M
401160	Headboard Cushion	36"W x 2"D x 13"H	Ea.
401170	Footboard Cushion	36"W x 2"D x 9"H	Ea.


MATTRESSES



558015



558051





558123




All Skil-Care Mattresses come with our support pod technology to optimally address the key issues of weight distribution, pressure relief, anti-shearing and provide a positive thermo-climate by eliminating heat build-up. Our Premium Pressure-Check mattress comes with a 4-way stretch fluid-proof wipe clean cover and a 7 year non-prorated warranty. Custom Bariatric sizes can be made per your specifications. The Psychiatric mattress comes in nylon or vinyl and are stitched closed without a zipper.

Premium Pressure-Check MATTRESSES



The Premium Pressure-Check 2-Layer, 3-zone, 6" Mattress is specifically designed to provide an industry leading pressure redistribution and pressure relief.

- 400 individual "Support Pods" provide weight distribution, pressure relief, anti-shearing and a positive thermo-climate by eliminating heat buildup.
- The top layer of the mattress is made out of a unique polyurethane foam material.
- The bottom layer of the mattress is high-density compression resistant foam to support weight distribution and extend the mattress life.
- 3-zones provide specific pressure redistribution between head, heel and trunk.
- The top portion of the cover is made out of Skil-Care's Low-Shear III 4-Way stretch super-durable, breathable, fluid-proof fabric that wipes clean for easy care while protecting sensitive skin against friction and shear damage.
- The bottom section of the cover made out of slip-resistant fabric.
- All of the mattresses come with a Seven (7) year non-prorated replacement warranty and comply with fire code 16 CFR 1633.
- Available with or without a 2" raised perimeter guard with an EZ cutout section.




Perimeter-Guard:
2" durable foam wall to prevent "roll-off"

ReOrder #	Description	ReOrder #	Description	Size	U/M
558015		558051		76"L x 35"W x 6"H	Ea.
558020	Pressure-Check Mattress	558052	Perimeter-Guard Mattress	80"L x 35"W x 6"H	Ea.
558025		558053		84"L x 35"W x 6"H	Ea.

PREMIUM PRESSURE-CHECK GEL-INFUSED VISCO MATTRESS

For residents who are looking for additional comfort the 3-Layer Gel-Infused visco 6" mattress has all the same attributes (properties) of the Premium Pressure-Check Mattress with the added benefit of a third layer of Gel-Infused Visco Foam.

- The Gel-Infused Visco combines the properties of softness and resiliency with gel infusion to keep mattress cool, fresh and dry.
- The Gel-Infused Visco foam keeps the weight off of the pressure points by cradling to the body.
- The gel in the foam helps maintain your body temperature by allowing the body heat to escape.
- Available with or without a 2" raised perimeter guard with EZ cutout section.



ReOrder #	Description	ReOrder #	Description	Size	U/M
558123		558133		76"L x 35"W x 6"H	Ea.
558125	Gel-Infused Visco Foam Mattress	558135	Perimeter-Guard Gel-Infused Visco Mattress	80"L x 35"W x 6"H	Ea.
558127		558137		84"L x 35"W x 6"H	Ea.


- All of the mattresses come with a Seven (7) year non-prorated replacement warranty and comply with fire code 16 CFR 1633.

PREMIUM PRESSURE-CHECK BARIATRIC MATRESSES

The Premium Pressure-Check Bariatric 2-Layer, 3-zone, 6" Mattress is made the same as the standard size Premium Pressure Check mattress with minor differences.

- The top portion of the cover is made out of Skil-Care's Low-Shear II super durable, breathable, fluid-proof stretch fabric that wipes clean for easy care while protecting sensitive skin against friction and shear damage.
- The mattresses come with a Five (5) year non-prorated replacement warranty and comply with fire code 16 CFR 1633.
- Available with Nylon or Vinyl cover

ReOrder #	Description	Size	U/M
558100		80"L x 42"W x 6"H	Ea.
558102	Bariatric Mattress	80"L x 48"W x 6"H	Ea.
558104		80"L x 54"W x 6"H	Ea.




* For other custom sizes call a Skil-Care customer service representative

PREMIUM PRESSURE-CHECK PSYCHIATRIC MATRESSES

The Premium Pressure-Check Psychiatric 2-Layer, 3-zone, 6" Mattress is made the same as the standard size Premium Pressure Check mattress with minor differences.

- The mattress cover is completely enclosed with a super strength stitching (No Zipper)
- The high strength, rip resistant, waterproof cover is easy to clean and comes with a non-slip bottom
- The mattresses come with a Five (5) year non-prorated replacement warranty and comply with fire code 16 CFR 1633.
- Available with or without a 2" raised perimeter guard with EZ cutout section




No Zipper
No Metal Parts

ReOrder #	Description	ReOrder #	Description	Size	U/M
912308		912324	Perimeter-Guard	76"L x 35"W x 6"H	Ea.
912310	Psych Mattress (Nylon Cover)	912328	Psych Mattress	80"L x 35"W x 6"H	Ea.
912314		912330	(Nylon Cover)	84"L x 35"W x 6"H	Ea.
912311	Psych Mattress (Vinyl Cover)	—	—	80"L x 35"W x 6"H	Ea.

- For Custom Sizes contract Customer Service or your Sale Representative.
- **No Zipper on cover.**
- Available with Vinyl cover (Model 912311).

REPLACEMENT COVERS



Waterproof comfortable stretch covers

- The Low-Shear III cover is a 4-Way Stretch moisture resistant, super durable, breathable, fluid-proof fabric that wipes clean for easy care while protecting sensitive skin against friction and shear damage.
- Heavy duty zipper with waterfall covering

ReOrder #	Description	Size	U/M
558091		72"L x 36"W x 6"H	Ea.
558093	Mattress Replacement Cover - 4-Way Stretch Low-Shear III Cover	76"L x 36"W x 6"H	Ea.
558092		80"L x 36"W x 6"H	Ea.
558094		84"L x 36"W x 6"H	Ea.
558089	Mattress Replacement Cover - Perimeter Guard 4-Way Low-Shear III Cover	80"L x 36"W x 6"H	Ea.

SYNTHETIC SHEEPSKIN PADS



Model 501021



Model 501010



Reduces pressure and friction

- Aids in the reduction of pressure sores
- Allows for air circulation beneath resident
- Remains soft and resilient even after repeated laundering

ReOrder #	Description	Size	U/M
501011		24"L x 30"W	Ea
501010	Sheepskin Pad - White	24"L x 30"W	Dz
501021		24"L x 30"W	Ea
501051		40"L x 30"W	Ea.
501090		60"L x 30"W	Ea.
501020	Sheepskin Pad - Beige	24"L x 30"W	Dz
501050		40"L x 30"W	Dz
501022		24"L x 30"W	6/pk
501210		30 inches x 30 feet	Roll

SAFETY STRAP



Ideal stabilizer for foam mattresses

- Conveniently affixes the foam mattress to the bed
- Easy side release buckles for convenient installation and release
- Supplied with two units, one each for headboard and footboard region
- Launderable
- Custom size available upon request



ReOrder #	Description	Size	U/M
558070	Mattress Safety Strap Fits 35"- 42"W Mattress	36"L x 8"W	Pair
558071	Mattress Safety Strap Fits 43"- 48"W Mattress	42"L x 8"W	Pair
558072	Mattress Safety Strap Fits 49"- 60"W Mattress	48"L x 8"W	Pair

GERI CHAIR



Geri Chairs are designed to provide a more substantial seating platform. Skil-Care understands comfortable seating is important and offers a selection of product specifically designed to complement the Geri Chair.

ADJUSTABLE LATERAL SUPPORT



original
Skil-Care
CORPORATION
Design



Easy adjustment to accommodate a range of body widths.

Lateral support is provided by width-adjustable side panels

- Minimizes lateral leaning
- Foam-padded sides adjust laterally for a custom fit

ReOrder #	Description	Size	U/M
706055	Adjustable Lateral Support - Medium	Adjusts from 12"W - 18"W	Ea.
706059	Adjustable Lateral Support Large/Bariatric	Adjusts from 16"W - 24"W	Ea.

UNIVERSAL CHAIR SIDE WINGS



original
Skil-Care
CORPORATION
Design



Provides additional comfort for lateral leaning

- Improves sitting posture when chair is in upright or reclined position
- Prevents sideways movement
- Foam padded for comfort
- Low-Shear II cover wipes clean for easy care
- Secures to all geri-chairs, wheelchairs and lounge chairs

ReOrder #	Description	Size	U/M
706319	Universal Chair Side Wings	18.5"W x 10"D x 10"H	Ea.

GERI CHAIR SMOKERS APRON



original
Skil-Care
CORPORATION
Design



906015 White



906016 Blue

Protects resident and geri-chair against cigarette burns

- Flame proof
- Full-sized to cover lap and chest
- Lightweight for resident comfort
- Resistant to cracking

ReOrder #	Description	Size	U/M
906015	White	39"W x 44"L	Ea.
906016	Blue		Ea.

GERI CHAIR ALARM SYSTEM

This magnetic Geri-Chair Alarm system comes with an alarm, bracket, and belt length accommodates most geri-chairs.

- Adjustable strap for easy installation
- Keeps alarm out of resident's reach

original
Skil-Care
CORPORATION
Design



ReOrder #	Description	Size	U/M
909212	Geri-Chair Alarm System	Universal	Ea.
909211	Personal Alarm	Universal	Ea.

See page 97 for more details



GERI CHAIR COZY SEAT



Original Skil-Care Corporation Design



Protects against skin tears and relieves pressure

- Ample foam padding improves geri-chair comfort
- Reversible with Cozy Cloth on both sides
- Covers metal chair components to reduce cuts and bruises
- Designed not to slip or bunch-up
- Two models: one with and one without leg-rest protection
- Covered with a soft launderable Cozy Cloth fabric

ReOrder #	Description	Size	U/M
703002	Backrest and Seat	52"L x 21"W	Ea.
703003	Backrest, Seat and Leg-rest (as shown)	70"L x 21"W	Ea.

GERI CHAIR GEL OVERLAY



Pressure relief for entire body

- Water-based gel cushion in the seating area relieves interface pressure at coccyx and ischials (20"W x 18"D)
- Foam padding for back and legs offers cushioned comfort
- Air-Flow Mesh fabric cover maintains comfort
- Reverse side is fluid-proof nylon
- Cover is launderable after removing the cushions
- Easily secures to all geri-chairs

ReOrder #	Description	Size	U/M
703001	Geri-Chair Gel Overlay	68"L x 19"W x 2"H	Ea.

GERI CHAIR FOOT CRADLE



Original Skil-Care Corporation Design



Positions and protects resident's feet

- Prevents feet from slipping off footrest
- Controls foot-drop
- Foam padded for comfort
- Vinyl cover wipes clean for easy care
- Easily secures to all geri-chairs

ReOrder #	Description	Size	U/M
703430	Geri-Chair Foot Cradle	14"W x 14"D x 18"H	Ea.

GERI CHAIR LEG POSITIONER



Original Skil-Care Corporation Design




Prevents residents legs from slipping off the geri-chair leg rest

- Comfortably positions legs on the leg rest
- Prevents legs from slipping off sides
- Assists in alignment of the hips
- Easily attaches to all geri chairs with a quick release strap
- Low-Shear 1 launderable outer cover that reduces skin damaging friction

ReOrder #	Description	Size	U/M
703420	Geri-Chair Leg Positioner	19"L x 5.25"H x 6.5"D	Ea.



Pro Series Alarms



911572

SAFETY



SkiL-Care's Fall Prevention line of safety alarms was developed to help your facility reduce the risk of fall-related injuries. Our alarms are designed to provide **Functionality, Flexibility and Durability** and come with a variety of features to alert caregivers when a patient attempts to arise unattended. The Skil-Care's Alarms are easy to use and can function in many situations including bed, chair, floor, door & window.

Functionality was achieved through effective design, the use of Auto-reset functions, easy mounting and attachment

and a choice of alarm tones, including a new "music" alarm tone, which supports noise reduction efforts in the facility.

Flexibility was achieved by designing the Pro Series of Alarms which includes the ChairPro, MultiPro and BedPro. All of the pressure sensors shown on the following pages can be utilized on any of the "Pro" alarm units. All Pro Series alarms provide an audio alarm at the electronics box and in addition the MultiPro and BedPro alarms have a choice of music or alarm tone, and nurses call capabilities. The BedPro alarm also provides the option to silence the audio alarm at the resident site while utilizing the nurse call system feature, and also comes with an optional a/c adapter.

Durability is achieved by embedding the solid state components in heavy duty circuit boards and also combining a heavy duty outside housing with rubber bumpers. All electronics units are guaranteed for one year. Sensor Pads are guaranteed from 45 days up to one year as indicated on the enclosed charts. ChairPro and BedPro alarms are sold as sets which include: sensor, electronics, batteries and accessories. The MultiPro items are sold as separate components.

Fall Mats

To help facilities reduce fall risk, Skil-Care developed an effective line of bedside safety products including bed rail pads, mattress extenders, positioning wedges and bed barrier surfaces. When these systems are not sufficient, Skil-care also provides a wide selection of popular floor mats that reduce the impact from fall-related injuries. Our mats come in many different sizes and styles. See page 102 for more details.

- Folding/Non-Folding
- Beveled/Non-Beveled
- Vinyl covered/Molded foam
- Various heights (1", 1.5", 2", 3.5" or 4")

CHAIRPRO SAFETY ALARM UNITS



The alarm unit provides sounds that alert the caregiver when joined with a sensor (provided on other pages).

The ChairPro alarm includes the features below:

- Auto-Reset
- High-Low Volume
- Low Battery Indicator
- Drop protection bumpers
- Green/Red Flashing Indicator Light
- Battery Test
- Buzzer Tone Only
- Cable Strain Relief
- Mountable to chair or bed

ReOrder #	Description	U/M
909367	ChairPro Alarm Unit w/ Accessories	Ea.
909366		10/pk

SENSOR PAD ALARM SETS

Monitors egress of seated residents

- Alarm activates if resident leaves/falls off the chair
- Automatically resets when resident sits down



7" x 15" Sensor Pad



16" x 16" Sensor Pad

Alarm Sets (Alarm Unit, Sensor Pad & Accessories)

ReOrder #	Description	Pad Size	U/M
909350	ChairPro Sensor Pad Alarm Set - Guarantee 45 Days	7" x 15"	Set
909355	ChairPro Sensor Pad Alarm Set - Guarantee 180 Days	7" x 15"	Set
909356	ChairPro Sensor Pad Alarm Set - Guarantee 1 Year	7" x 15"	Set
909430	ChairPro Sensor Pad Alarm Set - Guarantee 45 Days	16" x 16"	Set
909434	ChairPro Sensor Pad Alarm Set - Guarantee 180 Days	16" x 16"	Set

REPLACEMENT PADS

ReOrder #	Description	Pad Size	Guarantee	U/M
909360		7" x 15"	45 Days	Ea.
909362		7" x 15"	45 Days	10/pk
909363		7" x 15"	180 Days	Ea.
909365		7" x 15"	180 Days	10/pk
909431	Sensor Pad Only	16" x 16"	45 Days	Ea.
909432		16" x 16"	45 Days	10/pk
909435		16" x 16"	180 Days	Ea.
909436		16" x 16"	180 Days	10/pk

All ChairPro™ Alarm Sets Include Alarm Unit and a Unique Sensor




BI-FOLD SENSOR PAD FOR CHAIR & BED



Skil-Care developed a dual purpose waterproof sensor pad that can be used on the bed and/or the wheelchair.

- Sensor pad can be placed on the wheelchair (folded) to fit any size wheelchair greater than 15" x 15"
- Sensor pad can be used on the bed (unfolded - 15" x 30")
- Sensor pad can be used under or over the mattress/seat cushion
- Resistant to body fluids
- Coiled cord for easy connection
- Compatible with any Skil-Care Pro Series alarm units.



ReOrder #	Description (Sets)	Pad Size	Guarantee	U/M
909420	ChairPro Alarm Unit w/Bi-Fold Sensor Pad	15"W x 15"L (Folded)	1 Year	Set
909421	BedPro Alarm Unit w/Bi-Fold Sensor Pad	15"W x 30"L (Flat)	1 Year	Set

ReOrder #	Description (Pad Only)	Pad Size	Guarantee	U/M
909425		15"W x 15"L (Folded)	1 Year	Ea.
909426	Bi-Fold Sensor Pad for Chair or Bed	15"W x 30"L (Flat)	1 Year	5/pk
909427		15"W x 30"L (Flat)	1 Year	10/pk
909427P*			1 Year	10/pk

*All products with "P" suffix have 11.5" cord length.


BELT ALARM SETS



- Automatically resets when belt is secured
- Alarm sounds when closure is released
- Reinforced attachment holes and D-ring option for easy application
- Reminds resident to stay seated
- Alerts caregiver before resident leaves the chair
- Provides restraint-free wheelchair safety
- Universally sized belt secures to all wheelchairs

EZ RELEASE BELT ALARM SETS

Hook & Loop closure with red loop for easy release



Magnetic Closure



909372 - EZ Release Belt w/Adjustable Side Release Loop Attachment

ReOrder #	Description	Size	U/M
909371	ChairPro EZ Release Belt Alarm Set with Grommets	50"L	Set
909395		50"L	10/pk
909372	ChairPro EZ Release Belt Alarm Set w/Adjustable Loop Attachment	45"L	Set

SEAT BELT ALARM SETS

Nylon belt has metal press release buckle for restraint-free care



909370 - Metal Press Release Seat Belt with Grommets



909385 Seat Belt w/Adjustable Loop Attachment

ReOrder #	Description	Size	U/M
909370	Metal Press Release Seat Belt Alarm Set with Grommets	50"L	Set
909373		70"L	Set
909385	Metal Press Release Seat Belt Alarm Set with Adjustable Loop Attachment*	45"L	Set
909386		65"L	Set

***Note that when using the loop attachment for securing to the wheelchair no hardware is required.**

UNDERSEAT ALARM SET



Effective alarm not accessible to resident

- Under-seat sensor is protected against incontinence and wear
- Alarm immediately activates when resident rises
- Automatically resets when resident sits
- Sensor stays in place to reduce false alarms
- Equally effective with or without seat cushion
- Under-seat placement deters tampering
- Easily removed when power washing the chair



ReOrder #	Description	Size	U/M
909369	ChairPro UnderSeat Alarm Set	Fits 18"-22"W Wheelchairs	Set


All ChairPro™ Alarm Sets Include Alarm Unit and a Unique Sensor



GEL-FOAM CUSHION ALARM SET

Gel cushion provides pressure relief, comfort, stability & an alarm alert

- Alarm initiates when user rises and resets when he/she sits
- Fluid-proof, heat-sealed vinyl cover protects against incontinence
- Low-Shear I launderable outer cover reduces skin shear
- Alarm sensor is protected in vinyl pocket beneath cushion



ReOrder #	Description	Size	U/M
909380		16"W x 16"D x 2.5"H	Set
909382	ChairPro Gel-Foam Pad Alarm Set	18"W x 16"D x 2.5"H	Set
909384		20"W x 16"D x 2.5"H	Set
909389	Replacement Sensor Unit	-	Ea.

LEVELING PAD ALARM SET

Combines an alarm with a firm foundation for all cushions

- Alarm immediately activates when resident rises
- Automatically resets when resident sits
- Sensor is protected in a vinyl pocket beneath the leveling pad
- Vinyl-covered leveling pad eliminates sling-seat hammocking
- Wheelchair cushion required, sold separately



Eliminates Sling-Seat Hammocking



ReOrder #	Description	Size	U/M
909375	ChairPro Leveling Pad Alarm Set	Fits 16" - 18" Wheelchairs	Set
909535	Leveling Pad Sensor (Pad & Sensor Unit)	Fits 16"-18" Wheelchairs	Ea.
909389	Replacement Sensor Unit	-	Ea.

LIFT-OFF LAP CUSHION ALARM SET



Positioning plus a chair-exit alarm

- Automatically resets when Hook & Loop is secured
- Easy-open Hook & Loop closure has red loop to facilitate release by hand. Releases easily as resident rises
- Alarm sounds if closure is opened or resident stands
- Lap cushion is foam padded for comfort
- Covered in durable wipe-clean nylon
- Secures to all standard wheelchairs



All ChairPro™ Alarm Sets
Include Alarm Unit and a
Unique Sensor

ReOrder #	Description	Size	U/M
909374	ChairPro Lift-Off Lap Cushion Alarm Set	Fits 16"-18"W Wheelchairs	Set

GEL-FOAM TOILET SEAT CUSHION ALARM PAD



Layered with gel and foam to provide comfort and maximum pressure relief and pressure redistribution

- Built-In sensor pad
- Alarm alerts caregiver when the resident stands up
- Available with either ChairPro or MultiPro alarms
- Improves comfort on a toilet seat
- Prevents bruising and provides relief for residents with sores
- Convenient optional mounting bracket



ChairPro™
Alarm unit included with
model# 909376 ONLY



MultiPro™
Alarm unit included with
model# 909377 ONLY

Wall Mounting Bracket
Model 909398



909397

ReOrder #	Description	Size	U/M
909376	Gel-Foam Toilet Seat Cushion Alarm Set (Includes: ChairPro Alarm , Toilet Seat Cushion)		Set
909377	Gel-Foam Toilet Seat Cushion Alarm Set (Includes: MultiPro Alarm , Toilet Seat Cushion)	15"W x 18.5"D	Set
909378	Gel-Foam Toilet Seat Replacement Cushion w/ Built-In Sensor		Ea.
909397	Replacement Straps for Gel-Foam Toilet Seat Cushion	Universal	24 pieces
909367	ChairPro Alarm Unit with Accessories	—	Ea.
909510	MultiPro Alarm Unit with Accessories	—	Ea.
909398	Universal Alarm Wall Bracket	7"L x 2.5"W	Ea.

BEDPRO SAFETY ALARM UNITS



Protects residents against fall-related injuries resulting from unassisted bed exits.

The BedPro™ alarm includes the features below:

- Auto-Reset
- In-Room On-Off
- Battery Test
- Battery Included
- High-Low Volume
- Cable Strain Relief
- Mountable
- Nurse Call
- Low Battery Indicator
- AC Adaptor
- Two Buzzer/Two Music Tones
- Bumpers
- Indicator Light*

*Indicator light flashes GREEN when monitoring, RED when alarming

ReOrder #	Description	U/M
909335	BedPro Unit with accessories (Accessories: Nurse Call Cord, Y-Connector, AC Adapter & Batteries)	Ea.
909334	BedPro Unit without accessories	Ea.
909336	BedPro Unit with accessories	10/pk

Note: BedPro alarm can also be matched with any of the sensors shown in the ChairPro or MultiPro sections.

The BedPro™ alarm is available with Skil-Care's UnderMattress, OverMattress & Early Warning Sets.

UNDERMATTRESS SETS & SENSOR PAD

- Sensor is placed under mattress and protected against incontinence
- Weight of mattress keeps sensor pad in place to reduce false alarms
- No sensor pad on mattress to compromise pressure relief or comfort
- Universally connects to any bed, chair or multipro alarms
- Can be also used above the mattress



Alarm Sets (Alarm, Sensor Pad & Accessories)

ReOrder #	Description	Pad Size	Guarantee	U/M
909320		10" x 30"	90 Days	Set
909322		10" x 30"	180 Days	Set
909340	BedPro UnderMattress Alarm Set	10" x 30"	1 Year	Set
909321	w/Sensor Pad	20" x 30"	90 Days	Set
909323		20" x 30"	180 Days	Set
909341		20" x 30"	1 Year	Set



All BedPro™ Alarm Sets Include Alarm Unit and a Unique Sensor

Sensor Pad Only

ReOrder #	Pad Size	ReOrder #	Pad Size	Description	Guarantee	U/M
909325	10" x 30"	909328	20" x 30"	Sensor Pad Only	90 Days	Ea.
909327	10" x 30"	909330	20" x 30"		90 Days	10/pk
909331	10" x 30"	909337	20" x 30"		180 Days	Ea.
909333	10" x 30"	909339	20" x 30"		180 Days	10/pk
909401	10" x 30"	909403	20" x 30"		1 Year	Ea.
909401P*	10" x 30"	909404	20" x 30"		1 Year	10/pk
909402	10" x 30"					
909402P*	10" x 30"					

*All products with "P" suffix have 11.5" cord length.

OVERMATTRESS ALARM SETS & SENSOR PAD



All BedPro™ Alarm Sets Include Alarm Unit and a Unique Sensor

Alerts caregivers of unassisted bed exits

- Flexible wire connection increases reliability
- Sets comes with popular 10" x 30" sensor pad




Alarm Sets (Alarm Unit, Sensor Pad & Accessories)

ReOrder #	Description	U/M
909310	BedPro OverMattress Alarm System w/45 Day Sensor	Set
909311	BedPro OverMattress Alarm Set w/180 Day Sensor	Set

OVERMATTRESS SENSOR PADS

ReOrder #	Description	Pad Size	Guarantee	U/M
909314			45 Days	Ea.
909316	Sensor Pad Only	10" x 30"	45 Days	10/pk
909317			180 Days	Ea.
909319			180 Days	10/pk

MULTIPRO SAFETY ALARM UNIT



Simply connect the MultiPro alarm to any of the sensors shown in the alarm section to monitor all your fall-risk residents.

The MultiPro™ alarm includes the highlighted features below:

- Auto-Reset
- Battery Test
- Battery Included
- Cable Strain Relief
- Mountable
- Indicator Light*
- Nurse Call
- Low Battery Indicator
- High-Low Volume
- Bumpers
- Two Buzzer/Two Music Tones

*Indicator light flashes GREEN when monitoring, RED when alarming

ReOrder #	Description	U/M
909510	MultiPro Alarm Unit with Accessories (Accessories: Auxiliary Connector Cable, Battery and Hook & Loop)	Ea.
909515	Auxiliary Connector Cable, Battery and Hook & Loop	10/pk

The MultiPro™ alarm is sold as a unit and works with any of the safety sensors shown in this alarm section, which are sold separately

SENSOR PADS

Monitors movements of seated residents

- Alarm activates if resident leaves or falls from the chair
- Compatible with any ProSeries alarms
- Durable heat sealed construction
- Flexible wire connection increases reliability

ReOrder #	Description	Pad Size	Guarantee	U/M
909360			45 Days	Ea.
909362	Sensor Pad Only	7" x 15"	45 Days	10/pk
909363			180 Days	Ea.
909365			180 Days	10/pk
909431			45 Days	Ea.
909432	Sensor Pad Only	16" x 16"	45 Days	10/pk
909435			180 Days	Ea.
909436			180 Days	10/pk



16" x 16" Sensor Pad

UNDERMATTRESS SENSOR PADS

- Sensor is placed under mattress and protected against incontinence
- Weight of mattress keeps sensor pad in place to reduce false alarms
- No sensor pad on mattress to compromise pressure relief or comfort
- Universally connects to any bed, chair or multipro alarms
- Flexible wire connection increases reliability
- Can also be used above the mattress

***Note: This pad requires either Chair, Bed or MultiPro Alarm Unit**

Original
Skil-Care
CORPORATED
Design



Sensor Pad Only

ReOrder #	Pad Size	ReOrder #	Pad Size	Description	Guarantee	U/M
909325	10" x 30"	909328	20" x 30"	Sensor Pad Only	90 Days	Ea.
909327	10" x 30"	909330	20" x 30"		90 Days	10/pk
909331	10" x 30"	909337	20" x 30"		180 Days	Ea.
909333	10" x 30"	909339	20" x 30"		180 Days	10/pk
909401	10" x 30"	909403	20" x 30"		1 Year	Ea.
909401P*	10" x 30"	-	-		1 Year	Ea.
909402	10" x 30"	909404	20" x 30"		1 Year	10/pk
909402P*	10" x 30"	-	-		1 Year	10/pk


*All products with "P" suffix have 11.5" cord length.

OVERMATTRESS SENSOR PADS

Sensor Pads Only

ReOrder #	Pad Size	Guarantee	U/M
909314	10" x 30"	45 Days	Ea.
909316	10" x 30"		10/pk
909317	10" x 30"	180 Days	Ea.
909319	10" x 30"		10/pk

***Note: This pad requires either Chair, Bed or MultiPro alarm unit**



***Note: All the sensors in this section require either Chair, Bed or MultiPro alarm unit**

BI-FOLD SENSOR PAD FOR CHAIR & BED

Skil-Care developed a dual purpose no opening, waterproof sensor pad that can be used on the bed and/or the wheelchair.

- Sensor pad can be placed on the wheelchair (folded) to fit any size wheelchair greater than 15" x 15"
- Sensor pad can be used on the bed (unfolded - 15" x 30")
- Sensor pad can be used under or over the mattress/seat cushion
- Resistant to body fluids
- Coiled cord for easy connection
- Compatible with any Skil-Care Pro Series alarm units.



Original
Skil-Care
CORPORATED
Design

ReOrder #	Description (Pad Only)	Pad Size	Guarantee	U/M
909425		15"W x 15"L (Folded)	1 Year	Ea.
909426	Bi-Fold Sensor Pad for Chair or Bed	15"W x 30"L (Flat)	1 Year	5/pk
909427		15"W x 30"L (Flat)	1 Year	10/pk
909427P*			1 Year	10/pk

*All products with "P" suffix have 11.5" cord length.

UNDERSEAT ALARM SENSOR



Effective alarm sensor isn't accessible to resident

- Under-seat sensor is protected against incontinence and wear
- Sensor activates alarm when resident rises
- Automatically resets when resident sits
- Sensor stays in place to reduce false alarms
- Equally effective with or without seat cushion
- Under-seat placement deters tampering
- Easily removed when power washing the chair



ReOrder #	Description	Size	U/M
909391	UnderSeat Alarm Sensor	Fits 18"-22"W WC	Ea.

LEVELING PAD SENSOR



- Sensor activates alarm when resident rises
- Vinyl-covered leveling pad eliminates sling-seat hammocking
- Sensor is protected in vinyl pocket beneath leveling pad
- Wheelchair cushion required; sold separately
- Pad has built-in sensor



Eliminates Sling-Seat Hammocking



ReOrder #	Description	Pad Size	U/M
909535	Leveling Pad Sensor (Pad & Sensor Unit)	Fits 16"-18" Wheelchairs	Ea.
909389	Replacement Sensor Unit	-	Ea.

GEL-FOAM CUSHION AND SENSOR

- Alarm initiates when user rises and resets when he/she sits
- Cushion provides pressure relief and redistribution
- Fluid-proof, heat-sealed vinyl cover protects against incontinence
- Low-Shear I launderable outer cover reduces skin shear
- Alarm sensor is protected in vinyl pocket beneath cushion



Original Skil-Care CORPORATION Design

ReOrder #	Description	Size	U/M
909530		16"W x 16"D x 2"H	Ea.
909531	Gel-Foam Cushion and Sensor	18"W x 16"D x 2"H	Ea.
909532		20"W x 16"D x 2"H	Ea.
909389	Replacement Sensor Pad	-	

**Note: All the sensors in this section require either Chair, Bed or MultiPro alarm unit*

BELT ALARM SENSORS




- Automatically resets when belt is secured
- Alarm sounds when closure is released
- Reinforced attachment holes and D-ring option for easy application
- Reminds resident to stay seated
- Alerts caregiver before resident leaves the chair
- Provides restraint-free wheelchair safety
- Universally sized belt secures to all wheelchairs

**Note: These above features apply only when combined with ChairPro, BedPro or MultiPro alarm unit.*

EZ RELEASE BELT ALARM SENSORS

Hook & Loop closure with red loop for easy release



New Magnetic Closure



914119 - EZ Release Belt w/Adjustable Side Release Loop Attachment


ReOrder #	Description	Size	U/M
909393	EZ Release Belt Sensor with Grommets	50"L	Ea.
914119	EZ Release Belt Sensor w/Adjustable Loop Attachment	45"L	Ea.



909392 - Metal Press Release Seat Belt with Grommets

SEAT BELT ALARM SENSORS

Nylon belt has metal press release buckle for restraint-free care



909385 Seat Belt w/Adjustable Loop Attachment

ReOrder #	Description	Size	U/M
909392	Metal Press Release Seat Belt with Grommets	50"L	Ea.
909394	Belt with Grommets	70"L	Ea.
914117	Metal Press Release Seat Belt w/Adjustable Loop Attachment*	45"L	Ea.
914118	Belt w/Adjustable Loop Attachment*	65"L	Ea.

**Note that when using the loop attachment for securing to the wheelchair no hardware is required.*

LIFT-OFF LAP CUSHION WITH SENSOR



- Releases easily as resident rises
- Easy-open Hook & Loop closure has red loop to facilitate release by hand
- Lap cushion is foam padded for comfort
- Covered in durable wipe-clean nylon



ReOrder #	Description	Size	U/M
909525	Lift-Off Lap Cushion w/Sensor	Fits 16"-18" Wheelchairs	Ea.

**Note: All the sensors in this section require either Chair, Bed or MultiPro alarm unit*

FLOORPRO SAFETY ALARM UNIT



Alerts caregivers when residents exert pressure on the fall mat. The FloorPro™ alarm includes the features below:

- Auto-Reset
- Nurse Call
- Cable Strain Relief
- Indicator Light*
- Battery Test
- Bumpers
- In-Room On-Off
- AC Adaptor
- Mountable
- High-Low Volume
- Two Buzzer/Two Music Tones
- Battery Included

* Indicator light flashes GREEN when monitoring, RED when alarming

ReOrder #	Description	U/M
909292	FloorPro Electronics Unit w/Accessories: (Accessories: Nurse Call Cord, Y-Connector, AC Adapter, Sensor Pad Wire and Batteries)	Ea.



All FloorPro™ Alarm Sets
Include Alarm Unit


Note: All the sets to the right are provided with an alarm unit, a floor mat or pad, and the appropriate number of sensor pads

FLOORPRO FALL MAT SETS

Combines an impact-absorbing fall mat with an alarm*

- Weight-sensitive floor pad activates alarm at bedside and/or nurse's station (when connected)
- High-density impact foam provides absorbing protection against fall-related injuries
- Covered in durable wipe-clean vinyl

FLOORPRO SOFT FALL BEDSIDE MAT ALARM SETS



68"L x 26"W x 2"H mat folds for easy storage

Alarm Sets (Alarm Unit, Sensor Pad & Accessories)

ReOrder #	Description	U/M
909275	FloorPro Bedside Mat Alarm Set (Single-Sided) 180-Day Sensor	Set
909276	FloorPro Bedside Mat Alarm Set (Double-Sided) 180-Day Sensor	Set
909280	Soft Fall Mat Only with Pockets 1-Year Guarantee	Ea.

* Double-side system is recommended for residents who can exit or fall from either side of bed.

Note: Two sensor pads are required for the single-sided Soft-Fall Bedside Alarm System

ReOrder #	# Sensor Pads
909275	Two Sensor Pads
909276	Four Sensor Pads
909283	One Sensor Pad

Pads provided with initial purchase

FLOORPRO FLOOR PAD ALARM SETS

Floor Pad measures 36"L x 26"W x 1"H

Alarm Sets (Alarm Unit, Sensor Pad & Accessories)

ReOrder #	Description	Guarantee	U/M
909283	Floor Pad Alarm Set (Single-Sided)	180 Days	Set
909284	Floor Pad Alarm Set (Double-Sided)	180 Days	Set
909286	Floor Pad Only (no sensor pad/electronics)	6 Months	Ea.



FLOORPRO SENSOR PADS ONLY

909290		180 Days	Ea.
909293	Sensor Pad Only	180 Days	4/pk
909294		180 Days	10/pk
909403		1 Year	Ea.

Original
Skil-Care
CORPORATION
Design

ECONO ALARM

Affordable alarm

- Secures to wheelchairs, geri-chairs and lounge chairs
- Reliable multi-directional magnetic pull-switch
- Adjustable-length cord clips to resident's clothing
- Comes with a Hook & Loop strap for attachment to all chairs
- Battery operated (batteries included)
- Fixed alarm level (85 db)



ReOrder #	Description	U/M
909207		Ea.
909208	Econo-Alarm	3/pk
909206		10/pk
909209	Replacement Magnet & Cord	5/pk

PERSONAL ALARM



Effective, economical alarm for use in wheelchair or bed

- Integral mounting hole
- Hook & Loop bed attachment strip provided
- Selector switch permits choice of high (85 db) or low (70 db) alarm volume
- Reliable multi-directional magnetic pull-switch
- Adjustable-length cord clips to resident's clothing
- Battery operated (batteries included)



ReOrder #	Description	U/M
909211	Personal Alarm	Ea.
909209	Replacement Magnet & Cord	5/pk

DOOR/WINDOW ALARM



Sounds an alert when a door or window is opened

- Prevents resident elopements
- Installs in seconds with adhesive backing
- Three-position switch: off, chime, buzzer (model 909223)
- Resets itself when door or window is closed (model 909223)
- Perfect for supply closets, medicine cabinets, windows, and doors
- Model 909223 includes a three 1.5 volt batteries (#357, #228, AG13 or similar)
- Model 909222 includes a 9-Volt Battery



Door Window Alarm Model 909223



Door Window Alarm Model 909222

ReOrder #	Description	U/M
909222	Door/Window Alarm with Magnet and Cord	Ea.
909223	Door/Window Alarm	5/pk

AC REPLACEMENT ADAPTOR



All of the Skil-Care alarms are powered by batteries, but the BedPro and FloorPro alarms also have the capability to be powered by AC (plug-in) power using this adaptor. Please note that when the BedPro or FloorPro alarm is powered by AC and in the event power is lost, power is immediately diverted to the batteries. Also, note when initially purchased, this adaptor is included with all BedPro or FloorPro Sets.




ReOrder #	Description	U/M
909348	AC Replacement Adaptor	Ea.

NURSES CALL ACCESSORY CONNECTOR KIT



This is included with the BedPro and Floor Pro on initial purchase

- Permits the MultiPro, BedPro or FloorPro alarms to be connected to nurses call system.



Y-Connector (909247)

ReOrder #	Description	U/M
909560	Nurse Call Accessory Connector Kit (Cord and Y-Connector)	Set
909247	Y-Connector for BedPro and MultiPro Alarms	Ea.
909242	Nurses Call Cable for BedPro and MultiPro Alarms	Ea.

AUXILIARY CONNECTOR CABLE

This cable is required to connect any of the Skil-Care Pro Series (ChairPro, MultiPro & BedPro) alarms to the sensors listed here. Although this cable is included with the Chair and MultiPro alarms, it is shown separately in the event of any listed sensors are used with the BedPro alarm unit.

ReOrder #	Description	Size	U/M
909368	Auxiliary Connector Cable	3ft Long	5/pk



Required for sensors below

ReOrder #	Description
909391	UnderSeat Sensor
909392	Seat Belt w/Buckle Sensor
909394	Seat Belt w/Buckle Sensor
909530	Gel-Foam Cushion Sensor 16" x 16" x 2"
909531	Gel-Foam Cushion Sensor 18" x 16" x 2"
909532	Gel-Foam Cushion Sensor 20" x 16" x 2"
909535	Leveling Pad Sensor

ALARM BAG



- Convenient way to mount any of the Skil-Care ProSeries alarms out of the residents reach
- Attaches to any wheelchair or bed
- Laundryable
- Zipper closure



ReOrder #	Description	U/M
911366	Alarm Bag	6/pk

REPLACEMENT MAGNETIC/ CORD

The Skil-Care Econo (909206, 909207, 909208) and the Personal (909211) alarms are activated when a magnet is pulled from the alarm. This magnet is tethered to the resident and when the resident leaves or rises from the chair, the magnet is released from the alarm and an audio sound is activated. When the magnet is placed back on the alarm, the audio sound stops. The magnet and cord listed below is included in the initial purchase of models 909207 or 909211



ReOrder #	Description	U/M
909209	Replacement Magnet & Cord	5/pk

ALARM BRACKETS



- Adjustable strap for easy installation
 - Keeps alarm out of resident's reach
 - May be used with alarms from most other manufacturers
- Can accommodate all Pro Series and Personal Alarms

CHAIR ALARM MOUNTING BRACKET



Alarm Shown for Demo Only

ReOrder #	Description	Size	U/M
909201	Alarm Mounting Bracket	Bracket: 2.5"W x 7"H, Belt: 26"L	Ea.

BED ALARM MOUNTING BRACKET




Alarm Shown for Demo Only

ReOrder #	Description	Size	U/M
909399	Bed Alarm Bracket	Bracket: 2.5"W x 7"H for 36" Beds	Ea.

CLIP-ON ALARM BRACKET



- Permits most alarms to be mounted away from resident
- Easily clips to any wheelchair or bed frame with a 1" tubular diameter
- Includes Hook & Loop strip to hold alarm in place



1" Tubular Diameter

ReOrder #	Description	Size	U/M
909390	Clip-On Alarm Bracket	7"L x 2.5"W	Ea.

UNIVERSAL ALARM WALL BRACKET

- Permits the MultiPro, ChairPro or BedPro (w/bracket removed) alarms to be mounted on the wall
- Convenient Hook & Loop or screw attachment to the wall and Hook & Loop attachment of alarm to the bracket



ReOrder #	Description	Size	U/M
909398	Universal Alarm Wall Bracket	7"L x 2.5"W	Ea.



STOP STRIP



Visual barrier controls wandering

- Important part of restraint reduction program
- Directs confused residents away from off-limits areas
- Secures to door frames by means of Hook & Loop or Magnet



Personal Alarm for
Models 909217 & 909218



ReOrder #	Description	Size	U/M
909217	Stop Strip w/Alarm Unit	12"H x 50"L	Ea.
909218	Stop Strip w/Alarm Unit	18"H x 50"L	Ea.
909220	Stop Strip w/Hook & Loop Mounts w/Stop Sign	12"H x 50"L	Ea.
909221	Yellow Hook & Loop Stop Strip w/o Stop Sign	12"H x 50"L	Ea.
909225	Stop Strip w/Magnetic Mounts w/Stop Sign	12"H x 50"L	Ea.
909226	Stop Sign ONLY with Hook & Loop	12"H x 12"L	Ea.

HIP-EASE (PAD)



Helps prevent hip injury with shock-absorbing high-density foam hip pads

- Foam pads position over hip joints
- Completely machine launderable (after removing foam pads)
- Removable vinyl-covered foam hip pads



ReOrder #	Waist Size	Size	U/M
911450	26" - 28"	XS	Ea.
911452	28" - 30"	S	Ea.
911454	30" - 34"	M	Ea.

ReOrder #	Waist Size	Size	U/M
911456	34" - 38"	L	Ea.
911458	38" - 42"	XL	Ea.
911460	42" - 46"	XXL	Ea.

ReOrder #	Description	U/M
911470	Replacement Vinyl-Covered Foam Pads	Pair

FIRE BLANKET

- Effective covering to extinguish grease fires
- Fiberglass fabric is fire-proof
- Can be easily stored in/on kitchen cabinets
- Full 30"W x 30"L to cover all kitchen pans
- Reusable
- Use only on non-flammable containers (metal, ceramic)
- Not to be used as a substitute for a fire extinguisher



ReOrder #	Description	Size	U/M
915310	Fire Blanket	30"W x 30"L	Ea.

SMOKERS APRON

A must for residents who smoke
NO-FLAME™ fabric prevents burn holes

- Full sized to cover lap and chest
- Lightweight for resident comfort
- New fabric that resists cracking



Wheelchair Smokers Apron



ReOrder #	Description	Size	U/M
906010	Wheelchair Smokers Apron White	30"L x 32"W	Ea.
906011	Wheelchair Smokers Apron Blue	30"L x 32"W	Ea.

ReOrder #	Description	Size	U/M
906015	Geri-Chair Smokers Apron White	39"L x 44"W	Ea.
906016	Geri-Chair Smokers Apron Blue	39"L x 44"W	Ea.

SHOWER CHAIR SAFETY BELT



For use with the shower chair, toilet or wheelchair. Attaches to all shower/bath seats and commode safety frames.

- Makes showering/bathing resident safer and easier
- Attaches with a loop tie
- No metal or components that rust
- May also be used as a wheelchair safety belt
- Mildew resistant



Original Skil-Care
 CONSULTANT DESIGN



ReOrder #	Description	Size	U/M
909110	Safety Belt for Shower/Toilet or Wheelchair	60"L	Ea.

SLIP-GRIP NON-SLIP MATTING



Multi-purpose matting provides added slip control

- Prevents cushions from slipping off chairs
- Holds plates so that unsteady eaters can feed themselves
- Keeps things in place and cuts easily into custom sizes & shapes



Can be used under a cushion for added grip




ReOrder #	Description	Size	U/M
909250	Non-Slip Matting Blue	144"L x 12"W Rolls	Ea.
909251	Non-Slip Matting Green	144"L x 12"W Rolls	Ea.
909252	Non-Slip Matting Raspberry	144"L x 12"W Rolls	Ea.

EZ Landing II™ **NEW**


FALL MAT

The EZ Landing II fall mats are characterized by their ease of care, one piece elastomeric foam surface, beveled low profile edges

- Reduces impact energy and helps reduce injury due to falls from the bed
- Low profile beveled edges help reduce tripping hazard
- Low profile allows resident to stand, walk or wheel over the mat without difficulty
- Slip resistant bottom keeps mat securely in place
- Slip resistant top and bottom



Beveled Edge




ReOrder #	Description	Size	U/M
911585	EZ Landing II Fall Mat	72"L x 24"W x 1"H	Ea.
911590		72"L x 30"W x 1"H	Ea.

SAFE & SOUND FALL MAT **NEW**

A multi-purpose impact reduction floor mat. It can be placed next to a bed or chair of at-risk patients to cushion and reduce injuries from falls. The mat is engineered with a thin, floor hugging surface with a beveled edge and anti-slide bottom.

- The beveled edge permits carts, chairs and IV stands to roll over easily and will not impede walking activities or patient stability.
- The mat is covered with an easy to clean, water resistant PVC material with anti-bacterial protection.
- The inside contains a unique polyurethane foam material that provides energy absorption to reduce impact from accidental falls.
- The anti-slide bottom keeps the mat stable on the floor and the non-skid surface prevents the patient or caregiver from sliding.
- The mat is latex free.
- Folding model comes with straps that allows it to be attached to the foot board above the floor
- This mat has an **18-month** guarantee.
- Available in two styles; folding and non-folding



ReOrder #	Description	Size	U/M
911570	Safe and Sound Folding Fall Mat	70"L x 24"W x .75"H	Ea.
911575	Safe and Sound Non-Folding Fall Mat	70"L x 24"W x .75"H	Ea.
911539	Replacement Storage Strap	Universal	Ea.

SAFESIDE FALL MAT

Impact absorption plus step-on stability

- Shock-absorbing bonded foam cushions fall impact
- Step-on stability prevents falls when resident walks on mat
- Trip resistant beveled edges
- Slip-resistant bottom keeps mat securely in place
- Waterproof covering protects against moisture and bacteria
- Bi-fold design for easy storage
- Includes strap that secures mat to footboard




ReOrder #	Description	Size	U/M
911572	Folding	68"L x 24"W x 1"H	Ea.
911577		68"L x 36"W x 1"H	Ea.
911539	Replacement Storage Strap	Universal	Ea.

ROLL-ON BEDSIDE FALL MAT

Tapered edge for easy wheelchair roll-on

- High-density, closed-cell foam mat absorbs a fall from the bed
- Trip resistant beveled edges
- Slip resistant bottom to keep mat securely in place
- Durable vinyl heat-sealed cover wipes clean for easy care
- Folds in half for easy storage and includes straps that secures mat to footboard



ReOrder #	Description	Size	U/M
911548	Roll-On Folding Bedside Fall Mat	68"L x 26"W x 1.5"H	Ea.
911549		68"L x 36"W x 1.5"H	Ea.
911539	Replacement Storage Strap	Universal	Ea.

SOFT-FALL BEDSIDE MAT

Surround bed with mats when rails are lowered

- High density polyurethane foam for extra protection
- Two easy storage models: a non-folding version that slips under the bed and a folding version with a special strap that secures mat to footboard (strap included)



ReOrder #	Description	Size	U/M
911535	Folding	68"L x 36"W x 4"H	Ea.
911536		68"L x 24"W x 4"H	Ea.
911530	Non Folding	68"L x 36"W x 4"H	Ea.

ReOrder #	Description	Size	U/M
911545	Folding	68"L x 36"W x 2"H	Ea.
911546		68"L x 24"W x 2"H	Ea.
911540	Non Folding	68"L x 36"W x 2"H	Ea.
911539	Replacement Storage Strap	Universal	Ea.

TRI-FOLD BEDSIDE MAT

Designed for easy storage

- High-density closed-cell foam protects against fall-related injuries
- Mats provide protection against fall-related injuries for low height beds



ReOrder #	Description	Size	U/M
911547	Tri-Fold	68"L x 26"W x 1"H	Ea.
911539	Replacement Storage Strap	Universal	Ea.

Skil-Care™

CORPORATION

AMBULATION



252200



252010



914384



Transfer Aids help with the transferring of residents to/from hospital beds and chairs. Skil-Care offers a wide assortment of Gait Belts, Transfer Belts, and Transfer Lift devices to provide lifting and walking assistance. If you are helping a person stand from a seated position, try the cotton gait belt. If infection control is a concern and you need something that wipes clean after use, try the Bio-shield Gait Belts.

TRANSFER BELTS WITH ADJUSTABLE HANDLES



Comfortable transfer belt when assisting in resident ambulation

- Full 4" width for added comfort, pressure distribution and reduced friction
- Three fixed handles and two sliding handles that assure optimal handle location
- Available with a metal or Delrin® buckle
- Elongated 5' length fits most bariatric residents
- Available in a padded or unpadded models

ReOrder #	Description	Size	U/M
252200	Transfer Belt, Non-Padded, Metal Buckle	Universal	Ea.
252201	Transfer Belt, Non-Padded, Delrin® Buckle	Universal	Ea.
252202	Transfer Belt, Padded, Metal Buckle	Universal	Ea.
252203	Transfer Belt, Padded, Delrin® Buckle	Universal	Ea.

WALKER BELT



This walker belt is intended to be used as an ambulatory assistance device.

- Three convenient handles
- Available in two sizes
- Convenient plastic side release buckles
- Wipe clean or launderable 2" nylon

ReOrder #	Description	Size	U/M
251011	Walker Belt	Small/Medium	Ea.
251022	Walker Belt	Medium/Large	Ea.

COTTON GAIT BELTS



Provides caregivers with a surface to grab when helping patients ambulate or transfer

- Additional supports helps reduce the risk of caregiver back injuries
- Provides lifting and walking assistance
- Durable 2"-wide belt with a 60"L, 72"L or any customizable size (contact customer service 800-431-2972).
- Easy to hand clean or launder

ReOrder #	Description	Size	U/M
252055	Heavy-Duty Webbing, Delrin® Buckle, Stars&Stripes	60"L	Ea.
252015	Heavy-Duty Webbing, Metal Buckle, Stars&Stripes	60"L	Ea.
252060	Standard Webbing, Delrin® Buckle, Pastel Stripes	60"L	Ea.
252070	Standard Webbing, Metal Buckle, Pastel Stripes	60"L	Ea.
252061	Heavy-Duty Webbing, Delrin® Buckle, Pastel Stripes	60"L	Ea.
252071	Heavy-Duty Webbing, Metal Buckle, Pastel Stripes	60"L	Ea.
252050	Standard Webbing, Delrin® Buckle, Pinstripe	60"L	Ea.
252010	Standard Webbing, Metal Buckle, Pinstripe	60"L	Ea.
252051	Heavy-Duty Webbing, Delrin® Buckle, Pinstripe	60"L	Ea.
252011	Heavy-Duty Webbing, Metal Buckle, Pinstripe	60"L	Ea.
252073	Heavy-Duty Webbing, Metal Buckle, Stars&Stripes	72"L	Ea.
252072	Heavy-Duty Webbing, Metal Buckle, Pastel Stripes	72"L	Ea.
252075	Heavy-Duty Webbing, Delrin® Buckle, Pinstripe	72"L	Ea.
252074	Heavy-Duty Webbing, Metal Buckle, Pinstripe	72"L	Ea.

INFECTION CONTROL GAIT BELTS



These Bio-shield Gait Belts provide lifting and walking assistance while addressing infection control

- 2" wide plastic coated polyurethane webbing
- Smooth liquid proof wipeable surface
- Disinfects easily with spray on disinfectant or bleach type cloth
- Sealed construction with no stitching prevents penetration of bacteria
- Secures with durable Delrin® side-release buckle
- Comes in 60" or 72" lengths



ReOrder #	Description	Size	U/M
914376	Bio-shield Wipe-Clean - Yellow	60"L	Ea.
914378	Bio-shield Wipe-Clean - Green	60"L	Ea.
914380	Bio-shield Wipe-Clean - Blue	60"L	Ea.
914384	Bio-shield Wipe-Clean - Stars&Stripes	60"L	Ea.
914387	Bio-shield Wipe-Clean - Orange	60"L	Ea.
914386	Bio-shield Wipe-Clean - Stars&Stripes	72"L	Ea.
914377	Bio-shield Wipe-Clean - Yellow	72"L	Ea.
914379	Bio-shield Wipe-Clean - Green	72"L	Ea.
914382	Bio-shield Wipe-Clean - Blue	72"L	Ea.
914388	Bio-shield Wipe-Clean - Orange	72"L	Ea.
913150	Bio-shield Gait Belt Handles - See Page 107		Pair

GEO-PATTERN GAIT BELTS



Provides lifting and walking assistance

- Available in four (4) attractive Geo-Pattern designs
- Available in Delrin® or metal buckles
- Reinforced nylon for extra strength
- Easy to hand clean or launder
- All belts are 60”L

Delrin Buckle	Pattern	Metal Buckle	Size	U/M
252021	Geo-Pattern A	252020	60”L	Ea.
252024	Geo-Pattern B	252023	60”L	Ea.
252027	Geo-Pattern C	252026	60”L	Ea.
252030	Geo-Pattern D	252029	60”L	Ea.

EZ CLEAN GAIT BELTS



252000



252004



252250

Provides lifting and walking assistance

- The polypropylene webbing can be cleaned easily with detergents or anti-bacterial agents. Do not launder
- Available with an easy release Delrin® buckle or a heavy duty metal buckle
- Cost effective gait belt
- Soft 2” webbing is comfortable to the resident and caregiver.
- Available in two colors, blue or orange

Delrin Buckle	Color	Metal Buckle	Size	U/M
252000	Blue	252002	60”L	Ea.
252001	Blue	252003	60”L	24/pk
252004	Orange	252006	60”L	Ea.
252005	Orange	252007	60”L	24/pk
252250	Blue w/Two Handles	N/A	48”L	24/pk
252252	Blue w/Two Handles	N/A	60”L	24/pk

GAIT BELT HANDLES



Deluxe Rubber Gait Belt Handles



PathoShield Gait Belt Handle



Nylon Gait Belts Handles



PathoShield Gait Belt Handle (Belt sold separately)

The gait belt handles fits over 2” wide gait belts

- Provides a more secure and comfortable way to grip the gait belt
- Available in three webbing types.

ReOrder #	Description	U/M
913145	Nylon Gait Belt Handles	6/pk
914389	Deluxe Rubber Gait Belt Handles	6/pk
913150	Bio-shield Gait Belt Handles	Pair

SUPER-SLINGS FOR LIFTING & TRANSFER



Designed for Two Caregivers

2 Handles



Designed for One or Two Caregivers

4 Handles



Designed for use up to Four Caregivers

8 Handles



Designed for use up to Six (6) Caregivers

12 Handles



Facilitate resident transfer, lifting, and positioning while reducing exposure to back injury

All Super-Slings feature:

- Friction-reducing bottoms that ease positioning and transfer
- Reinforced stitching at handles
- Two-inch-wide nylon webbing for resident safety
- Durable nylon fabric

ReOrder #	Description	Handles	Size	U/M
914123	Supports up to 200lbs	2 Handles	20"L x 9"W	Ea.
253020	Supports up to 250lbs	4 Handles	40"L x 36"W	Ea.
253010	Supports up to 500lbs	8 Handles	46"L x 46"W	Ea.
253030	Supports up to 500lbs	12 Handles	72"L x 36"W	Ea.


Super-Slings may be used for emergency evacuations when wheelchairs or gurney's are either unavailable or not practical.

Note: The number of caregivers involved in positioning and transferring should be determined by a qualified caregiver

TLC POSITIONING PAD



555012 Low Friction Bottom



Reduces caregiver back-injury when residents are positioned in bed or chair

- Wicks urine away from sensitive skin
- Covered with a soft launderable Cozy Cloth fabric
- Water proof barrier protects bed from body fluids
- Stands up to repeated laundering
- Strong web handles facilitate positioning and transfer
- Supports up to 250lbs

ReOrder #	Description	Size	U/M
555012	TLC Positioning Pad w/Low-friction nylon bottom	48"L x 40"W	Ea.
555015	TLC Positioning Pad w/Hook & Loop compatible bottom (use w/In-Bed System 555036)	48"L x 40"W	Ea.
555018	TLC Positioning Pad w/Hook & Loop compatible bottom	40"L x 36"W	Ea.

BARIATRIC



705024



703475



503358



Skil-Care offers several products designed to accommodate additional weight bearing capabilities as well as comfort. These products are engineered to meet the needs of bariatric clients.

BARIATRIC HEEL-FLOATS

Heel-Float



Heel-Float Plus+




Heel-Float Wipe-clean

Off-loads the heel over an air cavity to eliminate pressure

- Floats the heel over an air cavity for zero pressure
- Provides air circulation to aid healing of existing ulcers
- Water-based gel insert lowers pressure beneath ankle
- Anti-foot-drop construction controls dorsal flexion
- Launderable and reusable (only when removal of gel pack)

ReOrder #	Description	Size	U/M
503036	Bariatric Heel-Float	5" wide	Ea.
503046	Bariatric Heel-Float Plus+	5" wide	Ea.
503049	Bariatric Heel-Float - Wipe-Clean	5" wide	Ea.
503099	Replacement Gel Pack for Heel Float Series	Large	Ea.

BARIATRIC HEELS-OFF



Suspends heels off the bed for total heel pressure relief

- Off-loads both heels
- Prevents formation of pressure sores
- Aids in healing existing ulcers
- Elevates legs to relax spine
- Low-Shear II cover protects against skin-damaging friction
- Wipes clean for easy care
- Optional gel adds superior comfort and cooling

ReOrder #	Description	Size	U/M
503150	Heels-Off - Medium/Large	16"L x 24"W x 7"H	Ea.
503165	Heels-Off w/Gel Pad Medium/Large	16"L x 24"W x 7"H	Ea.
503167	Bariatric Heels-Off w/Gel Pad	16"L x 26"W x 7"H	Ea.

BARIATRIC ECONO FOOTREST EXTENDER



Effective for lower extremity positioning

- Raises footrest for improved posture and seating alignment
- Prevents foot drop and feet from slipping behind footrests
- Cushioned back protects lower leg against contact with wheelchair
- Secures easily with quick-release straps
- Covered in durable, wipe-clean vinyl

ReOrder #	Description	Size	Foot Pad Thickness	U/M
703294		Fits	1 Inch	Ea.
703296	Econo-Footrest Extender	20" - 24"	2 Inch	Ea.
703297		W/C	3 Inch	Ea.

BARIATRIC WHEELCHAIR LEG PAD



Protects legs from contact with metal wheelchair parts

- New E-Z clip buckles and straps simplify attachment
- Covered in durable, wipe-clean vinyl
- Prevents resident's feet from slipping behind footrests
- Foam-padded, rigid-backed leg pad

ReOrder #	Description	Size	U/M
703072	Bariatric Wheelchair Leg Pad	Fits 20" - 24" Wheelchairs	Ea.

BARIATRIC SWING-AWAY FOOT SUPPORT



Swings away for easy wheelchair transfers

- Specially designed to fit wheelchairs 22" to 30" wide
- Cushioned side and back panels control foot positioning
- Covered in durable, wipe-clean nylon
- Easy installation
- Specify left, right, or pair when ordering

ReOrder #	Description	Size	U/M
703475	Bariatric Swing-Away Support - Left	Fits 22" - 30" Wheelchair	Ea.
703476	Bariatric Swing-Away Support - Right	Fits 22" - 30" Wheelchair	Ea.
703477	Bariatric Swing-Away Support - Pair (Left/Right)	Fits 22" - 30" Wheelchair	Pair

BARIATRIC GERI-SLEEVES

LIGHT TONE



PADDED



DARK TONE



New extra-heavy stockinette maintains skin integrity

- Breathable cotton-blend Lycra® protects against skin tears
- Thumb-hole cutouts keep arm sleeves in place
- Leg sleeves prevent skin shear at the heels and calves
- Protects I.V. and wound sites
- Reusable and launderable



ReOrder #	Description	U/M	Size	Stripe Color	Length
503357	Light Tone Geri-Sleeves - Arm	Pair	Large/Bariatric	Blue	20"
503360	Light Tone Geri-Sleeves - Leg	Pair	Universal	Green	16"
503366	Dark Tone Geri-Sleeves - Arm	Pair	Large/Bariatric	Blue	20"
503368	Dark Tone Geri-Sleeves - Leg	Pair	Universal	Green	16"
503358	Light Tone Geri-Sleeves - Arm	25 Pr/Pk	Large/Bariatric	Blue	20"
503361	Light Tone Geri-Sleeves - Leg	25 Pr/Pk	Universal	Green	16"
503367	Dark Tone Geri-Sleeves - Arm	25 Pr/Pk	Large/Bariatric	Blue	20"
503369	Dark Tone Geri-Sleeves - Leg	25 Pr/Pk	Universal	Green	16"

ReOrder #	Description	Size	U/M
503394	Padded Geri-Sleeves	Large/Bariatric	Ea.

BARIATRIC KNEE ELEVATOR



- Decreases leg and lower back pain
- Elevates the legs to relax the spine
- Covered in durable, wipe-clean anti-microbial vinyl
- Improves circulation
- Foam padded for comfort

ReOrder #	Description	Size	U/M
555046	Bariatric Knee Elevator Half Round	26"L x 8"D x 7"H	Ea.

BARIATRIC GAIT BELTS



Stars&Stripes - Heavy-Duty



Pastel - Heavy-Duty



Pinstripe - Heavy-Duty



Standard Gait Belt



Bio-shield Gait Belt

Extra strength Gait Belts that provides lifting and walking assistance for bariatric residents

- 72"-long belts accommodate your largest residents
- Strong 2"-wide webbing
- Secure non-slip buckles
- Fabric (Standard Gait Belts) are launderable
- Bio-shield Gait Belts assistance in lifting and walking while addressing infection control

ReOrder #	Description	Guarantee	U/M
252072	Heavy-Duty Webbing, Metal Buckle, Pastel Stripes - 72"L	6 months	Ea.
252073	Heavy-Duty Webbing, Metal Buckle, Stars&Stripes - 72"L	6 months	Ea.
252074	Heavy-Duty Webbing, Metal Buckle, Pinstripe Stripes - 72"L	6 months	Ea.

ReOrder #	Description	Guarantee	U/M
914382	Bio-shield Wipe-Clean - Blue - 72"L	1 Year	Ea.
914386	Bio-shield Wipe-Clean - Stars&Stripes - 72"L	1 Year	Ea.
914388	Bio-shield Wipe-Clean - Orange - 72"L	1 Year	Ea.
914377	Bio-shield Wipe-Clean - Yellow - 72"L	1 Year	Ea.

BARIATRIC SofTop™ WHEELCHAIR LAP TRAY w/STRAPS



Easy Release
Hook & Loop Strap

- Foam-padded surface is gentle on arms and elbows
- Helps achieve proper body position, balance, and alignment
- Padded to provide a comfortable, low-pressure surface for arms/elbows
- Around-the-chair strap prevents resident from sliding forward
- Easy-open Hook & Loop closures
- Specially designed to fit wheelchairs 24" to 28" wide

ReOrder #	Description	Size	U/M
705024	SofTop Lap Tray Nylon Cover	Fits 24" - 28" Wheelchairs	Ea.

BUDGET BARIATRIC FOAM CUSHION



Weight Tested to 500 lbs

This light weight high density foam cushion provides comfort and pressure relief

- High density foam provides comfort and pressure relief
- Light weight for easy carrying
- Stable surface
- Can be used under all temperature conditions
- Low-Shear II cover reduces skin-damaging friction and wipes clean for easy care
- Safety strap secures cushion to wheelchair
- Non-slip bottom on outer cover provides extra safety

ReOrder #	Description	Size	U/M
915060		20"W x 16"D x 4"H	Ea.
915061		20"W x 18"D x 4"H	Ea.
915062		22"W x 16"D x 4"H	Ea.
915063		22"W x 18"D x 4"H	Ea.
915064	Budget Bariatric Foam	24"W x 18"D x 4"H	Ea.
915065	Cushion w/LSII Cover	24"W x 20"D x 4"H	Ea.
915066		26"W x 18"D x 4"H	Ea.
915067		26"W x 20"D x 4"H	Ea.
915068		28"W x 18"D x 4"H	Ea.
915069		28"W x 20"D x 4"H	Ea.

BARI-GEL FOAM CUSHION



Original Skil-Care CORPORATION Design


Weight Tested to 500 lbs

Double-chamber construction relieves pressure for heavy residents

- High-density foam prevents bottoming-out
- Low-Shear I launderable outer cover reduces skin damaging friction
- Incontinent-proof vinyl inner cover

ReOrder #	Description	Size	U/M
915110		20"W x 18"D x 3"H	Ea.
915120		22"W x 18"D x 3"H	Ea.
915130		24"W x 18"D x 3"H	Ea.
915135	Bariatric Gel-Foam	26"W x 18"D x 3"H	Ea.
915133		26"W x 20"D x 3"H	Ea.
915134		28"W x 18"D x 3"H	Ea.
915137		28"W x 20"D x 3"H	Ea.
915139		30"W x 20"D x 3"H	Ea.

BARIATRIC FOAM WEDGE



Weight Tested to 550 lbs

This Gel-Infused Visco foam wedge controls forward leaning

- Designed to reduce sitting pressure
- Controls forward leaning
- Safety strap secures cushion to wheelchair
- Nylon Cover
- 3-layers of foam provides added comfort and pressure reduction

ReOrder #	Description	Size	U/M
910140		20"W x 18"D x 4.5" x 3"H	Ea.
910145	Bariatric Foam Wedge	22"W x 18"D x 4.5" x 3"H	Ea.
910150		24"W x 18"D x 4.5" x 3"H	Ea.
910155		26"W x 18"D x 4.5" x 3"H	Ea.

TRI-FOAM GEL-INFUSED VISCO BARIATRIC CUSHION



Gel-Infused Visco Foam

Weight Tested to 550 lbs

Maximum comfort and stability

- Top layer of Gel-Infused Visco foam provides pressure relief
- Middle layer of high-density foam adds comfort and cushioning
- A **Firm Foundation** controls sling-seat hammocking
- Three-part design prevents bottoming-out
- Low-Shear II wipe-clean outer cover that reduces skin damaging friction


Contoured Bottom - For Sling Seat Wheelchairs

ReOrder #	Description	U/M	ReOrder #	Description	U/M
754810	20"W x 16"D x 3"H	Ea.	754821	24"W x 20"D x 4"H	Ea.
754812	20"W x 18"D x 3"H	Ea.	754823	26"W x 18"D x 4"H	Ea.
754814	22"W x 16"D x 4"H	Ea.	754825	26"W x 20"D x 4"H	Ea.
754816	22"W x 18"D x 4"H	Ea.	754827	28"W x 18"D x 4"H	Ea.
754818	24"W x 18"D x 4"H	Ea.	754829	28"W x 20"D x 4"H	Ea.

Flat Bottom - For Solid Seat Wheelchairs

754833	20"W x 16"D x 3"H	Ea.	754843	24"W x 20"D x 4"H	Ea.
754835	20"W x 18"D x 3"H	Ea.	754845	26"W x 18"D x 4"H	Ea.
754837	22"W x 16"D x 4"H	Ea.	754847	26"W x 20"D x 4"H	Ea.
754839	22"W x 18"D x 4"H	Ea.	754849	28"W x 18"D x 4"H	Ea.
754841	24"W x 18"D x 4"H	Ea.	754851	28"W x 20"D x 4"H	Ea.


CUSTOMIZED BARIATRIC POSITION PLUS CUSHION




All custom cushions and wedges are made to user specification and measurement. Call Skil-Care customer service at **800-431-2972** for complete ordering information.

Weight Tested to 500 lbs


Available in Four Styles:



Level cushion for solid-seat wheelchairs



Level/contoured cushion for sling-seat wheelchairs



Wedge cushion for solid-seat wheelchairs



Wedge/contoured cushion for sling-seat wheelchairs

Four separate gel-filled chambers provides stability.

- Available in sizes up to 28" x 28" (in 2" increments)
- Custom-made to user's special seating requirements
- Water-based gel relieves pressure and cools the skin to prevent breakdown caused by perspiration and heat build-up

- Side gel pads improve lateral stability by controlling left-right leaning
- Rear gel pad provide pressure relief in the coccyx area
- Front gel pad controls forward sliding by eliminating posterior pelvic rotation
- High-density foam adds comfort and prevents bottoming out
- Incontinent-proof vinyl inner cover prevents urine from entering cushion wheelchairs

ReOrder #	Description	U/M
754650	Custom Sized Solid Seat Level Cushions	Ea.
754655	Custom Sized Solid Seat Wedge Cushions	Ea.
754660	Custom Sized Sling Seat Level Cushions - Contoured	Ea.
754665	Custom Sized Sling Seat Wedge Cushions - Contoured	Ea.

Note: Available for sling-seat and solid-seat wheelchairs

SkiL-Care™
CORPORATION



303210



306130



301300

PHYSICAL RESTRAINTS



Restraints provide numerous functions:

1. As a reminder or to prevent the resident from leaving the wheelchair or bed
2. Limits arm and leg mobility (mitts, limb holders etc)
3. Assists in achieving optimal wheelchair posturing.
4. Adds safety to caregiver and resident, especially those who are confused

Although there are regulations that call out that restraints should be used only if necessary, in many scenarios a restraint is necessary for a specific resident. Regulations also state that if a restraint is used, it should be used only when necessary. These restraints were developed and intended only to be used in an acute care or long-term care setting. Refer to our control point chart located in the literature section of our website (skil-care.com). SkiL-Care offers a wide selection of registered restraints.

CUSHION SLIDER BELTS



- Provides slide control and waist support
- Comfortable, six-inch-wide foam waist belt
 - Soft padded crotch piece prevents sliding
 - Completely launderable
 - Includes removable urine-proof crotch pad cover for incontinent residents
 - Universally sized

ReOrder #	Description	Size	U/M
303210	Cushion Slider Belt	Universal	Ea.

PELVIC HOLDER



- Prevents forward sliding
- Effective, easy-application slider control
 - Tapered, hourglass design fits comfortably between resident's legs
 - Waist straps secure behind wheelchair, out of resident's easy reach
 - Available in Poly Weave or Cooline Mesh

ReOrder #	Description	Size	U/M
303311	Poly Weave	Small/Medium	Ea.
303312		Medium/Large	Ea.
303321	Cooline Mesh	Small/Medium	Ea.
303322		Medium/Large	Ea.

CUSHION BELT



- Comfortable belt that reminds resident not to leave the wheelchair or bed
- Easy to apply in wheelchair or bed
 - Available with your choice of closures
 - For wheelchair use, select either ties or quick-release buckle
 - In bed, the closure allows resident to roll from side to side
 - Durable and completely launderable

ReOrder #	Description	Size	U/M
301200	Tie Closure - Wheelchair and Bed	Padded Foam - 23"L x 6"W, Tie Straps - 70"	Ea.
301201	Buckle Closure - Wheelchair Only	Padded Foam - 23"L x 6"W, Webbing min/max 41" / 70"	Ea.

ECONO-BELT




- Cost effective and easy to apply
- Two-inch-wide nylon belt
 - Attaches to wheelchair kick-spurs, out of resident's easy reach
 - Slider buckles allow full size adjustment for any resident

ReOrder #	Description	Size	U/M
701020	Econo-Belt	Fits up to 50" waist	Ea.

ADJUSTA-LOOP CUSHION BELT



301202



301205

Comfortable belt that reminds resident not to leave the wheelchair

- Eliminates the need for knots
- Convenient securing loops slip over wheelchair kick-spurs
- Slider buckles adjust webbing to resident size
- Removes quickly in an emergency
- Easy-on, easy-off application saves time and effort
- Comfortable and foam padded
- Durable and completely launderable
- Available in 1.5" polyester webbing or 2" heavy duty webbing

ReOrder #	Description	Size	U/M
301202	Adjusta-Loop Cushion Belt Standard 1.5" Webbing	Padded Foam 23"L x 6"W, Web Belt - 84" (7")	Ea.
301205	Adjusta-Loop Cushion Belt Heavy Duty 2" Webbing	Padded Foam 23"L x 6"W, Web Belt - 84" (7")	Ea.

ROLL BELT



- Allows resident to roll from side to side
- "Gentle reminder" not to leave the bed
- Extra-long straps secure under bed, out of resident's reach
- Durable and completely launderable
- Universally sized

ReOrder #	Description	Size	U/M
301300	Roll Belt	Universal	Ea.

SHEEPSKIN LIMB HOLDER



Limits mobility of limbs in bed

- Made of plush, high pile synthetic sheepskin
- Softly secures wrists or ankles
- Non-slip buckle fastens easily and stays in place
- Extra-long straps can be secured out of resident's reach
- Launderable synthetic sheepskin

ReOrder #	Description	Tie Strap	U/M
306010	Sheepskin Limb Holder	Double	Pair
306020	Sheepskin Limb Holder	Double	36 Pair/Cs

ECONOMY LIMB HOLDER

Limits mobility of limbs in bed

- Disposable foam or sheepskin limb holders
- Soft and secure wrist or ankle control
- Limb holders are easily positioned with Hook & Loop closures
- Non-slip buckle secures webbing around wrists or ankles
- Choose either single or double tie straps



Model 306044

ReOrder #	Description	ReOrder #	Description	U/M
306044	Foam - Single Strap	306045	Foam - Double Strap	Pair
306040	Foam - Single Strap	306048	Synthetic Sheepskin - Single Strap	36 Pair/Cs
306041	Foam - Double Strap	-	-	36 Pair/Cs

SAFETY VEST



Serves as a reminder not to leave the chair

- Designed for easy application
- Straps secure behind chair, out of resident's reach
- Bindings are color coded for proper size and fit
- Useful for non-ambulatory resident who does not slide
- Choice of two durable, launderable fabrics: Poly Weave or Cooline Mesh

ReOrder #	Description	Size	Binding Color	U/M
301101		Small	Red	Ea.
301102	Poly Weave	Medium	Green	Ea.
301103		Large	Yellow	Ea.
301111		Small	Red	Ea.
301112	Cooline Mesh	Medium	Green	Ea.
301113		Large	Yellow	Ea.
301121	Economy - Poly Weave	Small	Red	Ea.
301122	Single Patient Use	Medium	Green	Ea.
301123		Large	Yellow	Ea.

SLEEVED JACKET



The zipper-back closure makes this jacket easy to apply for the caregiver and very difficult to remove for the resident

- Waist straps provide wheelchair safety
- Posture control is provided by shoulder straps
- Available in Poly Weave or Cooline Mesh

ReOrder #	Description	Size	U/M
304203		Small	Ea.
304204	Poly Weave	Medium	Ea.
304205		Large	Ea.
304213		Small	Ea.
304214	Cooline Mesh	Medium	Ea.
304215		Large	Ea.

SLEEPER JACKET



Prevents the resident from leaving the bed

- Far more secure than vest-style restraints
- Sleeved jacket is difficult for resident to remove
- Zipper back makes application easy
- Wrap-around waist belt allows resident to sleep in any position
- All straps secure under the bed, safely out of resident's reach

ReOrder #	Description	Size	U/M
305103	Poly Weave	Small	Ea.
305104	Poly Weave	Medium	Ea.
305105	Poly Weave	Large	Ea.
305113	Cooline Mesh	Small	Ea.
305114	Cooline Mesh	Medium	Ea.
305115	Cooline Mesh	Large	Ea.

PADDED MITTS



Limits picking, scratching, and the ability to grab tubes

- Lightweight and comfortable
- Cushioned with 3/4-inch-thick foam padding
- CoolLine Mesh cover promotes air circulation and allows for inspection
- Easy to apply and completely launderable
- Convenient closure eliminates cumbersome shoelace ties
- Can be used in bed or wheelchair
- Packaged with tie straps to limit resident's arm movement
- Launderable
- Fits either hand

ReOrder #	Description	Size	U/M
306110	Padded Mitts	Universal	Pair

If there is doubt about use of a restraint, a Skil-Care™ representative may be helpful in suggesting a non-restrictive functional equivalent.

PADDED-PLUS MITTS



Limits picking, scratching, and the ability to grab tubes

- Amply cushioned with 1.5" thick padding
- Easy-to-secure closure eliminates cumbersome shoelace ties
- Added Hook & Loop for strap placement
- For use in wheelchair or bed
- Packaged with tie straps to limit resident's arm movement
- Launderable
- Fits either hand

ReOrder #	Description	Size	U/M
306120	Padded-Plus Mitts	Universal	Pair

RIGID PALM PADDED MITTS



Limits picking, scratching, and the ability to grab tubes

- Rigid insert restricts resident's ability to grab tubes
- Cushioned with 1" thick foam padding
- Separate finger channels control hand contracture
- Mesh top piece promotes air circulation and allows for visual inspection
- Packaged with separate straps to limit resident's arm movement*
- Launderable after removing the rigid insert
- Fits either hand

ReOrder #	Description	Size	U/M
306115	Rigid Palm Padded Mitts	Universal	Pair

***Note that the extra strap is used to immobilize the limb**

FINGER SEPARATOR MITTS



Limits picking, scratching, and the ability to grab tubes

- Cushioned with 2" thick padding
- Stitched finger separations prevent the user to rotate the mitt
- CoolLine Mesh top piece promotes air circulation and allows for visual inspection
- Easy to apply
- Packaged with separate straps to limit resident's arm movement*
- Launderable
- Fits either hand

ReOrder #	Description	Size	U/M
306130	Finger Separator Mitts	Universal	Pair

E-Z VIEW PADDED MITTS



Limits picking, scratching, and the ability to grab tubes

- Effective hand mitt with easy inspection capability
- Effective finger and hand control
- Easy inspection by caregiver
- Fits either hand
- One size fits most adults
- Packaged with separate straps to limit resident's arm movement*
- Launderable
- Fits either hand

ReOrder #	Description	Size	U/M
306125	E-Z View Padded Mitts	Universal	Pair

**Note that the extra strap is used to immobilize the limb*



902011



904150





Skil-Care has what you need to give your residents an extra layer of care.



910015

TRACH TUBE TIE/HOLDER



- Hook & Loop fastener to hold tracheotomy tube in place
- Soft foam collar is comfortable and non-irritating to resident's skin
 - Secures with convenient Hook & Loop fasteners
 - Disposable

ReOrder #	Description	Size	U/M
904150	One Piece Construction	Small/Medium (Neck 8"-18")	Dz
904151	Two Piece Construction	Medium/Large (Neck 16"-23")	Dz
904110	Two Piece Construction - Adjustable Holder	Universal	Dz

TELEMETRY UNIT POUCH



- Protects expensive telemetry units
- Assists in providing more accurate readings
 - Straps will adjust to resident's size
 - Limits movement when tied around the neck & waist

ReOrder #	Description	Size	U/M
901130	Telemetry Unit Pouch - Cloth	Universal	36/Cs

BIBS



- Protection from food and liquids landing on the lap
- Reversible absorbent microfiber on one side and vinyl on the other
 - Adjustable neck size with Hook & Loop attachment
 - Bound on the perimeter for added durability
 - Deep, pleated, 8" catch pocket for effective food pickup (910015)
 - Launderable and latex free

ReOrder #	Description	Size	U/M
910025	Bib with no pocket (as shown)	18"W x 30"L	12/Cs
910015	Bib with pocket	18"W x 30"L	12/Cs

POUCH ARM SLING



- Universal arm support for either arm
- Adjustable shoulder strap
 - Covered with a soft launderable Cozy Cloth fabric
 - Non-slip buckle

ReOrder #	Description	Size	U/M
902011	Pouch Arm Sling	Small-Medium	Ea.
902012		Medium-Large	Ea.

BATH MITT



- Dignified washing surface
- Soft sheepskin mitt calms as it cleanses
 - Mesh pocket holds soap and exfoliates the skin
 - Machine launderable for infection control
 - Durable and reusable

ReOrder #	Description	Size	U/M
913030	Bath Mitt	Universal	Ea.



INCONTINENT SUPPLIES



101080



102010



102015



The Skil-Care urinary devices help resident maintain dignity. The bag holders hang under the wheelchair and the straps help conceal and structure catheter tubes.

URINARY DRAINAGE BAG HOLDER



Standard Drainage Bag



Thru-View Drainage Bag

- Conceals urinary drainage bags from public view
- Maintains dignity of catheterized residents
 - Designed to hang under wheelchairs, geri-chairs, and from bed rails
 - Made from durable, easy-to-clean vinyl
 - Heat-sealed to contain spills
 - Universally sized to hold all drainage bags
 - Provides safe storage for tubing
 - Thru-View model features a clear window to allow fluid level monitoring
 - Attached with Hook & Loop straps

ReOrder #	Description	Size	U/M
102010		13"W x 8"H	Ea.
102012	Standard	13"W x 8"H	24/Cs
102011		13"W x 8"H	100/Cs
102015		13"W x 8"H	Ea.
102020	Thru-View*	13"W x 8"H	24/Cs.

* Note: Thru-View models have a window on one side only and can be reversed

CATHETER LEG STRAP



Reduces accidental dislocation of catheter tubes

- Comfortable, cloth-backed foam band wraps around resident's leg
- Catheter tube is held in place with a Hook & Loop fastener system
- Eliminates the need for hard-to-remove tape
- Non irritating to sensitive skin
- Completely adjustable and easy to apply
- Can be cut for optimizing the size

ReOrder #	Description	Size	U/M
101080	Catheter		Dz
101090	Leg Strap	30"L	36/Cs

RECREATION & EXERCISE



708030



914676



914154



Active resistive exercise improves strength, flexibility and the ability to perform normal activities such as feeding and walking. A moderate exercise program helps minimize falls and dependence on the caregivers. In a standard healthcare facility, an exercise program is conducted in an exercise room and during a scheduled period. The exercise equipment presented by Skil-Care Corp. can be used independently, if capable by the resident in his/her room and/or home environment.

Ex-Box™



Multi-Muscle training for legs, back, shoulders and arms

- Provides a wide range of upper and lower body resistance and range of motion exercises
- Resistance can be varied with tube diameter or length
- Builds arm-leg coordination and strength
- Can be used in any fixed position chair or in a locked wheelchair
- Includes a set of exercise handles
- Easy-attachment sockets permit quick adjustment
- Includes one set of resistance bands
- Is equally affective for upper and lower body

ReOrder #	Description	U/M
708030	ExBox w/Bands	Ea.
914539	Replacement Hand Grips (2)	Pair

Aqua Weight



Safe and convenient weights that are easy to use and transport

- Five convenient models
- Low shipping cost
- Easy to fill, hold and store
- Reduces injury to feet if dropped
- Adjustable model has 1-6 lb fill lines

ReOrder #	Description	Weights	U/M
914676		2lbs	Ea.
914677	Aqua Weight	4lbs	Ea.
914678		6lbs	Ea.
914679		8lbs	Ea.
914680	Aqua Adjustable Weight - Clear Face	1-6lbs	Ea.

ARCTIC THERMAL SLEEVE



Ankle

Instep/Sole




Elbow

Knee

- Provides heat or cold therapy
- Thermal sleeves provide for easy and stable positioning on lower leg, elbow, knee and shoulders
- Sleeves provide protection against direct contact of hot/cold pack with the skin.
- Sleeve is durable and launderable
- Universally sized
- Shipped with one set of gel packs

ReOrder #	Description	U/M
913100	Arctic Thermal Sleeve - Ankle	Ea.
913101	Arctic Thermal Sleeve - Instep/Sole	Ea.
913102	Arctic Thermal Sleeve - Elbow	Ea.
913103	Arctic Thermal Sleeve - Knee	Ea.
503092	Replacement Gel Pack for Arctic Sleeve	Ea.

Aqua Air Exerciser



These air exercisers are easy to store and are a perfect addition to your water workouts

- Can be filled with air from 3 to 10 lbs of resistance
- Convenient fill valve
- Easy to hold and store
- Low shipping cost
- Should not be used as a floatation device

ReOrder #	Description	Max Resistant Buoyant Force	U/M
914800		3 lbs	Ea.
914801	Aqua Air	5 lbs	Ea.
914802	Exerciser	8 lbs	Ea.
914803		10 lbs	Ea.

BED LADDER



This bed ladder is designed to enable residents to pull themselves upright in a bed and can also be used for upper body exercise.

- The bed ladder enables and facilitates the resident to sit up in a bed
- The ladder provides mild exercise for the upper body
- Easily attaches to any bed
- Can be grabbed easily
- Launderable

ReOrder #	Description	Size	U/M
914154	Bed Ladder	12' L	Ea.

CUSHION GRIP



Cloth Cushion Grip

Cozy Cloth Gel Grip



For residents with involuntary hand closure

- Prevents residents from digging their fingernails into their palms
- Provides a mild form of hand exercise
- Elastic band keeps cushion from being dropped
- Filled with resilient synthetic sheepskin (201020 & 201030) or gel (201170 & 201175)
- Cover is launderable

ReOrder #	Description	Size	U/M
201020	Polyester Fill	Universal	6/pk
201030	Cushion Grip		36/Cs
201170	Gel Grip	Small/Medium Hands	6/pk
201175		Medium/Large Hands	6/pk

FINGER CONTRACTURE CUSHION



Separates fingers and provides anti-contracture hand therapy

- Soft and filled with resilient polyester pile
- Prevents residents from digging fingernails into their palms
- Band keeps cushion from being dropped and separates fingers
- Covered with perspiration-absorbing fabric
- Machine launderable

ReOrder #	Description	Size	U/M
201140	Finger	Universal	6/pk
201150	Contracture		36/Cs

CONE GRIP



A cone-shaped anti-contracture device

- Elastic band prevents cone grip from being dropped
- Cloth cover is soft and perspiration absorbing
- Outer cover is removable for easy laundering
- Firm surface keeps hands open without stimulating finger contractions

ReOrder #	Description	Size	U/M
201060	Cone Grips	Universal	6/pk
201070			36/Cs


SENSORY INTEGRATION



914420



912422



914734



Sensory integration is a neurological process that organizes sensations from one's own body and from the environment. The five main senses are: Touch, Sound, Sight, Taste and Smell. To help Individuals that are deprived of their proprioception, Skil-Care developed a series of products to assist in improving sensory perception, eye-hand coordination and cognitive processing.

Typical products include:

- Cognitive training products provide an activity to support brain exercise.
- Weighted products provide proprioceptive input that helps increase body awareness, Improves attention span and concentration, and has calming benefits.

EASY-TO-USE SENSORY KEY



Creativity - Imagination & Arts



Fine Motor - Grasping, Eye-Hand Coordination



Visual - Seeing, Perceiving



Gross Motor - Large Muscle Movement



Thinking - Color Matching, Counting, Cause & Effect



Language - Reading, Writing



Social Skills - Relationships, Interaction



Tactile - Skin Awareness of Touch & Texture

GEL AQUARIUM



Sample Fish Inside Gel Aquarium




An exciting, stimulating pad that reduces anxiety and boredom

- Effective in reducing anxiety and boredom
- Serves as a tool in training for eye-hand coordination
- Helps improve finger strength and dexterity
- Can be placed on the lap to reduce anxiety

ReOrder #	Description	Size	U/M
912422	Gel Aquarium w/Four Fish	19"W x 12"D	Ea.

GEL MAZE



The gel-maze provides the resident the opportunity to move a marble through a maze pattern. This activity offers the sensory feeling of gel, the opportunity to project the path of the marble, and introduces the resident to linear geometry. All of these activities contribute to brain stimulation, improved eye-hand coordination, increased finger strength, dexterity, reduced boredom and anxiety.

ReOrder #	Description	Size	U/M
912425	Gel-Maze w/marbles - Green Gel	14"W x 14"D	Ea.
914517	Gel-Maze w/marbles - used with Light Box (Clear Background)	14"W x 14"D	Ea.

GEL SPIRAL MAZE





Original
Skil-Care
CORPORATION
Design



- Helps build eye-hand coordination skills
- Assists in developing finger strength and dexterity
- Offers a challenge in planning and anticipating outcomes from moving a marble through a maze

ReOrder #	Description	Diameter Size	U/M
914733B	Blue Gel Maze	17"	Ea.
914733R	Red Gel Maze		Ea.
914733Y	Yellow Gel Maze		Ea.
914733G	Green Gel Maze		Ea.

TIC-TAC-TOE GEL PAD



- Fun way to fine tune eye-hand coordination and increase attention span
- Allows two people to participate in an activity

ReOrder #	Description	Size	U/M
912435	Tic-Tac-Toe Gel Pad (Light Blue Gel)	10"W x 10"D	Ea.
914516	Tic-Tac-Toe Gel Pad - Clear for Light Box Use (no white background)	10"W x 10"D	Ea.

SPIRAL GEL PAD



Clear Background



White Background




Light Box sold separately (see pg 133) Item #914514

- Cushions tabletop to protect against injury caused by uncontrolled motor activity
- Improves finger strength and dexterity
- Provides sensory stimulation (tactile and visual)
- Clear vinyl background for use with light box

ReOrder #	Description	Size	U/M
912428	Spiral Gel Pad (White Background)	20"W x 22"D	Ea.
912441	Spiral Gel Pad used with Light Box (Clear Background)	15"W x 15"D	Ea.

LIGHT BOX



The light box features a 16" x 16" frosted light screen that is designed to be used either alone or in conjunction with the Skil-Care clear gel pads. The box also allows the user to trace line drawings, shapes and patterns as a daily activity. The touch responsive gel shapes encourage manipulation for tactile stimulation and hand/finger strengthening. Other expected outcomes are color recognition, improved eye hand coordination, reduction of boredom and increased attention span. **The light box is powered by 120VAC, with two 850 lumens light bulbs.**

ReOrder #	Description	Size	U/M
914514	Light Box	16"W x 16"D x 6"H	Ea.

Light Box Gel Pads

6 Spoke Snow Flake (15"x15")
 912447B - Blue
 912447G - Green
 912447C - Clear
 912447R - Red
 912447Y - Yellow



Quad Tree (15"x15")
 912442B - Blue
 912442G - Green
 912442C - Clear
 912442R - Red
 912442Y - Yellow




Other Light Box Gel Pads are listed in the Skil-Care Special Needs Catalog

PRESS N GLO LIGHT BOX



Model 912437






Model 912437

The unit can be placed in pressure sensitive mode to turn on and off intermittently when activated or it can be set to steady on mode to provide continuous illumination.


- Touch screen illuminates with LED lights.
- When combined with gel pads the unit provides a surface for tactile manipulation to help improve eye-hand coordination, reduce boredom and help to develop finger strength.
- Optional gel pads help develop learning skills and color and shape recognition.
- Unit operates under two power settings (pressure sensitive or continuous on).
- **Unit is portable and operates with "C" batteries (not included).**

ReOrder #	Description	Size	U/M
912436	Press N Glo Light Box	15"W x 15"D	Ea.
912437	Press N Glo Light Box w/Sensory Star	15"W x 15"D	Ea.

NUMBER-COLOR ASSOCIATION PAD



Model 914723



This pad can be used in several activities, namely:

- Color comparison and recognition
- Counting sequencing
- Number comparisons
- Eye-hand coordination
- Finger dexterity and strength exercises
- Supplied with a re-usable template for activity setup

ReOrder #	Description	Size	U/M
914720	Number-Color Association Pad	14"W x 8.5"D	Ea.
914723	Replacement Laminated Activity Sheets	14"W x 8.5"D	12/pk

BASEBALL GEL PAD



Original Skil-Care CORPORATION Design




It has been found that therapeutic touch is effective in decreasing boredom and anxiety. The baseball gel maze provides a familiar subject and media that will encourage users to move the marble throughout the gel. The process will also help in improving eye-hand coordination and finger strength and dexterity.

- When placed on the lap reduces anxiety
- Heat-sealed vinyl enclosure that is easy to clean

ReOrder #	Description	Size	U/M
914734	Baseball Gel Pad	15"W x 15"D	Ea.
914735	Baseball Gel Pad - used with Light Box (Clear Background)	15"W x 15"D	Ea.

HUMAN ANATOMY GEL PAD




Original Skil-Care CORPORATION Design

This gel pad acquaints the user with several body parts and provides tactile stimulation for developing eye-hand coordination and finger strength

- Serves as a useful tool in reading development
- Tactile stimulation decreases restlessness and helps improve eye-hand coordination

ReOrder #	Description	Size	U/M
914731	Human Anatomy Gel Pad	19"W x 12"D	Ea.

LETTER & NUMBER GEL PADS



A B C D E F G
H I J K L M
N O P Q R S
T U V W X Y Z

Model# 914630

E J L S H Y Q O
P C X A R A W O
D Z M T G N W

Model# 914631

a b c d e f g
h i j k l m
n o p q r s
t u v w x y z

Model# 914632


1 2 3 4 5
6 7 8 9 10
20 30 40 50 60
70 80 90 100

Model# 914633

Skil-Care's glitter gel pads are filled with a non-toxic gel and incorporates letters and numbers in the background. Research shows that mental effort can stimulate new brain cell growth. The best brain fitness activities are those that are new and challenging. These pads can be used to assist with visual perception, writing skills, letter and number recognition and tactile awareness.

ReOrder #	Description	Color	Size	U/M
914630	Alphabet Upper Case in Order	Blue	7"D x 16"W	Ea.
914632	Alphabet Lower Case in Order	Yellow	7"D x 16"W	Ea.
914631	Jumbled Alphabet with Numbers	Red	7"D x 16"W	Ea.
914633	Numbers in Order	Lt. Blue	7"D x 16"W	Ea.



GEL-FOAM PICTURE FRAME



There is no doubt that mental exercise nourishes the mind. One must continue to learn complex, new tasks in order to keep the brain fit and active. The brain has an amazing capacity for healing, whether the deterioration is due to aging, traumatic brain injury, or stroke. The picture frame is a heat sealed vinyl frame with a gel pocket in the front and a back pocket that allows one to insert any graphic that challenges or attracts the users interest. The product can be used as a learning tool or as a place to house familiar pictures.

ReOrder #	Description	Size	U/M
914722	Gel-Foam Picture Frame	17"W x 10"D	Ea.


WEIGHTED LAP PADS



- Provides proprioceptive input
- Helps reduce anxiety
 - Increases body awareness
 - Improves attention span
 - Vinyl cover wipes clean

ReOrder #	Description	Weight	U/M
914505	Rectangular Lap Pad Red Gel	3lbs	Ea.
914506	Rectangular Lap Pad Green Gel	5lbs	Ea.
914508	Rectangular Lap Pad Blue Gel	7lbs	Ea.
914682	12" Hexagon - Red Gel	3lbs	Ea.
914683	14" Hexagon - Green Gel	5lbs	Ea.
914684	16" Hexagon - Blue Gel	7lbs	Ea.
914689	Semi-Circle - Red Gel	3lbs	Ea.
914690	Half Circle - Green Gel	5lbs	Ea.
914691	Circle - Blue Gel	7lbs	Ea.
914692	Semi Circle Optional Cover	1lbs	Ea.

ACTIVITY AIDS






Activity Aids come in three styles: Apron, Vest and Overlay. These items are designed for patients with decreased cognitive function and the need for tactile stimulation.

- Helps redevelop fine motor skills
- Improves hand-eye coordination
- Helps improve dressing skills, dexterity and attention span
- Provides color and textures stimulation
- Attaches to tray, lap or worn as a vest depending on the model type.
- Pockets can hold small photos

ReOrder #	Description	Size	U/M
912461	Activity Apron	Small-Medium	Ea.
912460		Large	Ea.
912462	Activity Apron w/Hook & Loop Pocket	Large	Ea.
912464	Activity Vest	Small-Medium	Ea.
912465		Large	Ea.
912467	Activity Overlay	28"W x 11"D	Ea.


GEL WAVE PAD



- Recommended for improving finger dexterity and hand strength
- Reduces boredom and increases attention span
- Two tracks for use with two residents simultaneously

ReOrder #	Description	Size	U/M
912431	Gel-Wave Pad	21"W x 13"D	Ea.

SENSORY STIMULATION GEL PAD




Gel-filled with a kaleidoscope of colorful stars

- Soothing water-based gel calms agitated residents
- Bright stars move through gel for visual stimulation
- Tactile stimulation decreases restlessness and chair exits

ReOrder #	Description	Size	U/M
912420	Sensory Stimulation Gel Pad	16"W x 7"D	Ea.

GEL-TOP SENSORY TRAY



912410



914504

Bright stars in water-based gel calms agitated residents

- Convenient working surface for feeding & activities
- Water-based gel surface provides tactile stimulation
- Moving stars create visual appeal
- Decreases chair exits by reducing restlessness
- Gel surface relieves pressure on arms & elbows
- Covered in durable, wipe-clean vinyl

For other trays see page 46

ReOrder #	Description	Size	U/M
912410	Gel Top Sensory Tray, Black Vinyl, Lift-Away	Fits 16"-18" Wheelchair	Ea.
914504	Gel Top Sensory Tray, White Vinyl, Lift-Away	Fits 16"-18" Wheelchair	Ea.
912415	Gel Top Sensory Tray, Black Vinyl, Hook & Loop Strap	Fits 16"-18" Wheelchairs	Ea.

WAVY ACTIVITY TRAY

This therapy product helps develop motor skills while providing sensory sensations.

- Provides an entertaining way to improve motor skills
- Strengthens eye-hand coordination, visual perception and sensory stimulation
- All items attach with Hook & Loop and are easily removable and replacable for a puzzle type activity.



ReOrder #	Description	Size	U/M
914671	Wavy Activity Tray	22.5"W x 13.5"D	Ea.

ADJUSTABLE WEIGHTED LAP BLANKET



Gel pad partially removed from pocket



Weighted lap blanket provides therapeutic benefits.

- Provides calming comfort
- Provides gentle pressure similar to swaddling
- Helps with sleep issues and restlessness
- Reduces stress
- 5lbs. gel pack consists of eight 10oz gel weights
- There are no metal parts in this blanket and can be used in any hospital setting
- Individual pockets allow for customization of weight distribution
- Weights designed for children & adults
- Ultra-soft non-allergenic micro-fiber can be laundered (with gel weights removed)

ReOrder #	Description	Size	U/M
914000	Adjustable Weighted Blanket - 36 Pockets	54"L x 36"W	Ea.
914001	Adjustable Weighted Blanket w/ 5lbs of Weight - 36 Pockets	54"L x 36"W	Set
914006	Adjustable Weighted Blanket - 49 Pockets	63"L x 42"W	Ea.
914007	Adjustable Weighted Blanket w/5lbs of Weight - 49 Pockets	63"L x 42"W	Set
914008	Gel Weights, 8 pieces at 10oz each	5lbs	8 pcs.

VIBRATOR MITT




The vibrator mitt provides a soft vibrating massage for relaxation and sensory stimulation. This is especially helpful for listless residents prior to preparing for daily activities. The product offers a choice of two sensory materials; an ultra-soft fabric on one side and synthetic sheepskin on the other. Since the unit is pressure activated no on-off switch is needed. This product is launderable once the vibrator is removed.

Re Order #	Description	Size	U/M
914560	Vibrator Mitt	Universal	Ea.

SENSORY PILLOWS



Original Skil-Care Corporation Design



Provides gentle sensory stimulation with soothing textures.

- Reduces anxiety and boredom
- Multi-texture covered foam
- Use it as a pillow, place it on your lap or simply hold it
- Use at home or in clinical settings to provide a sensory experience
- Removable cover for easy cleaning
- Improves finger strength and dexterity

ReOrder #	Description	Size	U/M
914582	Small	14"W x 13.5"D x 5"H	Ea.
914584	Large	20"W x 16"D x 5"H	Ea.

Abduction/Contraction

- Abductor/Contracture Cushion 72
- Abduction Wedge 54
- E-Z Transfer Slider Pommel 24
- Gull Wing Leg Abductor 54
- Hip/Leg Abduction Pillow 71
- Leg Abductor Cushion 25,32
- Slide Guard Cushion and Wedge 25,32

Alarms 3-16

- Accessories 98-99
- BedPro Systems and Sensor Pads 91-92
- ChairPro Systems and Sensor Pads 87-90
- FloorPro Systems and Sensor Pads 96
- Magnetic Alarms 97
- MultiPro Alarms and Pads 92-95

Ambulation [see Gait Belts]

Anti-Rotational Gripper 5,6,8,9

Aprons

- Activity Apron 60,135
- Lap Blanket w/Hand Warmer 59
- Modesty Apron 60
- Smokers Apron 84,101

Aprons, Smokers 84,101

Arm Sling Pouch 12,123

Arm Supports

- Armrest Pads 43,58
- Armrest Wheelchair Cushion with Pouch .. 58
- Extended Full and Half Armrest Pads 58
- EZ On Lateral Support 51
- Lateral Body Support 50
- Lateral Stabilizer Armrest Bolster 51
- Lateral Stabilizer Arm Trough 51
- Mobile Arm Support 58
- Pillow Prop 70

Back Comfort & Support

- Adjustable Thickness Lumbar Support 56
- Air Lift Back Cushion 44
- Backrest Seat Combo 43
- Full Round Lumbar Support 56
- Lateral/Lumbar Support 49,56
- Reclining Wheelchair Backrest 43

Back Comfort & Support (cont.)

- Seat & Back Air Cushion 45
- Wheelchair Backrest Pad 45
- Wheelchair Cozy Seat 45

Bags/ Holders

- Armrest Wheelchair Cushion with Pouch .. 58
- CarPack 62
- Chart Holder 61
- ChairPack 62
- Footrest Bag for Wheelchairs 61
- Oxygen Cylinder Holder 61
- Universal Storage Bags 61
- Urinary Drainage Bag Holder 125
- Wheelchair/Walker Handy Bag 62

Bariatric Aids 109-115

- Cushions 36
- Gait Belts 113
- Geri-Sleeves 112
- Heel Floats 110
- Heels-Off 110

Leg Support

- Econo Footrest Extender, One-Piece 111
- Econo Footrest Extender, Two-Piece 53
- Knee Elevator 112
- Leg Pads 111
- Swing-Away Foot Support 111
- SoftTop Tray 113

Bath Mitt 123

Bed Ladder 128

Bed-Foot Support 72

Bed Positioning 65-73

- Abductor/Contractor Cushion 72
- Adjustable Anti-Reflux Wedge 72
- Bed Foot Support 72
- Bed Glide System 68
- Bed Support System w/Two Attached
- 30° Bolsters 67
- Bed System w/Slider Sheet
- & Two 16" or 24" Wedges 68
- Body Aligner 71
- Bolsters, Ultra Soft 73

Bed Positioning (cont.)

- Elevating Bed Wedge 73
- Hip/Leg Abduction Pillow 71
- In-Bed Resident Positioning System 66
- Infinity Pillow 38
- Knee Elevator 70
- Leg/Thigh Elevator 70
- Pillow Prop 70
- Positioning Wedge, 45° 66
- Positioning Wedge, 30° 66
- Roll Control Bolsters 67
- Roll Sheet 69
- Slide Sheet, Four Handle 69
- Slider Sheet with Continuous Handles 69
- Super Soft Head Pillow 38
- Universal Glide Sheet 68
- 30° Full Body Bed System w/4 Bolsters 67

Bed Rail Pads 74-77

- Bed Rail Wedge Pads 75
- Bed/Wall Protector 78
- Cushion for Foot and Headboard 78
- Cushion Top Vinyl Bed Rail Pads 76
- Gap Guard 76
- Half Rail Vinyl Bed Pads 76
- Split-Rail Vinyl Bed Rail Pads 75
- Synthetic Sheepskin Bed Rail Pads 77
- Thru-View Vinyl Bed Rail Pads 77
- Vinyl Bed Rail Pads (classic) 75

Belts

- Adjusta-Loop Cushion Belts 118
- Cushion Belt 117
- Cushion Slider Belt 117
- Econo Belt 117
- Economy Gait Belt 107
- E-Z Release Belt Sensor 88,95
- Gait Belts (see Gait Belts)

Non-Restrictive

- Resident-Release Econo Seat Belt 64
- Resident-Release Soft Belts 63
- Resident-Release Nylon Belts 63
- Resident-Release Slider Belts 64
- Resident-Release Slide Control Belts ... 64

Belts (cont.)

- Resident-Release Universal Belts 64
- Pelvic Holder 117
- Posture Support 48
- Resident-Release Slide Control Belt 64
- Roll Belt 118
- Seat Belt Alarm Sensor 88,95
- Shower Chair Safety Belt 101
- Transfer Belts w/Adjustable Handles 105
- Walker Belt 105

Bibs 123

BiFold Sensor Pad 88,93

Blankets

- Activity Apron 60,135
- Adjustable Weighted Lap Blanket 137
- Lap Blanket with Hand Warmer 59
- Shower Poncho 59
- Rain Cape for Wheelchair 59
- Modesty Apron 60

Bolsters

- Lateral Stabilizer Armrest Bolsters 51
- Lateral Stabilizer Armrest Trough 51
- Leg Bolsters 55
- Roll Control Bolsters 67
- 30° Full Body Bed System w/4 Bolsters 67
- Ultra-Soft Bolster 73

Catheter Leg Strap 125

Calf Pad Cover 55

Chart Holder 61

Commode Seat Cushion 19

Commode Seat Cushion w/Alarm 90

Cushions 14-41

Multi-Purpose

Gel-Foam

- Econo-Gel Cushion 16
- Gel-Foam Cushion (classic) 17
- Gel-Foam Toilet Seat Cushion 19
- Gel Lift Cushion 19
- Gel-Foam X-Cushion 18
- Starry Night Gel-Foam Cushion 18

Cushions (cont.)

Sittin Pretty Gel-Foam Cushion..... 17
 Super Soft Econo-Gel Cushion 16
Foam Only
 Comfort Foam Wheelchair Cushion 19
 Convoluted Foam Cushion..... 22
 EZ Dry Foam Cushion 21
 Gel-Infused Visco Cushion Topper 20
 Pressure-Check Foam Cushion 20
 Pressure-Check Light/Heavy Cushion .. 20
 Solid Foam Cushion w/Sheepskin..... 20
 Wedged Foam Wheelchair Cushion..... 21
Air
 Air Lift Seat Cushion..... 23
 Air Lock Seat Cushion 23
 Cushion Air 22
 Foam Air Cushion 22,30
Slide Control
 Anti-Thrust Cushion..... 24
 Anti-Thrust Visco Foam Cushion..... 26
 Contour Cushion..... 27
 Contour-Plus Cushion 27
 E-Z Transfer Slider Pommel 24
 Leg Abductor Cushion 25,32
 Position Plus - Leveled & Wedged ... 26,29
 Slide Guard Level and
 Wedged Cushion 25,32
 Stability Plus Gel-Foam Cushion..... 26
 Visco Conform Cushion..... 27
Thin-Line™
 Econo-Gel Cushion 29
 Foam Air Cushion 22,30
 Gel-Foam Cushion 28
Thin-Line™
 Gel-Foam X-Cushion..... 28
 Gel-Infused Visco Cushion Topper ... 20,28
 Position Plus Cushion..... 26,29
 Ultra Cushion 29,31
 Visco Conform Cushion..... 27,30
 X-Gel Overlay 30
Firm Foundation
 Foam Wheelchair Wedged Cushion..... 32
 J-Hook Drop Seat..... 34

Cushions (cont.)

Lateral Positioning Cushion..... 33
 Leg Abductor Cushion 25,32
 Slide Guard Leveled and
 Wedged Cushion 25,32
 Solid Seat Platform..... 34
 Tri-Foam Gel-Infused Visco Cushion 33,36
 Ultra Cushion Leveled 29,31
 Ultra Cushion Wedged 31
 Wheelchair Wedged Base..... 33
 Wheelchair Leveling Pad..... 34
Coccyx Cutout
 Foam Air Cushion 35
 Gel-Foam Cushion (classic) 35
 Gel-Foam X-Cushion..... 35
 Visco Conform Cushion..... 35
Bariatric
 Gel-Foam Cushion 36
 Position Plus Customized Cushion 36
 Tri-Foam Gel-Infused Visco Cushion ... 36
Reminder
 Air Lift Cushion 41
 Air Lock Cushion 41
 Comfort Foam Wheelchair Cushion 40
 Econo-Gel Cushion 39
 EZ Dry Foam Cushion..... 40
 Gel-Infused Visco Cushion Topper 39
 Super Soft Econo-Gel Cushion 39
 Wedged Foam Wheelchair Cushion..... 40
 X-Gel Overlay 41
Cushion Covers..... 37
Cushion Grips..... 129
Elbow Support
 Elbow Pad Protectors 12
 Gel-Infused Visco Foam Elbow Protector . 12
 Heel/Elbow Protector..... 12
 Padded Geri-Sleeves 13,112
 Pillow Prop 70
 Pouch Arm Sling 12,123
Elevating Bed Wedge..... 73
Exercise & Recreation
 Aqua Weights 127

Exercise & Recreation (cont.)

Aqua Air Exerciser 128
 Arctic Thermal Sleeve..... 128
 Bed Ladder 128
 Cone Grip 129
 Cushion Grip..... 129
 Ex-Box 127
 Finger Contractor Cushion 129
Fall Mats
 EZ Landing II 102
 Roll-On Bedside Fall Mat - 1.5" Beveled . 103
 Safe & Sound Fall Mat..... 102
 Safeside Fall Mat - 1" Beveled 103
 Soft-Fall Bedside Mat - 2" & 4" Beveled .. 103
 Tri-Fold Bedside Fall Mat - 1" Beveled ... 103
Finger Contractor Cushion..... 129
Fire Blanket..... 100
Floor Mats [see Fall Mats]
Foot Supports
 Abduction Wedge 54
 Calf Pad Cover 55
 Foot Cradle 52
 Foot/Heel Elevator 3
 Hip/Leg Abduction Pillow 71
 Geri-Chair Foot Cradle 85
 Geri-Chair Leg Positioner 85
 Gull Wing Leg Abductor 54
 Knee Elevator 70,112
 Leg Bolsters 55
 Leg/Thigh Foam Elevator 70
 One-Piece Econo-Footrest Extender 53
 Skin-Guard Leg Protector..... 55
 Side-Kick..... 53
 Swing-Away Foot Support 52,111
 Synthetic Sheepskin Relief Slippers..... 55
 Two-Piece Footrest Extender 53
 Ultra-Soft Bolsters..... 73
 Wheelchair Leg Pad 54
 Wheelchair Legrest with Padded Sides... 54
 Wheelchair Footrest Extender
 w/Leg Separation..... 52

Gait Belts

Bariatric..... 106,113
 Cotton, Standard..... 106,113
 Economy..... 107
 Gait Belt Handles..... 107
 Geo-Pattern 107
 Infection Control..... 106,113
 Transfer Belts..... 105
 Walker Belt..... 105
Gap Guard..... 76
Gel Pack, Replacement..... 6
Geri-Chair Accessories
 Adjustable Lateral Support 84
 Alarm System 84
 Cozy Seat 85
 Foot Cradle 85
 Gel Overlay..... 85
 Leg Positioner..... 85
 Side Wings..... 84
 Smokers Apron 84
Geri-Sleeves
 Bariatric..... 112
 Geri-Sleeve, Regular 13
Grip, Cushion..... 129
Gripper, Anti-Rotational..... 8
Head Control/Rests
 Adjustable Head Positioner w/Gel..... 57
 Anti-Reflux Wedge..... 72
 Elevating Bed Wedge 73
 Standard Headrest 57
 Triangular Head Positioner 57
Heel Protectors
 Bilateral Anti-Rotation Gripper 8
 Covered Foam Heel Cushion 11
 Econo Foam Heel Protector 11
 Econo Sheepskin Heel Protector 11
 Foam Pressure Relieving Heel Protector 9
 Foot Drop Boot 8
 Foot/Heel Elevator 3
 Gel-Foam Heel Cushion 10
 Heels-Off..... 3,110
 Heel/Elbow Protector..... 12

Heel Protectors (cont.)

Heel Floats

Adjustable 7
 Adjustable Walker 7
 Air Boot 8
 Bariatric 5-9,110

Heel Floats

Heel-Float II 5
 Original 5
 Plus 6
 Walker 7
 Wipe-Clean 6
 Pillow Prop 70
 Soft Boot 10
 Stay-On Heel Protector 12
 Super Soft Heel Protector 9
 Triple-Ply Heel Protector 11
 Tri-Function Heels Off-Loader 3
 Ultra-Soft Fiber-Fill Heel Cushion 10

Hip Protection 100

In-Bed Resident Positioning 66

Knee Elevator 70,98

Lateral Supports

Adjustable Lateral Support 48,84
 E-Z On Lateral Support 51
 Geri-Chair Side Wings 84
 Lateral Body Support 50
 Lateral/Lumbar Support 49,56
 Lateral Stabilizer Armrest Bolster 51
 Lateral Stabilizer Arm Trough 51
 Lateral Support Orthosis 44
 Mobile Arm Support 58
 One-Way Glide System 49
 Snug Supports 50
 Universal Chair Side Wings 50
 Wheelchair Side Wings 48
 Wheelchair Flip Tray 47
 Wheelchair Posture Support 48

Leg/Thigh Foam Elevator 70

Leg Positioning and Support

(see Foot Support)

Light Box 133

Limb Holders (see Restraints)

Lumbar Support

Adjustable Thickness Lumbar Support 56
 Full Round Lumbar Support 56
 Lateral/Lumbar Support 49,56

Non-Restrictive 63,64

Non-Slip Matting 101

Mattresses

Bariatric 81
 Premium Pressure Check 80
 Premium Gel-Infused Visco Foam 81
 Psychiatric 81
 Replacement Covers 82
 Safety Strap 82
 Synthetic Sheepskin Pads 82

Oxygen Cylinder Holder 61

Patient Care

Bath Mitt 123
 Bibs 123
 Pouch Arm Sling 12,123
 Trach Tube Tie/Holder 123
 Telemetry Unit Pouch 123

Pillows

Infinity Pillow 38
 Pillow Prop 70
 Super Soft Head Pillow 38

ProSeries (see Alarms)

Rain Cape 59

Reclining Wheelchair Backrest 43

Reminder Cushions

Air Lift Cushion 41
 Air Lock Cushion 41
 Comfort Foam Wheelchair Cushion 40
 Econo-Gel Cushion 39
 EZ Dry Foam Cushion 40
 Gel-Infused Visco Cushion Topper 39
 Super Soft Econo-Gel Cushion 39
 Wedged Foam Wheelchair Cushion 40
 X-Gel Overlay 41

Restraints, Physical

Belts

Adjusta-Loop Cushion Belt 84

Restraints, Physical (cont.)

Belts

Cushion Belt 117
 Cushion Slider Belt 117
 Econo Belt 117
 Pelvic Holder 117
 Roll Belt 118

Jackets & Vests

Safety Vest 119
 Sleeper Jacket 119
 Sleeved Jacket 119

Limb Holders

Economy 118
 Sheepskin 118

Padded Mitts

E-Z View Padded 121
 Finger Separator 121
 Padded 120
 Padded Plus 120
 Rigid Palm 120

Pelvic Control

Cushion Slider Belt 117
 Pelvic Holder 117

Roll Control, Bed Systems 66,68

Roll Control Bolster 67

Safety

Floor Fall Mats (see Fall Mats)

Resident

Hip-Ease 100
 Shower Chair Safety Belt 101
 Slip-Grip Non-Slip Matting 101
 Smokers Apron 101
 Stop Strip 100

Patient Handling

Roll Sheet 69
 Slide Sheet, Four Handle 69
 Slider Sheet with Continuous Handles 69

Sensory Integration

Activity Aprons 135
 Activity Vests 135
 Activity Overlay 135
 Adjustable Weighted Lap Blanket 137

Sensory Integration (cont.)

Baseball Gel Pad 134
 Gel Spiral Maze 132
 Gel Aquarium 131
 Gel-Foam Picture Frame 135
 Gel-Maze 131
 Gel-Wave Pad 136
 Gel-Top Sensory Tray 46,136
 Human Anatomy Gel Pad 134
 Letter & Number Gel Pads 134
 Light Box 133
 Light Box Gel Pads 133
 Number/Color Association Pad 133
 Press N Glo Light Box 133
 Sensory Pillows 137
 Sensory Stimulation Gel Pad 136
 Spiral Gel Pad 132
 Tic-Tac-Toe Gel Pad 132
 Vibrator Mitt 137
 Wavy Activity Tray 136
 Weighted Lap Pads 135

Sheepskin Coverings 59

Shower Poncho 59

Skin Protection

Armrest Pads 58
 Bath Mitt 123
 Calf Pad Cover 55
 Foot/Heel Elevator 81
 Geri-Sleeves 13
 Heel/Elbow Protectors 84
 Heel Protectors (see Heel Protectors)
 Leg Bolsters 55
 Pouch Arm Sling 12,123
 Skin-Guard Leg Protector 55
 Synthetic Sheepskin Pads 82
 Visco Foam Elbow Pad 12

Slip Grip Non-Slip Mat 101

Slippers, Synthetic Sheepskin 55

Slings 12,123

Smokers Apron 84,101

Snug Support 50

Stabilizers (see Heel Protector)

Stop Strip	100
Synthetic Sheepskin Pads	82
Telemetry Unit Pouch	123
Toilet Seat Cushion [see Cushions]	
Toilet Seat Cushion w/Alarm [see Commode Seat Cushion w/Alarm]	
Trach Tube Holder	123
Transfer Belts	105-107
Transfer Slings	
Super Slings	108
TLC Pad.....	108
Trays and Lap Cushions	
Flip Tray	47
Gel Top Sensory Stimulation Tray	46,136
Lap Top Cushion.....	47
Lap Cushion, Lift-Off.....	47
Lift Away Wheelchair Tray	46
Wheelchair Lap tray w/straps	46
Urinary	
Urinary Drainage Bag Holder	125
Catheter Leg Strap	125
Vibrator Mitt	137
Warmth [see Blanket]	
Wedges	
Adjustable Anti-Reflux	72
Body Aligner.....	71
Elevating	73
In-Bed Resident Positioning System	66
Elevating Arm Support.....	71
Roll-Control Bolsters.....	67
30° Full Body Bed Support System	67
30° Positioning	66
Bed Support System w/Two Attached 30° Bolsters	67
30° Bed Bolster System w/Slider Sheet & Two 16" or 24" Wedges.....	68
45° Positioning	66
Weights	
Aqua	127
Air.....	128
Weighted Lap Pads	135

Wheelchair Accessories

Aprons	
Activity Aprons	135
Lap Blanket w/Hand Warmer.....	59
Modesty Apron.....	60
Armrest	
Armrest Pads	43,58
Armrest Wheelchair Cushion with Storage Pouch.....	58
Mobile Arm Support	58
Sheepskin Coverings.....	59
Backrest	
Air Lift Backrest Cushion	44
Backrest Seat Combo	43
Reclining Wheelchair Backrest	43
Seat and Back Air Cushion.....	45
Wheelchair Backrest Pad	45
Wheelchair Cozy Seat	45
Bags	
Car Pack	62
Chair Pack	62
Chart Holder	61
Oxygen Cylinder Holder	61
Wheelchair Footrest	61
Wheelchair Walker Handy	62
Universal Storage	61
Belts	
Resident-Release Nylon.....	63
Resident-Release Slider.....	64
Resident-Release Slide Control	64
Resident-Release Universal Safety.....	64
Headrest	
Adjustable Head Positioner w/Gel.....	57
Standard	57
Triangular Head Positioner	57
Lateral Support	
Adjustable Lateral Support	48,84
Chair Glide System.....	49
E-Z on Lateral Support	51
Lateral/Lumbar Support.....	49,56
Lateral Body Support.....	50
Lateral Stabilizer Armrest Bolster	51

Wheelchair Accessories (cont.)

Lateral Stabilizer Arm Trough	51
One-Way Glide System	49
Snug Support.....	50
Universal Chair Side Wings.....	50
Wheelchair Posture Support.....	48
Wheelchair Side Wings	48
Leg Support	
Abduction Wedge	54
Calf Pad	55
Foot Cradle	52
Gull Wing Leg Abductor.....	54
Leg Bolsters.....	55
One-Piece Footrest Extender	53
Side-Kick.....	53
Skin-Guard Leg Protector.....	55
Swing-Away Foot Support.....	52,111
Synthetic Sheepskin Slippers	55
Two-Piece Footrest Extender	53
Universal Chair Side Wings.....	50
Wheelchair Footrest Extender w/Leg Separation.....	52
Wheelchair Leg Pad	54
Wheelchair Legrest w/Padded Sides	54
Lumbar Support	
Adjustable Thickness.....	56
Full Round	56
Lateral/Lumbar.....	49,56
Rain Cape & Poncho	
Rain Cape.....	59
Shower Poncho	59
Wheelchair Posture Support	48
Wheelchair Cushions [see Cushions]	
Wound Care [see Skin Protection or Heel Protectors]	

MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#
101080	125	252250	107	306130	121	503034	5	503180	8	553061	73	556155	67,69	703006	45
101090	125	252252	107	307010	47	503035	5	503181	8	554010	66	556157	68	703010	58,59
102010	125	253010	108	307012	47	503036	5,110	503182	8	554020	66	556158	68,69	703015	58
102011	125	253020	108	307014	47	503038	8	503185	3	554022	66	558015	80	703016	58
102012	125	253030	108	307020	47	503039	8	503186	3	554025	66	558020	80	703018	58
102015	125	301101	119	307030	47	503040	11	503187	3	554026	66	558025	80	703019	58
102020	125	301102	119	307032	47	503044	6	503188	3	554028	66	558051	80	703020	59
201020	129	301103	119	307034	47	503045	6	503189	3	554030	66	558052	80	703040	59
201030	129	301111	119	307040	47	503046	6,110	503198	3	554040	72	558053	80	703050	59
201060	129	301112	119	307042	47	503047	6	503199	3	554110	73	558070	82	703060	55
201070	129	301113	119	307044	47	503048	6	503203	12	554120	73	558071	82	703061	55
201140	129	301121	119	307050	47	503049	6,110	503204	12	555012	67,108	558072	82	703063	45
201150	129	301122	119	307051	47	503050	12	503205	12	555015	66,108	558089	82	703066	56
201170	129	301123	119	401010	75	503060	12	503213	12	555018	66,108	558091	82	703068	55
201175	129	301200	117	401015	77	503070	11	503214	12	555025	66	558092	82	703070	54
251011	105	301201	117	401017	77	503077	11	503215	12	555035	66	558093	82	703071	54
251022	105	301202	118	401040	75	503080	11	503230	12	555036	66	558094	82	703072	54,111
252000	107	301205	118	401045	77	503087	11	503240	12	555040	70	558100	81	703073	55
252001	107	301250	63	401047	77	503088	10	503310	12	555041	70	558102	81	703075	54
252002	107	301255	63	401050	75	503092	128	503320	12	555043	70	558104	81	703077	54
252003	107	301270	63	401080	75	503095	5,6	503350	13	555045	70	558123	81	703080	59
252004	107	301275	63	401090	76	503097	5,6	503352	13	555046	70,112	558125	81	703082	34
252005	107	301300	118	401100	75	503098	5,6,8,9	503355	13	555055	72	558127	81	703101	43
252006	107	303210	117	401140	75	503099	5,6,57,110	503356	13	555056	72	558133	81	703102	43
252007	107	303311	117	401150	75	503100	3	503357	13,112	555060	72	558135	81	703103	43
252010	106	303312	117	401160	78	503101	3	503358	13,112	556010	67	558137	81	703105	43
252011	106	303321	117	401165	78	503134	5	503360	13,112	556012	67	610111	48	703106	43
252015	106	303322	117	401170	78	503135	5	503361	13,112	556020	67	610112	48	703107	43
252020	107	304203	119	401210	76	503136	5	503362	13	556021	67	610115	48	703111	43
252021	107	304204	119	401220	76	503140	3	503363	13	556025	67	610116	48	703112	43
252023	107	304205	119	401230	76	503141	3	503364	13	556030	67	610220	64	703113	43
252024	107	304213	119	401235	75	503142	7	503365	13	556038	67	610230	64	703115	57
252026	107	304214	119	401240	75	503143	7	503366	13,112	556065	68	610310	64	703116	57
252027	107	304215	119	401310	76	503144	7	503367	13,112	556066	68	610315	64	703117	57
252029	107	305103	119	401320	76	503147	7	503368	13,112	556067	68	620010	64	703121	57
252030	107	305104	119	402010	77	503148	7	503369	13,112	556068	69	620020	64	703130	43,58
252050	106	305105	119	501010	82	503149	7	503390	13	556069	69	701010	63	703131	43,58
252051	106	305113	119	501011	82	503150	3,110	503392	13	556070	69	701011	63	703199	43
252055	106	305114	119	501020	82	503151	3	503394	13	556080	49	701015	63	703210	53
252060	106	305115	119	501021	82	503152	7	503400	9	556081	49	701020	117	703220	53
252061	106	306010	118	501022	82	503153	7	503410	9	556084	49	701030	63	703230	53
252070	106	306020	118	501050	82	503154	7	503420	9	556085	49	701031	63	703240	53
252071	106	306040	118	501051	82	503160	10	503430	5,6,8,9	556086	49	701035	63	703250	53
252072	106,113	306041	118	501090	82	503165	3	503450	9	556087	68	701055	64	703251	53
252073	106,113	306044	118	501210	82	503165	110	503510	8	556088	68	702210	44	703252	53
252074	106,113	306045	118	503005	11	503166	3	553010	73	556089	68	702220	44	703253	53
252075	106	306048	118	503006	11	503167	3,110	553011	73	556150	68	702230	44	703256	53
252200	105	306110	120	503015	11	503170	8	553030	73	556151	68	703001	85	703257	53
252201	105	306115	120	503016	11	503171	8	553031	73	556152	68,69	703002	85	703258	53
252202	105	306120	120	503025	10	503172	8	553050	73	556153	68	703003	85	703259	53
252203	105	306125	121	503030	10	503174	8	553051	73	556154	69	703005	45	703262	52

MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#
703264	52	705019	46	706362	24	751050	17	751350	32	753272	19,40	754528	17,35	757109	24
703267	52	705020	46	706363	24	751051	17	751370	25	753280	19	754650	115	757110	24
703269	52	705021	46	706364	24	751052	17	751420	26,29	753285	19,40	754655	115	757112	24
703271	53	705022	46	706365	56	751100	19	751435	26,29	753290	19,40	754660	115	757115	24
703272	53	705024	46,113	707010	61	751101	19	751450	26,29	753295	19	754665	115	757117	24
703273	53	705025	46	707026	61	751148	16	751465	26	753400	21,40	754750	33	757120	24
703274	53	705026	46	707027	61	751149	16	751480	26	753402	21,40	754755	33	757122	24
703280	53	705027	46	707028	61	751150	16,39	751495	26	753404	21,40	754757	33	757125	24
703281	53	705031	47	707029	61	751151	16,39	751540	22	753406	21,40	754810	33,36,115	757127	24
703282	53	705032	47	707030	60	751153	16,39	751545	22	753408	21,40	754812	33,36,115	757130	24
703283	53	705091	46	707035	62	751154	16,39	751550	22	753410	21,40	754814	36,115	758010	25,32
703286	53	706005	51	707038	62	751155	16,39	751560	22	754010	21,40	754816	36,115	758015	25,32
703288	53	706010	50	707052	19,40	751156	16	751562	22	754016	21,40	754818	36,115	758020	25,32
703291	53	706011	50	708030	127	751157	16,39	751564	22	754020	21,40	754821	36,115	758025	25,32
703292	53	706012	50	710010	71	751162	16	751566	22	754021	21,40	754823	36,115	781014	37
703293	53	706015	50	710015	71	751163	16	751567	22	754032	21	754825	36,115	781015	37
703294	53,111	706016	50	751010	17	751164	16	751568	22	754035	21	754827	36,115	781016	37
703295	53	706017	50	751011	17	751166	16,39	751569	22	754036	21	754829	36,115	781017	37
703296	53,111	706020	50	751012	17	751167	16,39	751580	8,22,23,44	754037	21	754833	36,115	781018	37
703297	53,111	706021	50	751013	17	751168	16,39	751610	27	754039	21,40	754835	36,115	781019	37
703300	45	706050	48	751014	17	751170	29	751615	27	754040	21,40	754837	36,115	781028	37
703350	55	706055	48,84	751018	17	751171	29	751620	27	754041	21,40	754839	36,115	781029	37
703351	55	706059	48,84	751019	17	751172	29	751625	27	754043	21,40	754841	36,115	781030	19,37
703352	55	706110	51	751020	17	751173	29	751630	27	754044	21	754843	36,115	781031	37
703353	55	706111	51	751021	17	751174	29	751635	27	754048	21,40	754845	36,115	781032	37
703410	52	706112	51	751024	17,35	751175	29	751636	27	754050	21,32	754847	36,115	781033	37
703416	52	706114	51	751025	17,35	751176	29	751637	27	754051	21,32	754849	36,115	781034	37
703420	85	706115	51	751026	18,28	751177	29	751675	27	754054	21	754851	36,115	781035	37
703430	85	706116	51	751027	18,28	751178	29	751680	27	754110	31	756010	33	781036	37
703455	53	706120	51	751028	18	751179	29	753150	27,30	754114	31	756210	23,44	781037	37
703460	54	706121	51	751029	18	751180	16	753155	27,30,35	754120	31	756211	23	781038	37
703470	52	706122	51	751030	18,28	751181	16	753157	27	754124	31	756215	23,41,44	781039	37
703471	52	706124	51	751031	18,28	751188	19,40	753160	26	754130	31	756219	23,41,44	781040	37
703472	52	706125	51	751032	18,28	751189	19,40	753178	20	754132	31	756225	23,44	781041	37
703475	52,111	706126	51	751033	18,28	751196	18	753179	20	754134	31	756226	23	781042	37
703476	52,111	706140	51	751034	18,28,35	751197	18	753180	20	754160	31	756230	23,41,44	781043	37
703477	52,111	706141	51	751035	18,28,35	751198	18	753181	20	754164	31	756234	23,41,44	781044	37
704310	34	706142	51	751036	18,28,35	751200	17	753182	20	754170	31	756235	23	781045	37
704311	34	706210	58	751037	18	751201	17	753183	20	754174	31	756236	23	781046	37
704312	34	706215	58	751038	18,28	751202	17	753184	20	754180	31	756237	23,41	781047	37
704315	34	706220	58	751039	18,28	751210	25,32	753186	20	754184	31	756239	23,41	781048	37
704320	34	706225	58	751040	18,28	751230	25	753210	19,40	754310	31	756240	23	781049	37
704321	34	706309	50	751041	18,28	751260	25	753220	20	754312	31	756241	23	781054	37
704325	34	706310	50	751042	18,28	751265	25	753230	19,40	754320	31	756245	23,41	781055	37
704330	34	706319	50,84	751043	18,28	751269	25	753232	19,40	754322	31	756250	23	781056	37
704331	34	706320	48	751044	18	751280	32	753240	19	754324	31	756251	23	781057	37
704348	34	706322	48	751045	18	751310	32	753250	19,40	754330	31	756255	23,41	781058	37
705014	46	706350	49,56	751046	26	751314	32	753251	19,40	754332	31	757010	34	781059	37
705015	46	706351	49,56	751047	26	751316	25	753260	20	754334	31	757015	34	781065	30,41
705016	46	706352	49,56	751048	26	751320	25	753270	19,40	754339	29,31	757017	34	781066	30,41
705018	46	706359	49,56	751049	18	751324	25	753271	40	754340	29,31	757108	24	901130	123

MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#	MODEL#	PG#
902011	12,123	909293	96	909389	89,94	911460	100	914001	137	914452	66	915065	114
902012	12,123	909294	96	909390	99	911470	100	914006	137	914471	33	915066	114
904110	123	909310	92	909391	94	911530	103	914007	137	914472	33	915067	114
904150	123	909311	92	909392	95	911535	103	914008	137	914473	33	915068	114
904151	123	909314	92,93	909393	95	911536	103	914117	95	914481	18,28	915069	114
905010	70	909316	92,93	909394	95	911539	102,103	914118	95	914482	18,28	915105	36
905015	70	909317	92,93	909395	88	911540	103	914119	95	914483	18,28	915110	36,114
905020	38	909319	92,93	909397	19,90	911545	103	914123	108	914504	46,136	915115	36
905022	38	909320	91	909398	90,99	911546	103	914154	128	914505	135	915120	36,114
905023	38	909321	91	909399	99	911547	103	914161	61	914506	135	915125	36
905025	38	909322	91	909401	91,93	911548	103	914210	50	914508	135	915130	36,114
905030	71	909323	91	909402	91,93	911549	103	914220	51	914514	133	915131	36
905031	71	909325	91,93	909403	91,93,96	911570	102	914230	58	914516	132	915133	36,114
905032	71	909327	91,93	909404	91,93	911572	103	914235	58	914517	131	915134	36,114
906010	101	909328	91,93	909420	88	911575	102	914247	22,30,35	914539	127	915135	36,114
906011	101	909330	91,93	909421	88	911577	103	914248	22,30,35	914560	137	915137	36,114
906015	84,101	909331	91,93	909425	88,93	911585	102	914249	22,30,35	914582	137	915139	36,114
906016	84,101	909333	91,93	909426	88,93	911590	102	914250	22,30	914584	137	915310	100
909110	101	909334	91	909427	88,93	912220	77	914251	22,30	914630	134	915325	19
909120	59	909335	91	909430	87	912308	81	914252	22,30	914631	134	915330	19
909130	59	909336	91	909431	87,92	912310	81	914341	61	914632	134	909401P	91,93
909140	59	909337	91,93	909432	87,92	912311	81	914362	61	914633	134	909402P	91,93
909150	59	909339	91,93	909434	87	912314	81	914363	61	914671	136	909427P	88,93
909201	99	909340	91	909435	87,92	912324	81	914364	57	914676	127	912442B	133
909206	97	909341	91	909436	87,92	912328	81	914366	43	914677	127	912442C	133
909207	97	909348	98	909510	90,92	912330	81	914368	43	914678	127	912442G	133
909208	97	909350	87	909515	92	912410	46,136	914372	43	914679	127	912442R	133
909209	97,98	909355	87	909525	95	912415	46,136	914373	43	914680	127	912442Y	133
909211	97	909356	87	909530	94	912420	136	914374	43	914682	135	912447B	133
909212	84	909360	87,92	909531	94	912422	131	914375	43	914683	135	912447C	133
909217	100	909362	87,92	909532	94	912425	131	914376	106	914684	135	912447G	133
909218	100	909363	87,92	909535	89,94	912428	132	914377	106,113	914689	135	912447R	133
909220	100	909365	87,92	909560	98	912431	136	914378	106	914690	135	912447Y	133
909221	100	909366	87	910015	123	912435	132	914379	106	914691	135	914733B	132
909222	97	909367	87,90	910025	123	912436	133	914380	106	914692	135	914733G	132
909223	97	909368	98	910110	22	912437	133	914382	106,113	914720	133	914733R	132
909225	100	909369	89	910112	22	912441	132	914384	106	914722	135	914733Y	132
909226	100	909370	88	910115	22	912460	60,135	914386	106,113	914723	133		
909242	98	909371	88	910117	22	912461	60,135	914387	106	914731	134		
909247	98	909372	88	910140	114	912462	60,135	914388	106,113	914734	134		
909250	101	909373	88	910145	114	912464	60,135	914389	107	914735	134		
909251	101	909374	90	910150	114	912465	60,135	914393	61	914761	59		
909252	101	909375	89	910155	114	912467	135	914394	50	914800	128		
909275	96	909376	90	911182	71	913030	123	914395	50	914801	128		
909276	96	909377	90	911184	71	913100	128	914396	62	914802	128		
909280	96	909378	90	911366	98	913101	128	914429	30,41	914803	128		
909283	96	909380	89	911450	100	913102	128	914430	30,41	915060	114		
909284	96	909382	89	911452	100	913103	128	914441	20,28,39	915061	114		
909286	96	909384	89	911454	100	913145	107	914443	43	915062	114		
909290	96	909385	88	911456	100	913150	106,107	914444	43	915063	114		
909292	96	909386	88	911458	100	914000	137	914451	66	915064	114		

Reusable Face Shield



REUSABLE FACE SHIELD

ReOrder #	Description	Size	U/M
914332	Reusable Face Shield	Universal	12 Per Case
914333	Replacement Foam & Elastic Banding	Universal	24 Per Pack
914334	Reusable Face Shield (2) w/Extra Foam (8) & Banding (8)	N/A	Set



Visit Our Website skil-care.com for more information

Reusable Face Shield



REUSABLE FACE SHIELD

ReOrder #	Description	Size	U/M
914332	Reusable Face Shield	Universal	12 Per Case
914333	Replacement Foam & Elastic Banding	Universal	24 Per Pack
914334	Reusable Face Shield (2) w/Extra Foam (8) & Banding (8)	N/A	Set



Visit Our Website skil-care.com for more information

Skil-Care™ CORPORATION

29 Wells Ave., Yonkers, NY 10701

Please Route to:

- Administrator
- Director of Nursing
- Rehab Director
- Activities Director

Founded in 1978, Skil-Care Corporation has been developing and designing products for healthcare facilities worldwide.

29 Wells Ave
Yonkers, NY 10701

toll-free: 800.431.2972,

local: 914.963.2040

fax: 914.963.2567

online: www.skil-care.com

