



Re-Order Number	Descriptions	Size	Units
914481	Ultra Thin Gel Foam Wheelchair Cushion w/Low Shear II Cover	16" x 16"	Each
914482	Ultra Thin Gel Foam Wheelchair Cushion w/Low Shear II Cover	18" x 16"	Each
914483	Ultra Thin Gel Foam Wheelchair Cushion w/Low Shear II Cover	20" x 16"	Each

Purpose

This very low gel-foam cushion is especially effective for leg-propelling residents. It provides pressure relief comparable to any of Skil-Care's thin-line gel cushions.

Important Safety Feature

The Half Gel/Half Foam Cushion with a cloth outer cover has straps with buckles that must be secured to the wheelchair to prevent the cushion from sliding off the seat.

Easy Care

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. The Low-Shear outer cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.

This product is guaranteed for 1 year.