


Anti-Entrapment Bed Rail Pads

Re-Order #	Size	Unit
☐401090	Fits rails up to 37" x 15"	Pair



PURPOSE

Skil-Care's Anti-Entrapment Bed Rail Pads are designed to fit the new bed rails that have narrower spaces between the metal bars. Foam padded, they protect patients against contact with metal rails and guard against limb entrapment that may occur, even with the narrower bar spacing. A full 15 inches wide, these pads extend below the mattress to guard against entrapment between the rail and the mattress.

EASY TO USE

Place the pads over the rail and press the Velcro[®] strips together to secure the pads.

EASY CARE

The durable vinyl cover cleans easily with spray-type cleanser and a damp cloth. A mild disinfectant may be used.