


## ECONO-GEL PAD

Reorder #	Cover	Size
□751180	Bimini Blue	16" x 16" x 2"
□751181	Bimini Blue	18" x 16" x 2"
□751184	Bimini Pink	16" x 16" x 2"
□751185	Bimini Pink	18" x 16" x 2"



### Purpose

Colorful Sittin' Pretty Econo-Gel Pads are designed to relieve pressure and increase comfort for wheelchair users. High-density foam layers on top and bottom provide stability and cushioning. The gel-filled center pod adds pressure-relieving flotation.

### Important notice

The Sittin' Pretty Econo-Gel Pad has safety ties that must be secured to the wheelchair to prevent the cushion (and possibly the user) from slipping off the wheelchair seat.

### Easy care

Skil-Care's Sittin' Pretty Econo-Gel Pad is heat sealed to prevent urine and other fluids from seeping into the cushion. The vinyl cover may be cleaned with a spray-type cleanser. A mild disinfectant may be used.