

ELBOW PROTECTORS



Re-Order #	Previous #	Unit
<input type="checkbox"/> 503310	8010	Pair
<input type="checkbox"/> 503320	8010BP	36 Pr./Cs.

Effective protection against friction burns and decubitus ulcers on sensitive elbows. Designed to flex naturally with your patient's arm, this elbow protector provides the comfortable cushioning of plush synthetic sheepskin. Two hook-and-loop closures - one at the forearm and one above the elbow-keep the pad securely in place.

LAUNDERING INSTRUCTIONS

Skil-Care's Elbow Protectors may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended. To prolong the useful life of the hook-and-loop closures, secure them before laundering.