


FOOT ELEVATOR

Re-Order #	Previous #	Size	Unit
<input type="checkbox"/> 503100	8020	4" Width	Each
<input type="checkbox"/> 503140	8022	6" Width	Each



PURPOSE

Designed for effective d.u. prevention and management. The Skil-Care Foot Elevator raises your patient's heel off the bed to eliminate pressure and shearing and provide the air circulation necessary for accelerated tissue granulation and healing of existing decubiti.

EASY APPLICATION

Wrap Foot Elevator around patient's ankle and adjust to a snug, comfortable fit. Secure hook-and-loop closure. There's no need to wrap a liner around the patient's ankle; the soft flannelette material will not irritate patient's skin.

LAUNDERING

This device may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.