




## HIGH-DENSITY FOAM CUSHION

Reorder #	Cover Color	Size
□751188	Bimini Blue	18" x 16" x 2"
□751192	Bimini Pink	18" x 16" x 2"
□751189	Bimini Blue	18" x 16" x 3"
□751193	Bimini Pink	18" x 16" x 3"



### PURPOSE

Colorful Sittin' Pretty Foam Cushions reduce pressure and increase comfort for residents in wheelchairs, geri-chairs, lounge chairs and just about any other type of chair.

### IMPORTANT NOTICE

All Sittin' Pretty Cushion have safety ties that must be secured to the chair in order to prevent the cushion (and possibly the user) from slipping off the chair seat.

### EASY CARE

Skil-Care's Sittin' Pretty Foam Cushions are heat sealed to prevent urine and other fluids from seeping into the cushion. The vinyl covers may be cleaned with a spray-type cleanser. A mild disinfectant may be used.