


HIGH-DENSITY FOAM CUSHIONS

Reorder #	Cover	Cushion Size (inches)
<input type="checkbox"/> 753210	Nylon	16" x 18" x 2"
<input type="checkbox"/> 753230	Vinyl	16" x 18" x 2"
<input type="checkbox"/> 753240	No Cover	16" x 18" x 2"
<input type="checkbox"/> 753250	Nylon	16" x 18" x 3"
<input type="checkbox"/> 753270	Vinyl	16" x 18" x 3"
<input type="checkbox"/> 753280	No Cover	16" x 18" x 3"
<input type="checkbox"/> 753285	Nylon	16" x 18" x 4"
<input type="checkbox"/> 753290	Vinyl	16" x 18" x 4"
<input type="checkbox"/> 753295	No Cover	16" x 18" x 4"
<input type="checkbox"/>		

Maximum cushioning with minimum compression make Skil-Care High Density Foam Cushions your best choice for wheelchair and geri-chair comfort.



COVER CARE

Nylon and vinyl covers may be cleaned with a spray-and-wipe cleaner and damp cloth. Laundering the nylon cover will cause it to lose its water repellent quality. If the nylon cover is laundered, do not wash or dry it at temperatures over 180°F (82°C). Bleaching is not recommended.