


## PILLOW PROP



Reorder #	Size	Unit
<input type="checkbox"/> 905010	13 1/2" x 13 1/2"	Each
<input type="checkbox"/> 905015	15" x 20"	Each

### Purpose

The Pillow Prop may be used to elevate the arm of a patient with edema. It can also serve as an effective positioning and support device for a side-lying patient. Filled with high-density foam, the Pillow Prop will not "cave-in" under the weight on an arm or the patient's body weight.

### Note

For patient comfort, it is suggested that the Pillow Prop be covered with a pillowcase.

### Easy Care

The fluid-proof, anti-bacterial nylon cover wipes clean with a damp cloth and spray-type cleanser. Do not launder.