


PRESSURE-CHECK[™] FOAM CUSHION

Reorder #	Size
□753180	18" x 16" x 3"



Effective Design

Separate foam sections adjust individually to eliminate pressure-causing high and low spots. The base are below the sections prevents bottoming out and provides maximum weight distribution for relief reduction and comfort.

New Low-Shear II Cover

The cover is two-way stretch Low-Shear II fabric. It moves with the wheelchair user to reduce skin-damaging shear and friction.

Important Safety Features

The Pressure-Check[™] Cushion has safety strap that secures to the wheelchair to prevent the cushion from sliding off the wheelchair. It also features an anti-slip bottom for added safety.

Easy Care

The fluid-proof Low-Shear II cover wipes clean with spray-and-wipe cleanser. It may be disinfected with a mild, non-chlorine disinfectant.