


Visit us online  
at Skil-Care.com

**Skil-Care**<sup>™</sup>  
CORPORATION

29 Wells Avenue  
Yonkers, NY 10701

914-963-2040  
800-431-2972



## STAY-ON HEEL PROTECTOR

Re-Order #	Previous #	Unit
<input type="checkbox"/> 503050	8002	Pair
<input type="checkbox"/> 503060	_____	Dz./Pr.

When keeping heel protectors on your patient is a problem, the Stay-On Heel Protector is the solution. The oversized design plus the slide-lock buckle closure make this pad difficult to kick off. Made from a double layer of plush synthetic sheepskin for superior pressure relief and anti-decubitus protection.

### LAUNDERING

The Stay-On Heel Protects may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommend.

**Skil-Care**<sup>™</sup>  
CORPORATION

*Advanced design for nursing and therapy.*<sup>™</sup>