

# Thin-Line™


## ECONO-GEL PAD

Reorder #	Cover	Size
□751170	Vinyl w/Polyester	16" x 16" x 1½"
□751171	Vinyl w/Polyester	18" x 16" x 1½"
□751172	Vinyl Only	16" x 16" x 1½"
□751173	Vinyl Only	18" x 16" x 1½"

**Skil-Care's new Thin-Line™ Cushions feature a low-profile design that enables most smaller wheelchair users to reach the floor or footrests.**

• **Polyester outer cover** reduces skin shear. Two-tone fabric serves as a reminder to reposition patient. Removes for laundering.

• **Fluid-proof vinyl cover** protects against incontinence. Wipes clean for easy care.



• **Safety ties** keep the cushion in place.

• **High-density foam layers** on top and bottom increase comfort and stability.

• **Gel-filled pod** between foam layers provides pressure-relieving flotation

### Important Safety Reminder

Each Econo-Gel Pad has safety ties that must be secured to the wheelchair frame to prevent the cushion from slipping off the seat.

### Easy care

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. The polyester cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.