


TRIPLE-PLY HEEL PROTECTORS



<u>Re-Order #</u>	<u>Unit</u>
-------------------	-------------

<input type="checkbox"/> 503010	Pair
---------------------------------	------

<input type="checkbox"/> 503020	36 Pr./Cs.
---------------------------------	------------

A uniquely designed heel protector that places three layers of plush synthetic sheepskin directly under the critical pressure point of your patient's heel. Maintains maximum resiliency for constant pressure relief.

LAUNDERING INSTRUCTIONS

Wash and dry at temperatures up to 180°F (82°C). Bleaching is not recommended. To prolong the useful life of the hook-and-loop secure closure before laundering.